



# Agilent 1290 Infinity Binary Pump VL

User Manual



**Agilent Technologies**

# Notices

© Agilent Technologies, Inc. 2009,  
2010-2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

## Manual Part Number

G4220-90005 Rev. C

## Edition

08/2012

Printed in Germany

Agilent Technologies  
Hewlett-Packard-Strasse 8  
76337 Waldbronn

**This product may be used as a component of an in vitro diagnostic system if the system is registered with the appropriate authorities and complies with the relevant regulations. Otherwise, it is intended only for general laboratory use.**

## Warranty

**The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.**

## Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

## Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will

receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

## Safety Notices

### CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

### WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

# In This Guide...

This manual covers the Agilent 1290 Infinity Binary Pump (G4220B).

## **1 Introduction**

This chapter gives an introduction to the pump, instrument overview and internal connectors.

## **2 Site Requirements and Specifications**

This chapter provides information on environmental requirements, physical and performance specifications.

## **3 Installing the Module**

This chapter gives information about the preferred stack setup for your system and the installation of your Agilent 1290 Infinity Binary Pump.

## **4 Using the Pump**

This chapter explains the operational parameters of the Agilent 1290 Infinity Binary Pump.

## **5 How to Optimize the Performance of Your Module**

This chapter gives hints on how to optimize the performance or use additional devices.

## **6 Troubleshooting and Diagnostics**

Overview about the troubleshooting and diagnostic features.

## **7 Error Information**

This chapter describes the meaning of error messages, and provides information on probable causes and suggested actions how to recover from error conditions.

**8 Test Functions and Calibrations**

This chapter describes the tests for the module.

**9 Maintenance**

This chapter describes the maintenance of the Agilent 1290 Infinity Binary Pump.

**10 Parts and Materials for Maintenance**

This chapter provides information on parts for maintenance.

**11 Identifying Cables**

This chapter summarizes information on all cables.

**12 Hardware Information**

This chapter describes the pump in more detail on hardware and electronics.

**13 LAN Configuration**

This chapter provides information on connecting the detector to the Agilent ChemStation PC.

**14 Appendix**

This chapter provides addition information on safety, legal and web.

# Contents

<b>1</b>	<b>Introduction</b>	<b>9</b>
	Features	10
	Overview of the Binary Pump	11
	Pump Principle	12
<b>2</b>	<b>Site Requirements and Specifications</b>	<b>15</b>
	Site Requirements	16
	Physical Specifications	19
	Specifications	20
<b>3</b>	<b>Installing the Module</b>	<b>23</b>
	Unpacking the Module	24
	Optimizing the Stack Configuration	26
	Removing the Transport Foam	31
	Installing the Pump	32
	Flow connections to the pump	35
	Installation of seal wash option	39
<b>4</b>	<b>Using the Pump</b>	<b>41</b>
	Preparing the Binary Pump	42
	Algae Growth in HPLC Systems	43
	Setting up the Pump with the Instrument Control Interface	44
	Priming the Pump	54
	Solvent Information	56
<b>5</b>	<b>How to Optimize the Performance of Your Module</b>	<b>61</b>
	Delay Volume and Extra-Column Volume	62
	How to Configure the Optimum Delay Volume	63
	How to Achieve Higher Resolution	65
	Using Solvent Calibration Tables	68

<b>6</b>	<b>Troubleshooting and Diagnostics</b>	<b>69</b>
	Overview of the Module's Indicators and Test Functions	70
	Status indicators	71
	Available Tests vs User Interfaces	73
	Agilent Lab Advisor Software	74
<b>7</b>	<b>Error Information</b>	<b>75</b>
	What Are Error Messages	77
	General Error Messages	78
	Pump Error Messages	85
<b>8</b>	<b>Test Functions and Calibrations</b>	<b>101</b>
	Introduction	102
	System Pressure Test	103
	Pump Head Leak Test	105
<b>9</b>	<b>Maintenance</b>	<b>109</b>
	Introduction to Maintenance	110
	Warnings and Cautions	112
	Overview of Maintenance	113
	Cleaning the Module	114
	Installing Fittings and Capillaries	115
	Replacing the Shutoff Valve Panel	116
	Replacing the Pressure Sensor	118
	Replacing the Inlet valve	120
	Replacing the Outlet valve	122
	Replacing the Solvent Selection Valve (SSV)	125
	Changing configuration or replacing the Jet Weaver	127
	Replacing the Seal Wash Pump	130
	Releasing a Stuck Inlet Valve	131
	Replacing the Pump Heads	134
	Disassembling the Pump Head	139
	Disassembling the Primary Pump Head	140
	Disassembling the Secondary Pump Head	143
	Assembling the Pump Head	146
	Replacing the Purge Valve Head	154
	Replacing Parts of the High Pressure Filter Assembly	157

Installing the Valve Rail Kit	159
Replacing the Main Power Fuses	160
Replacing Module Firmware	162
Preparing the Pump Module for Transport	163

## **10 Parts and Materials for Maintenance 167**

Overview of Maintenance Parts	168
Capillaries	169
Solvent Cabinet Kit	170
Seal Wash Option	171
Pump Head Assembly Parts	172
Primary Pump Head Parts	174
Secondary Pump Head Parts	178
Purge Valve	182
Cover Parts	183
Leak Parts	184
Fuses	184
Accessory Kit	185
Others	186

## **11 Identifying Cables 189**

Cable Overview	190
Analog cables	192
Remote Cables	194
BCD Cables	197
CAN/LAN Cable	199
RS-232 Cable Kit	200
Agilent 1200 Module to Printer	201

## **12 Hardware Information 203**

Firmware Description	204
Electrical Connections	207
Interfaces	209
Setting the 8-bit Configuration Switch	216
Instrument Layout	220
Early Maintenance Feedback	221

### **13 LAN Configuration 223**

- What You Have to Do First 224
- TCP/IP parameter configuration 226
- Configuration Switch 227
- Initialization mode selection 228
- Dynamic Host Configuration Protocol (DHCP) 232
- Link configuration selection 236
- Automatic configuration with Bootp 237
- Manual Configuration 247
- PC and User Interface Software Setup Setup 253

### **14 Appendix 255**

- General Safety Information 256
- The Waste Electrical and Electronic Equipment (WEEE) Directive (2002-96-EC) 259
- Radio Interference 260
- Sound Emission 261
- Agilent Technologies on Internet 262





# 1

## Introduction

Features 10

Overview of the Binary Pump 11

Pump Principle 12

This chapter gives an introduction to the pump, instrument overview and internal connectors.



## Features

The Binary pump is designed for highest performance, GLP compliance and easy maintenance. It includes the following features:

- Seal wash for continued high lifetime of pump seals using high buffer concentrations.
- Solvent selection valve for method flexibility.
- Jet Weaver for optimum mixing performance with a minimum of delay volume.
- Automatic purge function for ease of use and unattended preparation of the system.
- Auto tuning of the delivery cycle for compensation of elasticity and dead volume effects.
- Solvent selection for optimum density correction.
- Fast defill function to increase robustness of the pump.
- Two pistons in series design for increased reliability.
- Smooth motion control to prevent shock movements.

For specifications, see [Table 2](#) on page 20.

### NOTE

This Binary pump has been introduced together with the Agilent 1290 Infinity Liquid Chromatograph.

---

## Overview of the Binary Pump


The Agilent 1290 Infinity Binary Pump comprises two identical pumps integrated into one housing. Binary gradients are created by high-pressure mixing. A degassing unit is included for applications that require best flow stability, especially at low flow rates, for maximum detector sensitivity. The flow path of the pump has been optimized for minimal delay of gradients. Typical applications are high throughput methods with fast gradients on high resolution 2.1 mm columns. The pump is capable of delivering flow in the range of 0.05 -2 mL/min against up to 1200 bar. A solvent selection valve allows forming binary mixtures (isocratic or gradient) from one of two solvents per channel. Active seal wash (optional) is available for use with concentrated buffer solutions.

## Pump Principle

The Binary Pump is based on a two-channel, dual-piston in-series design which comprises all essential functions that a solvent delivery system has to fulfill. Metering of solvent and delivery to the high-pressure side are performed by two pump assemblies which can generate pressure up to 1200 bar.

Each channel comprises a pump assembly including independent pump drive for each piston, pump head, inlet valve, outlet valve, solvent heat exchanger and an outlet filter. The two channels are fed into a low-volume mixing groove in an automatic purge valve and a Jet Weaver mixer, with 35 or 100 µl volume can be added downstream for optimum mixing performance.

A system pressure sensor, for monitoring the pump pressure, is attached to the purge valve, normally connected in the B-channel of the pump, before the mixing groove, in order to minimize delay volumes.




**Figure 1** The hydraulic path

In the user interface the solvent in use for each channel can be selected in order to get optimum flow and composition accuracy. Although selecting the right solvent is not required for low ripple, especially for composition accuracy across pressure it is vital for the compensation algorithm to use the right solvent properties. The density of the solvents is increased under the influence of pressure and a certain displaced solvent will expand again when the pressure is released, for example across the column. In order to achieve

the correct volumetric flow while sample passes through the detector it is necessary to correct for density related flow inaccuracies in the pump module.

In order to always deliver the best possible pressure stability, the pump constantly tunes the delivery cycle for elasticity and dead volume effects. With this feature the pump is able to deliver a stable and accurate flow without requiring individual calibration settings. A further feature of the control and compensation algorithm is leak correction. With this it is even possible to compensate for minor leaks in primary pump chamber (inlet valve and seal), without the performance of the pump being affected.

To increase the robustness of the pump it uses a fast defill function which reduces the delivery time of the primary piston, thereby reducing the net effects of leaks considerably. Furthermore each pump channel has only two valves on its two pump heads which also reduce the potential of failures. In order to reduce stresses on the pump hardware, it uses a smooth motion control, which slowly increases or decreases the speed of the pistons to prevent shock movements. In order to be able to control these movements a high resolution encoder unit is attached to the pump drives which resolves a revolution into 65000 steps, and each step translates to a volume of about 300 pl.



## 2

# Site Requirements and Specifications

Site Requirements 16

Physical Specifications 19

Specifications 20

This chapter provides information on environmental requirements, physical and performance specifications.



## Site Requirements

A suitable environment is important to ensure optimal performance of the instrument.

### Power Consideration

The module power supply has wide ranging capabilities and accepts any line voltage in the range mentioned in [Table 1](#) on page 19. Consequently, there is no voltage selector in the rear of the module. There are also no externally accessible fuses, because automatic electronic fuses are implemented in the power supply.

#### **WARNING**

**Module is partially energized when switched off, as long as the power cord is plugged in.**

**Repair work at the module can lead to personal injuries, e.g. shock hazard, when the cover is opened and the module is connected to power.**

- Make sure that it is always possible to access the power plug.
- Remove the power cable from the instrument before opening the cover.
- Do not connect the power cable to the Instrument while the covers are removed.

#### **WARNING**

**Incorrect line voltage at the module**

**Shock hazard or damage of your instrument can result if the devices are connected to line voltage higher than specified.**

- Connect your module to the specified line voltage.



**CAUTION**

Inaccessible power plug.

In case of emergency it must be possible to disconnect the instrument from the power line at any time.

- Make sure the power connector of the instrument can be easily reached and unplugged.
  - Provide sufficient space behind the power socket of the instrument to unplug the cable.
- 

## Power Cords

Different power cords are offered as options with the module. The female end of all power cords is identical. It plugs into the power-input socket at the rear. The male end of each power cord is different and designed to match the wall socket of a particular country or region.

**WARNING**

**Absence of ground connection or use of unspecified power cord**

**The absence of ground connection or the use of unspecified power cord can lead to electric shock or short circuit.**

- Never operate your instrumentation from a power outlet that has no ground connection.
  - Never use a power cord other than the Agilent Technologies power cord designed for your region.
- 

**WARNING**

**Use of unsupplied cables**

**Using cables not supplied by Agilent Technologies can lead to damage of the electronic components or personal injury.**

- Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.
-

#### **WARNING**

##### **Unintended use of supplied power cords**

**Using power cords for unintended purposes can lead to personal injury or damage of electronic equipment.**

- Never use the power cords that Agilent Technologies supplies with this instrument for any other equipment.
- 

## **Bench Space**

The module dimensions and weight (see [Table 1](#) on page 19) allow you to place the module on almost any desk or laboratory bench. It needs an additional 2.5 cm (1.0 inches) of space on either side and approximately 8 cm (3.1 inches) in the rear for air circulation and electric connections.

If the bench shall carry a complete HPLC system, make sure that the bench is designed to bear the weight of all modules.

The module should be operated in a horizontal position.

## **Condensation**

#### **CAUTION**

Condensation within the module

Condensation will damage the system electronics.

- Do not store, ship or use your module under conditions where temperature fluctuations could cause condensation within the module.
  - If your module was shipped in cold weather, leave it in its box and allow it to warm slowly to room temperature to avoid condensation.
-

## Physical Specifications

**Table 1** Physical Specifications

Type	Specification	Comments
Weight	21.8 kg (48 lbs)	
Dimensions (height × width × depth)	240 x 345 x 435 mm (9.3 x 13.5 x 17 inches)	
Line voltage	100 – 240 VAC, ± 10 %	Wide-ranging capability
Line frequency	50 or 60 Hz, ± 5 %	
Power consumption	350 VA / 270 W / 922 BTU	Maximum
Ambient operating temperature	4–55 °C (41–131 °F)	
Ambient non-operating temperature	-40 – 70 °C (-4 – 158 °F)	
Humidity	< 95 % r.h. at 40 °C (104 °F)	Non-condensing
Operating altitude	Up to 2000 m (6562 ft)	
Non-operating altitude	Up to 4600 m (15091 ft)	For storing the module
Safety standards: IEC, CSA, UL	Installation category II, Pollution degree 2	For indoor use only.

## Specifications

**Table 2** Performance specifications

Type	Specification	Comments
Hydraulic system	Two dual pistons in series pumps with proprietary servo-controlled variable stroke design and smooth motion control.	
Settable flow range	Setpoints 0.001 — 2 ml/min, in 0.001 ml/min increments.	Executed in 300 pl/step increments
Flow precision	≤0.07 % RSD or 0.005 min SD, whatever is greater (0.2 — 2.0 ml/min).	Based on retention time at constant room temperature.
Flow accuracy	±1 % or 10 µl/min, whatever is greater.	Pumping degassed H <sub>2</sub> O at 10 MPa (100 bar)
Pressure operating range	Operating pressure range 120 MPa (1200 bar), up to 2 ml/min.	
Pressure pulsation	<1 % amplitude or < 0.5 MPa (5 bar), whatever is greater.	At 1 ml/min water
Compressibility compensation	Automatic, pre-defined, based on mobile phase selection.	
Gradient formation	High pressure binary mixing.	
Delay volume	Jet Weaver V35: <45 µl Jet Weaver V100: <75 µl	JetWeaver generally recommended. For applications requiring lowest delay volumes, JetWeaver can be removed out of the flow path. Delay volume is then solely determined by the volume of the connection capillary.
Composition range	Settable range: 0 — 100 %	Recommended range: 1 — 99 % or 5 µl/min per channel, whatever is greater.

**Table 2** Performance specifications

Type	Specification	Comments
Composition precision	<0.15 % RSD, or 0.01 min SD, whatever is greater.	0.2 —5.0 ml/min; based on retention time at constant room temperature
Composition accuracy	±0.35 % absolute (5 – 95 %, 0.2 – 5.0 ml/min)	Water/caffeine tracer
Solvent selection valve	Default	Standard part of the pump
Integrated degassing unit	Number of channels: 2 Internal volume per channel: 1.5 mL	
Control	Agilent ChemStation for LC (B.04.02 or above) EZChrom Elite (3.3.2 SP1 or above) OpenLAB (3.3.2 SP3) Masshunter (B.02.01 SP1 or above)	
Local control	Agilent Instant Pilot (G4208A) (B.02.08 or above)	
Communications	Controller-area network (CAN), RS232C, APG remote: ready, start, stop and shutdown signals, LAN	
Safety and maintenance	Extensive diagnostics, error detection and display (through Agilent Lab Advisor), leak detection, safe leak handling, leak output signal for shutdown of the pumping system. Low voltage in major maintenance areas.	
GLP features	Early maintenance feedback (EMF) for continuous tracking of instrument usage in terms of purge valve switches and volume of pumped mobile phase with pre-defined and user settable limits and feedback messages. Electronic records of maintenance and errors.	
Housing	All materials recyclable.	

## **2 Site Requirements and Specifications**

### **Specifications**



## 3 Installing the Module

Unpacking the Module	24
Optimizing the Stack Configuration	26
One Stack Configuration	26
Two Stack Configuration	29
Removing the Transport Foam	31
Installing the Pump	32
Flow connections to the pump	35
Installation of seal wash option	39

This chapter gives information about the preferred stack setup for your system and the installation of your Agilent 1290 Infinity Binary Pump.



## Unpacking the Module

### Damaged Packaging

If the delivery packaging shows signs of external damage, please call your Agilent Technologies sales and service office immediately. Inform your service representative that the instrument may have been damaged during shipment.

#### CAUTION

"Defective on arrival" problems

If there are signs of damage, please do not attempt to install the module. Inspection by Agilent is required to evaluate if the instrument is in good condition or damaged.

- Notify your Agilent sales and service office about the damage.
  - An Agilent service representative will inspect the instrument at your site and initiate appropriate actions.
-



## Delivery Checklist

Ensure that all parts and materials have been delivered with your module. The delivery checklist is shown below. For parts identification please check the illustrated part lists in [“Parts and Materials for Maintenance”](#) on page 167. Please report any missing or damaged parts to your local Agilent Technologies sales and service office.

**Table 3** Binary Pump Checklist

Description	Quantity
1290 Infinity Bin Pump User Manual	1
1290 Infinity System Manual	1
Agilent Lab Advisor	1
LC HW User Information & Utilities DVD	1
Power Cord	1
Solvent Cabinet Kit with 4 bottles	1
Accessory Kit (see <a href="#">“Accessory Kit”</a> on page 185)	1
RRHD Eclipse Plus C18, 2.1x50 mm, 1.8 $\mu$	1

## Optimizing the Stack Configuration


If your module is part of a complete Agilent 1290 Infinity Liquid Chromatograph, you can ensure optimum performance by installing the following configurations. These configurations optimize the system flow path, ensuring minimum delay volume.

For other possible configurations, please refer to the Agilent 1290 Infinity System Manual.

### One Stack Configuration

Ensure optimum performance by installing the modules of the Agilent 1290 Infinity Binary LC System in the following configuration (See [Figure 2](#) on page 27 and [Figure 3](#) on page 28). This configuration optimizes the flow path for minimum delay volume and minimizes the bench space required.


The Agilent 1290 Infinity Binary Pump should always be installed at the bottom of the stack.



**Figure 2** Recommended stack configuration for 1290 Infinity with binary pump (front view)

### 3 Installing the Module


#### Optimizing the Stack Configuration



**Figure 3** Recommended stack configuration 1290 Infinity with binary pump (rear view)

## Two Stack Configuration


In case the autosampler thermostat is added to the system, a two-stack configuration is recommended, which places both heavy modules (1290 Infinity pump and thermostat) at the bottom of each stack and avoids high stacks. Some users prefer the lower height of this arrangement even without the autosampler thermostat. A slightly longer capillary is required between the pump and autosampler. (See [Figure 4](#) on page 29 and [Figure 5](#) on page 30).



**Figure 4** Recommended two stack configuration for 1290 Infinity with binary pump (front view)

### 3 Installing the Module


#### Optimizing the Stack Configuration




**Figure 5** Recommended two stack configuration for 1290 Infinity with binary pump (rear view)

## Removing the Transport Foam


1 Open the front cover of the module.



2 Carefully remove the protective foam.



3 Close the front cover.




# Installing the Pump

Parts required	#	Description
	1	Pump
	1	Power cord
	1	Agilent Control Software and/or Instant Pilot G4208

**Preparations**

- Locate bench space
- Provide power connections
- Unpack the pump


- 1 Place the module on the bench in a horizontal position.
- 2 Ensure the power switch on the front of the module is OFF (switch stands out).



**Figure 6**    Front of Binary Pump



- 3 Connect the power cable to the power connector at the back of the 1290 Infinity Binary Pump.
- 4 Connect the required interface cables to the rear of the 1290 Infinity Binary Pump.



**Figure 7** Rear of Binary Pump

**NOTE**

In an Agilent 1290 Infinity System, the individual modules are connected by CAN cables. An Agilent 1200 Series Instant Pilot can be connected to the CAN bus of any module. Connection to an Agilent data system is established through the built-in LAN port of the detector. The LAN port of the detector must be used as the detector generates the highest data rate of all modules. For more information about connecting the Instant Pilot or Agilent Data System, please refer to the respective user manual. For setting up the LAN access, see “LAN Configuration” on page 223.

### 3 Installing the Module

#### Installing the Pump

- 5 Turn on the power by pushing the button at the lower left hand side of the module.

The power button stays pressed in and the status LED should be green.

---

**NOTE**

When the line power button stands out and the green light is off, the module is turned off.


---

**NOTE**

The module was shipped with default configuration settings. For changing these settings, refer to section *Setting the 8-bit configuration switch*.

---

## Flow connections to the pump



The pump is shipped with tubing and capillary connections installed between solvent selection valve, degassing unit, pump heads, pressure sensor, purge valve and Jet Weaver. This section describes the installation of additional flow connections.

### 3 Installing the Module

Flow connections to the pump

Parts required	#	p/n	Description
	1		Other modules
	1	G4220-68705	Accessory kit
	1	G5067-1531	Solvent Cabinet Kit

**Preparations** Pump is installed in the LC system.


#### **WARNING**

**Toxic, flammable and hazardous solvents, samples and reagents**

**The handling of solvents, samples and reagents can hold health and safety risks.**


- When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.
- The volume of substances should be reduced to the minimum required for the analysis.
- Do not operate the instrument in an explosive atmosphere.

- 1 Remove the front cover by pressing the snap fasteners on both sides.




**Figure 8** Removing the Front Cover

- 2 Place the solvent cabinet on top of the module stack that includes the 1290 Infinity Binary Pump.
- 3 Put the four bottles into the solvent cabinet and screw a bottle head assembly onto each bottle.
- 4 Install the shutoff valve panel at the top left corner of the instrument.



- 5 Connect the solvent tubes from the bottle head assemblies to the inlet connectors A1, A2, B1 and B2 of the solvent selection valves. Use the brown bottle for the aqueous solvent (usually channel A1).




### **3** **Installing the Module**

**Flow connections to the pump**

- 6** Label the tubes accordingly using the supplied stickers and fix the tubes in the clips of solvent cabinet and 1290 Infinity Binary Pump.
- 7** Connect the outlet of the Jet Weaver to the autosampler.
- 8** Connect Waste tubes (G4220-67000) to the purge valve outlets at ports 5 and 6.
- 9** Connect the corrugated waste tube to the outlet of the leak panel.
- 10** Route the corrugated waste tube to a waste container.
- 11** Route drain tubes coming from modules on top of the pump through the pump.
- 12** Purge your system prior to the first use (see [“Priming the Pump”](#) on page 54).

## Installation of seal wash option



**Figure 9** Binary Pump with Seal Wash Function

The 1290 Infinity Binary Pump is optionally available with a seal wash function. This option is recommended when using buffers or other non-volatile solvents or additives that could deposit on pistons and seals. It is used for regularly cleaning these parts automatically.

- 1 Place a wash solvent reservoir into the solvent cabinet. A mixture of distilled water and isopropanol (90/10) is a good choice for many applications.
- 2 Put the solvent inlet tube into the solvent reservoir, close it and connect the tube to the seal wash pump.
- 3 Route the outlet of the wash tube into a waste container.

### **3 Installing the Module**

#### Installation of seal wash option





## 4 Using the Pump

Preparing the Binary Pump	42
Algae Growth in HPLC Systems	43
How to Prevent and-or Reduce the Algae Problem	43
Setting up the Pump with the Instrument Control Interface	44
Overview	44
Instrument Configuration	45
The Pump User Interface (Dashboard Panel)	46
Control Settings	48
Method Parameter Settings	50
Priming the Pump	54
Solvent Information	56

This chapter explains the operational parameters of the Agilent 1290 Infinity Binary Pump.



## Preparing the Binary Pump

For best performance of the pump:

- Place solvent cabinet with the solvent bottles always on top (or at a higher level) of the pump.
- For optimum performance, use the built-in degasser. This is mandatory for flow rates below 0.5 mL/min and for configurations without Jet Weaver.
- When using the pump with vacuum degassing unit, flush the degassing unit with at least 5 mL per channel before operating the pump, especially when the pumping system had been turned off for a certain length of time (for example, overnight) and volatile solvent mixtures are used in the channels.
- Prevent blocking of solvent inlet filters (never use the pump without solvent inlet filters). Growth of algae should be avoided, see [“Algae Growth in HPLC Systems”](#) on page 43.
- Check pump outlet filters and column frit in regular time intervals. A blocked pump outlet filter can be identified by black, yellow or greenish layers on its surface.
- Whenever possible use a minimum flow rate of 5  $\mu$ L/min per solvent channel to avoid crossflow of solvent into the unused pump channel.
- When using buffer solutions, flush the system with water before switching it off.
- The optional seal wash function should be used when buffer solutions with concentrations of 0.1 M or higher are being pumped for long periods of time.
- Never leave an unused pump with water in a channel for an extended period of time (2-3 days). Always flush with organic solvent or add 10 % isopropanol to water.

## Algae Growth in HPLC Systems

The presence of algae in HPLC systems can cause a variety of problems that may be incorrectly diagnosed as instrument or application problems. Algae grow in aqueous media, preferably in a pH range of 4-8. Their growth is accelerated by buffers, for example phosphate or acetate. Since algae grow through photosynthesis, light will also stimulate their growth. Even in distilled water small-sized algae grow after some time.

### Instrumental Problems Associated With Algae

Algae deposit and grow everywhere within the HPLC system causing:

- Blocked solvent filters or deposits on inlet or outlet valves resulting in unstable flow, composition or gradient problems or a complete failure of the pump.
- Small pore high pressure solvent filters, usually placed before the injector to plug resulting in high system pressure.
- PTFE frits blockage leading to increased system pressure.
- Column filters to plug giving high system pressure.
- Flow cell windows of detectors to become dirty resulting in higher noise levels (since the detector is the last module in the flow path, this problem is less common).

### How to Prevent and-or Reduce the Algae Problem

- Always use freshly prepared solvents, especially use demineralized water which was filtered through about 0.2  $\mu\text{m}$  filters.
- Never leave mobile phase in the instrument for several days without flow.
- Always discard old mobile phase.
- Use the amber solvent bottle (Solvent bottle, amber (9301-1450)) supplied with the instrument for your aqueous mobile phase.
- If possible add a few mg/l sodium azide or a few percent organic solvent to the aqueous mobile phase.

## **Setting up the Pump with the Instrument Control Interface**

### **Overview**

Parameters described in following sections is offered by the instrument control interface and can usually be accessed through Agilent instrument control software. For details, please refer to manuals and online help of respective user interfaces.

## Instrument Configuration

Use the **Instrument Configuration** dialog box to examine and, if necessary, modify your instrument configuration. The **Configurable Modules** panel contains a list of all modules available for configuration. The **Selected Modules** panel contains the list of configured modules.

**Auto Configuration:** Under **Communication settings**, select either the **Host Name** option or the **IP address** option and enter the appropriate value for the host computer to enable automatic detection of the hardware configuration. The system configures the instrument automatically with no further manual configuration necessary.

The Binary Pump configuration parameters are in two sections:

- **Communication**
- **Options**

**Communication:** The parameters in this dialog box are detected automatically during autoconfiguration.

- **Device name,**
- **Type ID,**
- **Serial number,**
- **Firmware revision,**
- Button **Connection settings**

**Options:**

- **Pressure Unit:**  
select the pressure units from the drop-down list (bar, psi or MPa).
- **Seal wash installed:**  
This check box is marked to indicate that an optional seal wash has been detected during autoconfiguration.
- **ISET installed:**  
This check box is marked to indicate that ISET is installed. Click **ISET Configurations** to open the **ISET Configuration** dialog box, which allows you to configure a sampler for the ISET emulation.

**Configure Solvent Type Catalogs:** Displays the **Solvent Type Catalogs** dialog box, which allows you to import and export solvent calibration data. See [“Importing Solvent Calibration Tables”](#) on page 68.

Please refer to the online help of your user interface for more detailed information.

## The Pump User Interface (Dashboard Panel)

### Module Graphic

The items in the pump graphic have the following meaning and function:



Indicates that an External Contacts board is installed.




The level of solvent in the bottle is denoted by the green area; when the solvent level falls below the specified volume, the area turns yellow; when the bottle is empty, the area turns red. Clicking on the solvent bottle displays the **Bottle Fillings** dialog box. The tooltip for the bottle shows the solvent name.



Indicates that the ISET option is installed but with no active method (gray) or installed and active (orange). (G4220A only)



Indicates the presence of a solvent selection valve. Click the graphic to switch the valve; the animation shows when the valve is switched.



The pressure setpoints. The red line shows the current maximum pressure limit; the green area shows the current pressure (also shown as text).

The current solvent flow rate (in mL/min) is displayed above the pressure display.

## Instrument Signals

The following pump signals are displayed:

<b>Flow</b>	The current solvent flow rate (in mL/min).
<b>Pressure</b>	The current pump pressure (in bar, psi or MPa, see <a href="#">“Instrument Configuration”</a> on page 45).
<b>Pressure Limit</b>	The current maximum pressure limit.
<b>Composition A:B</b>	The current solvent composition. When a solvent selection valve is fitted, the channels are shown in the graphic.

## Context Menu

The context menu of the dashboard panel contains the following commands:



**Control** Displays the pump's **Control** dialog box.



**Method** Displays the pump's **Method Setup** dialog box.

**Set Error Method** Sets the method that is loaded if an error occurs to the method that is currently available in the hardware.

**Identify Device** Causes the LED on the front of the module to blink for a few seconds.



**Switch Pump On/Off** Toggles the status of the pump, on or off.

**Bottle Fillings** Displays the **Bottle Fillings** dialog box.

**Purge On/Off** Allows you to control the purging of the system.

**Prime On/Off** Allows you to prime the pump heads for initially drawing solvent.

**Conditioning On/Off** Allows you to switch pump conditioning on and off. The conditioning function is useful for removing small air bubbles inside the pump flow path.

# Control Settings

The Binary Pump control parameters are in six sections:

- **Pump**
- **Seal Wash**
- **Automatic Turn On**
- **Purge**
- **Prime**
- **Conditioning**

**Table 4**    Pump control parameters

Parameter	Limits	Description
<b>Pump</b>		Enables you to switch the pump <b>On</b> , <b>Off</b> or to a <b>Standby</b> condition. In the <b>Standby</b> condition, the pump motor is still active, and when the pump is switched on again, does not need to be re-initialized.
<b>Seal Wash</b>		<p>The seal wash can be set up to be run once (<b>Single wash</b>) or periodically (<b>Periodic</b>).</p> <ul style="list-style-type: none"><li>• <b>Off</b>: no seal wash is used.</li><li>• <b>Single wash</b>: the seal will be purged for a specified time.</li><li>• <b>Periodic</b>: a periodic wash will be applied for a defined period in minutes.</li></ul> <p>The option is available only when the pump has seal wash capability. The seal wash capability is detected by the module itself. If seal wash is installed, it is recommended to use it in order to increase the primary seal lifetime.</p>
<b>Seal Wash Run Mode</b>		<p>Allows you to define when to use the seal wash:</p> <ul style="list-style-type: none"><li>• <b>Off</b>: The seal wash is inactive.</li><li>• <b>On when pump is on</b>: The seal wash is active only when the pump is on.</li><li>• <b>On all the time</b>: The seal wash is active when the pump is on or in standby mode.</li></ul>
<b>Automatic Turn On</b>		Module can be turned on at a specified date/time. This feature can only be used if the module power switch is turned on.



**Table 4** Pump control parameters

Parameter	Limits	Description
<b>Purge</b>	<b>Time:</b> 0 – 100.00 min in steps of 0.01 . <b>Flow:</b> 0.000 – 5.000 mL/min for each channel, in steps of 0.001 (10.000 mL/min maximum).	<p>Setup and activation of <b>Purge</b> parameters. The automatic purge valve can be used for purging the system. The process has been automated for ease of use.</p> <ul style="list-style-type: none"> <li>• <b>Off:</b> Turns off the purge.</li> <li>• <b>On:</b> The device is purged.</li> <li>• <b>Purge Flow, Time</b> and <b>Composition</b> during purge have to be defined.</li> </ul> <p>As soon as the duration time of the purge ends, the module automatically switches to analytical conditions again.</p>
<b>Prime</b>		<p>Select <b>On</b> to start priming, <b>Off</b> to turn priming off.</p> <p>The <b>Prime</b> function is helpful for filling empty solvent lines or if air has entered the pump heads. The module draws solvent, at high speed with all four pump drives simultaneously, and dispenses it against the waste position of the automatic purge valve. This is done 20 times, before the process comes to an end.</p>
<b>Conditioning</b>	at least 200 bar (> 500 bar is better).	<p>Use this function if you see excessive pressure or composition ripple, and you are sure that the solvent type (aqueous/organic or specific solvent/solvent mix) is correctly set, and there is no evidence of leakage in the pump.</p> <p>Conditioning may be necessary if the pump may contain air, for example after running out of solvent, after a long period of standby or after service or repair.</p>

## Method Parameter Settings

The Binary Pump method setup parameters are in eight sections:

- **Flow**
- **Solvents A and B**
- **Stoptime**
- **Posttime**
- **Pressure Limits**
- **Timetable**
- **Advanced**
- **External Contacts**

**Table 5** Method parameters

Parameter	Limits	Description
<b>Flow</b>	0.00 – 5.00 mL/min in steps of 0.001 . Recommended flow range: 0.05 – 5.00 mL/min .	The flow is the rate of movement of eluent along the column. It is important that the flow rate is kept constant to ensure precise retention time, and peak measurements. Variations in flow rate can occur as a result of the partial failure of the pumping system, air in the pumping system, a change in the mobile phase viscosity or a temperature change.
<b>Solvents A and B</b>		For each channel, you can select which of the two solvents to deliver. You can set the percentage of solvent B to any value from 0 through 100 % . Solvent A always delivers the remaining volume: 100 - %B. The solvent B check boxes allow you to turn the solvent B channels on (checked) or off (cleared).  When the <b>Use solvent types</b> check box in the <b>Compressibility</b> section is checked (see “ <a href="#">Advanced Settings</a> ” on page 52), you click the down arrow and select either a <b>Generic</b> solvent or a calibrated <b>Solvent</b> .
<b>Stoptime</b>	0.01 – 99999 min or <b>As Injector/No Limit</b> (an infinite run time).	The stoptime sets a time limit for your analysis. After the stoptime, all gradients are stopped and the pump parameters return to their initial values. The pump can be used as a stoptime master for the complete analytical system. The pump also stops the detectors if they have a <b>No Limit</b> stoptime setting. If no limit is given, a method will have to be stopped manually.
<b>Posttime</b>	0.01 – 99999 min or <b>Off</b> (0.0 min ).	Your instrument remains in a not ready state during the posttime to delay the start of the next analysis. You can use the <b>Posttime</b> to allow your column to equilibrate after changes in solvent composition (for example after gradient elution).

**Table 5** Method parameters

Parameter	Limits	Description
<b>Pressure Limits</b>	<p><b>Max:</b> 1200 bar (17400 psi ) for flow rates up to 2 mL/min . For flow rates between 2 mL/min and 5 mL/min , the maximum pressure ramps down to 800 bar (11600 psi).</p> <p><b>Min:</b> any value between 0 and the upper pressure limit setting.</p>	<p>Sets the maximum and minimum pressure limits for the pump.</p> <ul style="list-style-type: none"> <li>• <b>Max</b> is the maximum pressure limit at which the pump will switch itself off, protecting the analytical system against over-pressure.</li> <li>• <b>Min</b> is the minimum limit at which the pump will switch itself off, for example, if any solvent reservoir is empty, this prevents system damage by pumping air.</li> </ul>
<b>Timetable</b>		See <a href="#">“Timetable Settings”</a> on page 53
<b>Advanced</b>		See <a href="#">“Advanced Settings”</a> on page 52
<b>External Contacts</b>		<p>The <b>External Contacts</b> section enables you to set up the switching of the external contacts.</p> <div style="background-color: #cccccc; padding: 5px; text-align: center;"><b>NOTE</b></div> <p>The <b>External Contacts</b> section is present only when a BCD/external contacts board is installed.</p>

## Advanced Settings

The Binary Pump advanced method setup parameters are in three sections:

- **Minimum Stroke**
- **Compressibility**
- **Maximum Flow Gradient**

**Table 6** Advanced method parameters

Parameter	Limits	Description
<b>Minimum Stroke</b>	20 – 100 µL	<p>The Stroke Volume is used for optimizing between performance of the module and seal life time. For performance a low stroke volume is beneficial, as it divides disturbances into smaller packages, but a larger volume is extending the life time of the pump seals.</p> <p>If <b>Automatic</b> is activated, the pump tries to achieve an optimized stroke volume for the Jet Weaver geometry.</p> <p><b>Synchronized:</b> Select this option to synchronize the strokes for both channels; the values for Channel B are set to the same as those for Channel A. This is done to avoid floating disturbances affecting instrument performance.</p>
<b>Compressibility</b>		<p>The compressibility of the mobile phase has an effect on the performance of the pump. For best flow accuracy and mixing performance, you can set the parameter according to the mobile phase being used.</p> <p><b>Use solvent types:</b></p> <ul style="list-style-type: none"> <li>• Clear this check box to display the compressibility fields, which allow you to enter compressibility values.</li> <li>• When the check box is selected, the compressibility fields are not displayed, and the enhanced compressibility calibration is enabled. Select the required calibrated solvents from the drop-down lists using the combo boxes in the <b>Solvents</b> section.</li> </ul>
<b>Maximum Flow Gradient</b>	1.000 – 1000.000 mL/min/min in steps of 0.001 mL/min/min Default value: 100.000 mL/min/min	<p>You can set a limit on the rate of change of the solvent flow to protect your analytical column. You can set individual values for <b>Flow ramp up</b> and <b>Flow ramp down</b>.</p>

## Timetable Settings

Use the **Timetable** to program changes in the pump parameters during the analysis by entering a time in the **Time** field and appropriate values in the following fields of the timetable. Changes in flow rate occur linearly from either time zero or the time of the last defined change; other parameters change instantaneously at the time defined in the timetable.

Show **Advanced Timetable** toggles the timetable display between standard mode and advanced mode.

The following parameters can be changed:

- **Change Contacts**
- **Change Flow**
- **Change Max. Pressure Limit**
- **Change Solvent Composition** - You can only use solvents, which have been enabled in the solvents section.
- **Function centric view** - This checkbox allows you displaying parameter changes instead of a time table.

## Priming the Pump

When the solvents have been exchanged or the pumping system has been turned off for a certain time (for example, overnight) oxygen will re-diffuse into the solvent channel between the solvent reservoir, vacuum degassing unit (when available in the system) and the pump. Solvents containing volatile ingredients will slightly lose these. Therefore priming of the pumping system is required before starting an application.

- 1 Initiate a purge in the controlling software with a Purge flow set to 3 – 5 ml/min per channel.
- 2 Flush all tubes with at least 30 ml of solvent.

**Table 7** Choice of Priming Solvents for Different Purposes

Activity	Solvent	Comments
After an installation	Isopropanol	Best solvent to flush air out of the system
When switching between reverse phase and normal phase (both times)	Isopropanol	Isopropanol is miscible with both normal phase and reverse phase solvents.
After an installation	Ethanol or Methanol	Alternative to Isopropanol (second choice) if no Isopropanol is available
To clean the system when using buffers	Bidistilled water	Best solvent to re-dissolve buffer crystals
After a solvent change	Bidistilled water	Best solvent to re-dissolve buffer crystals
Before turning off system for an extended period of time	Organic or 10 % isopropanol in water	

### NOTE

The pump should never be used for priming empty tubings (never let the pump run dry). Use a syringe to draw enough solvent for completely filling the tubings to the pump inlet before continuing to prime with the pump.

**If the system has been run dry or air has diffused into the pump it might require additional steps to get rid of the air again. Following the procedure below will give the best and fastest results.**

- 1** Change solvent to isopropanol on both channels.
- 2** Turn on the Prime function.
- 3** Purge the system with 10 ml, composition 50/50 and for 10 min.
- 4** Attach a column suitable for isopropanol and set the Max. pressure limit to the limit of the column.
- 5** Run the system at composition 50/50 and a flow rate that gives a pressure close to the limit of the column.
- 6** Observe the pressure fluctuations. The system is air free as soon as the pressure is stable.
- 7** Change solvents and column according to the analytical conditions and purge the system to change solvents.

## Solvent Information

### Introduction

Observe the following recommendations on the use of solvents.

- Follow recommendations for avoiding the growth of algae, see [“Algae Growth in HPLC Systems”](#) on page 43.
- Small particles can permanently block capillaries and valves. Therefore, always filter solvents through 0.4 µm filters.
- Avoid or minimize the use of solvents that may corrode parts in the flow path. Consider specifications for the pH range given for different materials like flow cells, valve materials etc. and recommendations in subsequent sections.

### Materials in Flow Path

Following materials are used in the flow path of this module:

Part	Materials
Degasser chamber	TFE/PDD Copolymer, FEP, PEEK, PPS
Microfluidic structures <sup>1</sup>	SST
SSV	PEEK, FEP, PFA, Al <sub>2</sub> O <sub>3</sub> -based ceramic, ruby, sapphire, SST
Passive inlet valve	SST, gold, ruby, ZrO <sub>2</sub> -based ceramic, tantalum
Outlet valve	SST, gold, ruby, ZrO <sub>2</sub> -based ceramic, tantalum
Pump head	SST
Pistons	ZrO <sub>2</sub> -based ceramic
Piston/wash seals	UHMW-PE, SST
Pressure sensor	SST
Automatic purge valve	Polyimide, SST, DLC

<sup>1</sup> Jet Weaver, Heat Exchanger



## Material Information

Materials in the flow path are carefully selected based on Agilent's experiences in developing highest quality instruments for HPLC analysis over several decades. These materials exhibit excellent robustness under typical HPLC conditions. For any special conditions please consult the material information section or contact Agilent.

### Disclaimer

Subsequent data were collected from external resources and are meant as a reference. Agilent cannot guarantee the correctness and completeness of such information. Data is based on compatibility libraries, which are not specific for estimating the long-term life time under specific but highly variable conditions of UHPLC systems, solvents, solvent mixtures and samples. Information can also not be generalized due to catalytic effects of impurities like metal ions, complexing agents, oxygen etc. Apart from pure chemical corrosion, other effects like electro corrosion, electrostatic charging (especially for non-conductive organic solvents), swelling of polymer parts etc. need to be considered. Most data available refers to room temperature (typically 20 – 25 °C, 68 – 77 °F). If corrosion is possible, it usually accelerates at higher temperatures. If in doubt, please consult technical literature on chemical compatibility of materials.

### PEEK

PEEK (Polyether-Ether Ketones) combines excellent properties with regard to chemical resistance, mechanical and thermal stability. It is stable in a pH range of 1 to 12.5 and inert to many common solvents. There are several known incompatibilities with chemicals such as chloroform, methylene chloride, THF, DMSO > 1 %, strong acids (nitric acid > 10 %, sulphuric acid > 10 %, trichloroacetic acid, sulfonic acids), halogenes or aqueous halogene solutions, phenol and derivatives (cresols, salicylic acid etc.). When used above room temperature, PEEK is sensitive to bases and various organic solvents, which can cause it to swell.

## **Polyimide**

Agilent uses semi-crystalline polyimide for rotor seals in valves and needle seats in autosamplers. One supplier of polyimide is DuPont, which brands polyimide as Vespel, which is also used by Agilent.

Polyimide is stable in a pH range between 1 and 10 and in most organic solvents. It is incompatible with concentrated mineral acids (e.g. sulphuric acid), glacial acetic acid, DMSO and THF. It is also degraded by nucleophilic substances like ammonia (e.g. ammonium salts in basic conditions) or acetates.

## **Polyethylene (PE)**

Agilent uses UHMW (ultra-high molecular weight)-PE/PTFE blends for yellow piston and wash seals, which are used in 1290 Infinity pumps and for normal phase applications in 1260 Infinity pumps.

Polyethylene has a good stability for most common inorganic solvents including acids and bases in a pH range of 1 to 12.5 . It is compatible to many organic solvents used in chromatographic systems like methanol, acetonitrile and isopropanol. It has limited stability with aliphatic, aromatic and halogenated hydrocarbons, THF, phenol and derivatives, concentrated acids and bases. For normal phase applications, the maximum pressure should be limited to 200 bar.


## **Tantalum (Ta)**

Tantalum is inert to most common HPLC solvents and almost all acids except fluoric acid and acids with free sulfur trioxide. It can be corroded by strong bases (e.g. hydroxide solutions > 10 %, diethylamine). It is not recommended for the use with fluoric acid and fluorides.

## Stainless Steel (SST)

Stainless steel is inert against many common solvents. It is stable in the presence of acids and bases in a pH range of 1 to 12.5. It can be corroded by acids below pH 2.3. It can also corrode in following solvents:

- Solutions of alkali halides, their respective acids (for example, lithium iodide, potassium chloride, and so on) and aqueous solutions of halogens.
- High concentrations of inorganic acids like nitric acid, sulfuric acid and organic solvents especially at higher temperatures (replace, if your chromatography method allows, by phosphoric acid or phosphate buffer which are less corrosive against stainless steel).
- Halogenated solvents or mixtures which form radicals and/or acids, for example:



This reaction, in which stainless steel probably acts as a catalyst, occurs quickly with dried chloroform if the drying process removes the stabilizing alcohol.

- Chromatographic grade ethers, which can contain peroxides (for example, THF, dioxane, di-isopropylether). Such ethers should be filtered through dry aluminium oxide which adsorbs the peroxides.
- Solutions of organic acids (acetic acid, formic acid, and so on) in organic solvents. For example, a 1 % solution of acetic acid in methanol will attack steel.
- Solutions containing strong complexing agents (for example, EDTA, ethylene diamine tetra-acetic acid).
- Mixtures of carbon tetrachloride with 2-propanol or THF.

## Diamond-Like Carbon (DLC)

Diamond-Like Carbon is inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.

### **Fused silica and Quartz (SiO<sub>2</sub>)**

Fused silica is used in 1290 Infinity Flow Cells and capillaries. Quartz is used for classical flow cell windows. It is inert against all common solvents and acids except hydrofluoric acid and acidic solvents containing fluorides. It is corroded by strong bases and should not be used above pH 12 at room temperature. The corrosion of flow cell windows can negatively affect measurement results. For a pH greater than 12, the use of flow cells with sapphire windows is recommended.

### **Gold**

Gold is inert to all common HPLC solvents, acids and bases within the specified pH range. It can be corroded by complexing cyanides and concentrated acids like aqua regia.

### **Zirconium Oxide (ZrO<sub>2</sub>)**

Zirconium Oxide is inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.

### **Platinum/Iridium**

Platinum/Iridium is inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.

### **Fluorinated polymers (PTFE, PFA, FEP, FFKM)**

Fluorinated polymers like PTFE (polytetrafluorethen), PFA (perfluoroalkoxy) and FEP (fluorinated ethylene propylene) are inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications. FFKM is perfluorinated rubber, which is also resistant to most chemicals. As an elastomer, it may swell in some organic solvents like halogenated hydrocarbons.

### **Sapphire, Ruby and Al<sub>2</sub>O<sub>3</sub>-based ceramics**

Sapphire, ruby and ceramics based on aluminum oxide Al<sub>2</sub>O<sub>3</sub> are inert to almost all common acids, bases and solvents. There are no documented incompatibilities for HPLC applications.



## 5 How to Optimize the Performance of Your Module

Delay Volume and Extra-Column Volume 62

Delay Volume 62

How to Configure the Optimum Delay Volume 63

How to Achieve Higher Resolution 65

Using Solvent Calibration Tables 68

This chapter gives hints on how to optimize the performance or use additional devices.



## Delay Volume and Extra-Column Volume

The *delay volume* is defined as the system volume between the point of mixing in the pump and the top of the column.

The *extra-column volume* is defined as the volume between the injection point and the detection point, excluding the volume in the column.

### Delay Volume

In gradient separations, this volume causes a delay between the mixture changing in the pump and that change reaching the column. The delay depends on the flow rate and the delay volume of the system. In effect, this means that in every HPLC system there is an additional isocratic segment in the gradient profile at the start of every run. Usually the gradient profile is reported in terms of the mixture settings at the pump and the delay volume is not quoted even though this will have an effect on the chromatography. This effect becomes more significant at low flow rates and small column volumes and can have a large impact on the transferability of gradient methods. It is important, therefore, for fast gradient separations to have small delay volumes, especially with narrow bore columns (e.g., 2.1 mm i.d.) as often used with mass spectrometric detection.

The delay volume in a system includes the volume in the pump from the point of mixing, connections between pump and autosampler, volume of the flow path through the autosampler and connections between autosampler and column.

## How to Configure the Optimum Delay Volume

The physical delay volume of the pump depends primarily on the use of the Jet Weaver mixer. For UV detection the Jet Weaver should always be used but for mass spectrometric detection the user can decide to bypass the Jet Weaver in order to reduce the delay volume. This only makes sense for ultra-fast gradient operation (less than 0.5 min) or for use with very small volume columns. If the Jet Weaver is bypassed the connection tubing to the autosampler is routed directly from the purge valve.

### NOTE

Before disconnecting a Jet Weaver from the flow path, flush it with organic solvent. Avoid leaving water or buffers inside the Jet Weaver, which may cause the growth of microorganisms like algae or bacteria.

Sometimes it may be advisable to increase the delay volume in the pump. Specifically this can be the case when UV detection is employed and a strongly UV-absorbing compound has been added to the mobile phase. This can have the effect of emphasizing any pump noise and the most common example is the use of trifluoroacetic acid (TFA) in the analysis of proteins and peptides. The effect can be mitigated by increasing the mixer volume.

The following different Jet Weaver configurations are available:

- The Jet Weaver 35 µL/100 µL(G4220-60006) has two alternative volumes in the same unit.

The switch from the lower volume, 35 µl, to the higher volume, 100 µl, is done by uninstalling it, turning it around from front to back and re-installing it, see [“Changing configuration or replacing the Jet Weaver”](#) on page 127. The mixing volume (and hence delay volume) is increased by 65 µl and the baseline performance with additives like TFA will be improved. The configuration of the Jet Weaver is logged automatically by an attached RFID tag.

- The 380 µL Jet Weaver high performance mixer is optionally available for demanding applications, which use solvents in different channels (for example A versus B), that differ strongly in their UV/Vis absorption, for example by using trifluoroacetic acid (TFA) as a modifier, which has a high absorbance.

## 5 How to Optimize the Performance of Your Module

### How to Configure the Optimum Delay Volume

Solvent packages created by the pump may persist until the solvent reaches the detector flow cell. Absorption fluctuations can then show up as baseline noise, also referred to as mixing noise. Applications like impurity quantitation or lowest level compound detection require minimizing this noise. The 380  $\mu\text{L}$  Jet Weaver strongly improves mixing and therefore reduces baseline noise and improves sensitivity in detection. Patented Agilent microfluidic technology offers high mixing performance at a low internal volume of 380  $\mu\text{L}$ .



## How to Achieve Higher Resolution

Increased resolution in a separation will improve the qualitative and quantitative data analysis, allow more peaks to be separated or offer further scope for speeding up the separation. This section considers how resolution can be increased by examining the following points:

- Optimize selectivity
- Smaller particle-size packing
- Longer Columns
- Shallower gradients, faster flow

Resolution between two peaks is described by the resolution equation:

$$R_s = \frac{1}{4} \sqrt{N} \frac{(\alpha - 1)}{\alpha} \frac{(k_2 + 1)}{k_2}$$

where

- $R_s$ =resolution,
- $N$ =plate count (measure of column efficiency),
- $\alpha$ =selectivity (between two peaks),
- $k_2$ =retention factor of second peak (formerly called capacity factor).

The term that has the most significant effect on resolution is the selectivity,  $\alpha$ , and practically varying this term involves changing the type of stationary phase (C18, C8, phenyl, nitrile etc.), the mobile phase and temperature to maximize the selectivity differences between the solutes to be separated. This is a substantial piece of work which is best done with an automated method development system which allows a wide range of conditions on different columns and mobile phases to be assessed in an ordered scouting protocol. This section considers how to get higher resolution with any chosen stationary and mobile phases. If an automated method development system was used in the decision on phases it is likely that short columns were used for fast analysis in each step of the scouting.

The resolution equation shows that the next most significant term is the plate count or efficiency,  $N$ , and this can be optimized in a number of ways.  $N$  is inversely proportional to the particle size and directly proportional to the

length of a column and so smaller particle size and a longer column will give a higher plate number. The pressure rises with the inverse square of the particle size and proportionally with the length of the column. This is the reason that the 1290 Infinity LC system was designed to go to higher pressures so that it can run sub-two-micron particles and column length can be increased to 100 mm or 150 mm. There are even examples of 100 mm and 150 mm columns linked to give 250 mm length. Resolution increases with the square root of N so doubling the length of the column will increase resolution by a factor of 1.4. What is achievable depends on the viscosity of the mobile phase as this relates directly to the pressure. Methanol mixtures will generate more back pressure than acetonitrile mixtures. Acetonitrile is often preferred because peak shapes are better and narrower in addition to the lower viscosity but methanol generally yields better selectivity (certainly for small molecules less than about 500 Da). The viscosity can be reduced by increasing the temperature but it should be remembered that this can change the selectivity of the separation. Experiment will show if this leads to increase or decrease in selectivity. As flow and pressure are increased it should be remembered that frictional heating inside the column will increase and that can lead to slightly increased dispersion and possibly a small selectivity change both of which could be seen as a reduction in resolution. The latter case might be offset by reducing the temperature of the thermostat by a few degrees and again experiment will reveal the answer.

The van Deemter curve shows that the optimum flow rate through an STM column is higher than for larger particles and is fairly flat as the flow rate increases. Typical, close to optimum, flow rates for STM columns are: 2 ml/min for 4.6 mm i.d.; and 0.4 ml/min for 2.1 mm i.d. columns.

In isocratic separations, increasing the retention factor, k, results in better resolution because the solute is retained longer. In gradient separations the retention is described by  $k^*$  in the following equation:

$$k^* = \frac{t_G}{\Delta\%B} \cdot \frac{F}{V_m} \cdot \frac{100}{S}$$

where:

- $k^*$  = mean k value,
- $t_G$  = time length of gradient (or segment of gradient) (min),
- F = flow (ml/min),
- $V_m$  = column delay volume,

- $\Delta\%B$  = change in fraction of solvent B during the gradient,
- $S$  = constant (ca. 4-5 for small molecules).

This shows that  $k$  and hence resolution can be increased by having a shallower gradient (2 to 5 %/min change is a guideline), higher flow rate and a smaller volume column. This equation also shows how to speed up an existing gradient – if the flow is doubled but the gradient time is halved,  $k^*$  remains constant and the separation looks the same but happens in half the time. Recently published research has shown how a shorter STM column (at temperatures above 40 °C) can generate higher peak capacity than a longer STM column by virtue of running it faster. (Refer to *Petersson et al., J.Sep.Sci, 31, 2346-2357, 2008, Maximizing peak capacity and separation speed in liquid chromatography*).

## Using Solvent Calibration Tables

### Importing Solvent Calibration Tables

RC.NET based Agilent graphical user interfaces (ChemStation, EZChrom Elite, OpenLab etc.) include data for most commonly used solvents in HPLC. This data contains solvent properties and is used for optimum pump control in order to ensure best flow and composition accuracy.

If your solvent is not included to the software, please check the Agilent web site

<http://www.chem.agilent.com/en-US/Support/Downloads/firmware/Pages/LC.aspx> for additional libraries (registration required), which also provides updates and optimized data.

If your solvent is neither available in the user interface nor in the library, please use generic solvents. "Generic aqueous" gives good results for most solvent mixtures with at least 50 % water, which have similar properties as pure water. For other solvents with high organic percentage, "Generic organic" gives a good approximation.

### Importing Solvent Calibration in ChemStation

- 1 Go to menu **Instrument > Instrument configuration**.
- 2 In the **Instrument Configuration** screen choose your module and click **Configure**.
- 3 Click **Configure Solvent Type Catalogs**.
- 4 In **Solvent Type Catalogs** click **Import**.
- 5 Navigate to the location of the solvent calibration table and click **Open**.
- 6 The new solvent will now appear in the **Solvent Type Catalogs**.
- 7 The imported solvent is now available for selection as a solvent type, see [Table 5](#) on page 50.



## 6 Troubleshooting and Diagnostics

Overview of the Module's Indicators and Test Functions 70

Status indicators 71

Power Supply Indicator 71

Module Status Indicator 72

Available Tests vs User Interfaces 73

Agilent Lab Advisor Software 74

Overview about the troubleshooting and diagnostic features.



## Overview of the Module's Indicators and Test Functions

### Status Indicators

The module is provided with two status indicators which indicate the operational state (prerun, run, and error states) of the module. The status indicators provide a quick visual check of the operation of the module.

### Error Messages

In the event of an electronic, mechanical or hydraulic failure, the module generates an error message in the user interface. For each message, a short description of the failure, a list of probable causes of the problem, and a list of suggested actions to fix the problem are provided (see chapter Error Information).

### Test Functions


A series of test functions are available for troubleshooting and operational verification after exchanging internal components (see Tests and Calibrations).

### Diagnostic Signals

The module has several signals (internal temperatures, voltages and currents of lamps) that can be used for diagnosing baseline problems. These can be added like normal signals in the Agilent ChemStation software.

## Status indicators

Two status indicators are located on the front of the module. The lower left indicates the power supply status, the upper right indicates the module status.



**Figure 10** Location of Status Indicators

### Power Supply Indicator

The power supply indicator is integrated into the main power switch. When the indicator is illuminated (*green*) the power is *ON*.

## Module Status Indicator

The module status indicator indicates one of six possible module conditions:

- When the status indicator is *OFF* (and power switch light is on), the module is in a *prerun* condition, and is ready to begin an analysis.
- A *green* status indicator, indicates the module is performing an analysis (*run* mode).
- A *yellow* indicator indicates a *not-ready* condition. The module is in a not-ready state when it is waiting for a specific condition to be reached or completed (for example, immediately after changing a set point), or while a self-test procedure is running.
- An *error* condition is indicated when the status indicator is *red*. An error condition indicates the module has detected an internal problem which affects correct operation of the module. Usually, an error condition requires attention (e.g. leak, defective internal components). An error condition always interrupts the analysis.

If the error occurs during analysis, it is propagated within the LC system, i.e. a red LED may indicate a problem of a different module. Use the status display of your user interface for finding the root cause/module of the error.

- A *blinking* indicator indicates that the module is in resident mode (e.g. during update of main firmware).
- A *fast blinking* indicator indicates that the module is in a low-level error mode. In such a case try to re-boot the module or try a cold-start (see [“Special Settings”](#) on page 218). Then try a firmware update (see [“Replacing Module Firmware”](#) on page 162). If this does not help, a main board replacement is required.



## Available Tests vs User Interfaces

- Depending on the user interface, the available tests and the screens/reports may vary (see chapter *"Test Functions and Calibrations"*).
- Preferred tool should be the Agilent Lab Advisor software, see ["Agilent Lab Advisor Software"](#) on page 74.
- The Agilent ChemStation B.04.02 and above may not include any maintenance/test functions.
- Screenshots used within these procedures are based on the Agilent Lab Advisor software.

## Agilent Lab Advisor Software

The Agilent Lab Advisor software is a standalone product that can be used with or without data system. Agilent Lab Advisor software helps to manage the lab for high quality chromatographic results and can monitor in real time a single Agilent LC or all the Agilent GCs and LCs configured on the lab intranet.

Agilent Lab Advisor software provides diagnostic capabilities for all Agilent 1200 Infinity Series modules. This includes diagnostic capabilities, calibration procedures and maintenance routines for all the maintenance routines.

The Agilent Lab Advisor software also allows users to monitor the status of their LC instruments. The Early Maintenance Feedback (EMF) feature helps to carry out preventive maintenance. In addition, users can generate a status report for each individual LC instrument. The tests and diagnostic features as provided by the Agilent Lab Advisor software may differ from the descriptions in this manual. For details refer to the Agilent Lab Advisor software help files.

The Instrument Utilities is a basic version of the Lab Advisor with limited functionality required for installation, use and maintenance. No advanced repair, troubleshooting and monitoring functionality is included.



## 7 Error Information

What Are Error Messages	77
General Error Messages	78
Timeout	78
Shutdown	79
Remote Timeout	79
Lost CAN Partner	80
Leak Sensor Short	80
Leak Sensor Open	81
Compensation Sensor Open	81
Compensation Sensor Short	82
Fan Failed	82
Leak	83
Open Cover	83
Cover Violation	84
Pump Error Messages	85
Pressure of binary pump above upper limit	85
Pressure below lower limit	86
Target pressure not reached for binary pump degasser	86
Degasser's pressure limit violation	87
Solvent counter exceeded limit	87
Waste counter limit exceeded	88
Flow rate limit exceeded	88
Binary pump shutdown during analysis	89
Reading the pump encoder tag failed	89
Writing the pump encoder tag failed	89
Pump drive blocked or encoder failed	90
Drive current too low	90
Drive Encoder failed	90



Drive current too high	91
Drive timeout	91
Overcurrent of pump drive	91
Overcurrent of solvent selection valve (SSV)	92
Deliver underrun	92
Defect connection between main board and pump drive encoder	92
Pump drive encoder defect	93
Purge valve failed	93
Reading of purge valve tag failed	93
Pump drive encoder rollover	94
Drive position limit	94
Insufficient power of drive encoder LED	94
Drive encoder error	95
Writing the purge valve tag failed	95
Current of primary pump drive too high	96
Current of secondary pump drive too high	96
Unknown purge valve type	97
Pump drive encoder error	97
Pump drive error	97
Pump drive stroke blocked	98
Pump drive stop not found	98
Pressure sensor calibration wrong or missing	98
Seal wash pump was missing when tried to turn on	99

This chapter describes the meaning of error messages, and provides information on probable causes and suggested actions how to recover from error conditions.

## What Are Error Messages

Error messages are displayed in the user interface when an electronic, mechanical, or hydraulic (flow path) failure occurs which requires attention before the analysis can be continued (for example, repair, or exchange of consumables is necessary). In the event of such a failure, the red status indicator at the front of the module is switched on, and an entry is written into the module logbook.

## General Error Messages

General error messages are generic to all Agilent series HPLC modules and may show up on other modules as well.

### Timeout

**Error ID: 0062**

The timeout threshold was exceeded.

**Probable cause**

- 1** The analysis was completed successfully, and the timeout function switched off the module as requested.
- 2** A not-ready condition was present during a sequence or multiple-injection run for a period longer than the timeout threshold.

**Suggested actions**

Check the logbook for the occurrence and source of a not-ready condition. Restart the analysis where required.

Check the logbook for the occurrence and source of a not-ready condition. Restart the analysis where required.

## Shutdown

### Error ID: 0063

An external instrument has generated a shutdown signal on the remote line.

The module continually monitors the remote input connectors for status signals. A LOW signal input on pin 4 of the remote connector generates the error message.

Probable cause	Suggested actions
1 Leak detected in another module with a CAN connection to the system.	Fix the leak in the external instrument before restarting the module.
2 Leak detected in an external instrument with a remote connection to the system.	Fix the leak in the external instrument before restarting the module.
3 Shut-down in an external instrument with a remote connection to the system.	Check external instruments for a shut-down condition.

## Remote Timeout

### Error ID: 0070

A not-ready condition is still present on the remote input. When an analysis is started, the system expects all not-ready conditions (for example, a not-ready condition during detector balance) to switch to run conditions within one minute of starting the analysis. If a not-ready condition is still present on the remote line after one minute the error message is generated.

Probable cause	Suggested actions
1 Not-ready condition in one of the instruments connected to the remote line.	Ensure the instrument showing the not-ready condition is installed correctly, and is set up correctly for analysis.
2 Defective remote cable.	Exchange the remote cable.
3 Defective components in the instrument showing the not-ready condition.	Check the instrument for defects (refer to the instrument's documentation).

## Lost CAN Partner

### Error ID: 0071

During an analysis, the internal synchronization or communication between one or more of the modules in the system has failed.

The system processors continually monitor the system configuration. If one or more of the modules is no longer recognized as being connected to the system, the error message is generated.

#### Probable cause

- 1**   CAN cable disconnected.
- 2**   Defective CAN cable.
- 3**   Defective main board in another module.

#### Suggested actions

- Ensure all the CAN cables are connected correctly.
  - Ensure all CAN cables are installed correctly.
- Exchange the CAN cable.
- Switch off the system. Restart the system, and determine which module or modules are not recognized by the system.

## Leak Sensor Short

### Error ID: 0082

The leak sensor in the module has failed (short circuit).

The current through the leak sensor is dependent on temperature. A leak is detected when solvent cools the leak sensor, causing the leak sensor current to change within defined limits. If the current increases above the upper limit, the error message is generated.

#### Probable cause

- 1**   Defective leak sensor.
- 2**   Leak sensor incorrectly routed, being pinched by a metal component.

#### Suggested actions

- Please contact your Agilent service representative.
- Please contact your Agilent service representative.



## Leak Sensor Open

### Error ID: 0083

The leak sensor in the module has failed (open circuit).

The current through the leak sensor is dependent on temperature. A leak is detected when solvent cools the leak sensor, causing the leak-sensor current to change within defined limits. If the current falls outside the lower limit, the error message is generated.

#### Probable cause

- 1 Leak sensor not connected to the main board.
- 2 Defective leak sensor.
- 3 Leak sensor incorrectly routed, being pinched by a metal component.

#### Suggested actions

- Please contact your Agilent service representative.
- Please contact your Agilent service representative.
- Please contact your Agilent service representative.

## Compensation Sensor Open

### Error ID: 0081

The ambient-compensation sensor (NTC) on the main board in the module has failed (open circuit).

The resistance across the temperature compensation sensor (NTC) on the main board is dependent on ambient temperature. The change in resistance is used by the leak circuit to compensate for ambient temperature changes. If the resistance across the sensor increases above the upper limit, the error message is generated.

#### Probable cause

- 1 Defective main board.

#### Suggested actions

- Please contact your Agilent service representative.

## Compensation Sensor Short

### Error ID: 0080

The ambient-compensation sensor (NTC) on the main board in the module has failed (short circuit).

The resistance across the temperature compensation sensor (NTC) on the main board is dependent on ambient temperature. The change in resistance is used by the leak circuit to compensate for ambient temperature changes. If the resistance across the sensor falls below the lower limit, the error message is generated.

#### Probable cause

- 1 Defective main board.

#### Suggested actions

Please contact your Agilent service representative.

## Fan Failed

### Error ID: 0068

The cooling fan in the module has failed.

The hall sensor on the fan shaft is used by the main board to monitor the fan speed. If the fan speed falls below a certain limit for a certain length of time, the error message is generated.

Depending on the module, assemblies (e.g. the lamp in the detector) are turned off to assure that the module does not overheat inside.

#### Probable cause

- 1 Fan cable disconnected.
- 2 Defective fan.
- 3 Defective main board.

#### Suggested actions

Please contact your Agilent service representative.

Please contact your Agilent service representative.

Please contact your Agilent service representative.

## Leak

### Error ID: 0064

A leak was detected in the module.

The signals from the two temperature sensors (leak sensor and board-mounted temperature-compensation sensor) are used by the leak algorithm to determine whether a leak is present. When a leak occurs, the leak sensor is cooled by the solvent. This changes the resistance of the leak sensor which is sensed by the leak-sensor circuit on the main board.

#### Probable cause

- 1 Loose fittings.
- 2 Broken capillary.

#### Suggested actions

- Ensure all fittings are tight.
- Exchange defective capillaries.

## Open Cover

### Error ID: 0205

The top foam has been removed.

#### Probable cause

- 1 Foam not activating the sensor.
- 2 Dirty or defective sensor.

#### Suggested actions

- Please contact your Agilent service representative.
- Please contact your Agilent service representative.

## Cover Violation

### Error ID: 7461

The top foam has been removed.

The sensor on the main board detects when the top foam is in place. If the foam is removed while the lamps are on (or if an attempt is made to switch on for example the lamps with the foam removed), the lamps are switched off, and the error message is generated.

#### Probable cause

- 1 The top foam was removed during operation.
- 2 Foam not activating the sensor.

#### Suggested actions

Please contact your Agilent service representative.

Please contact your Agilent service representative.

## Pump Error Messages

These errors are pump specific.

### Pressure of binary pump above upper limit

**Error ID: 22014**

The pressure has exceeded the upper pressure limit.

- Parameter: Measured pressure

**Probable cause**

- 1** Blockage in flow path after the pressure sensor.
- 2** Inappropriate settings (pressure limit, flow rate).

**Suggested actions**

- Check for blockages in the LC system, e.g. purge valve, Jet Weaver, degraded column, column frits, needle, needle seat, capillaries etc.
- Check for particles in the solvent.
- Decrease flow rate.
- Increase pressure limit.

## Pressure below lower limit

### Error ID: 22015

The pressure has dropped below the lower limit.

- Parameter: None

#### Probable cause

- 1 Leak
- 2 Bottle empty
- 3 Wrong solvent (viscosity)
- 4 Inappropriate setting
- 5 Column degradation

#### Suggested actions

- Check for leaks.
- Check bottle filling.
- Check solvent.
- Check flow rate and lower pressure limit.
- Replace column.

## Target pressure not reached for binary pump degasser

### Error ID: 22031

The target pressure of the binary pump degasser has not been reached within the expected time.

- Parameter: Pressure in mbar

#### Probable cause

- 1 Condensation in degasser chamber due to temperature fluctuation.
- 2 Degasser is defect.

#### Suggested actions

- Equilibrate and restart module.
- Please contact your Agilent service representative.

## Degasser's pressure limit violation

### Error ID: 22032

Pressure too far above the limit.

#### Probable cause

- 1 Leak in degasser chamber or degasser tubing.
- 2 Defect vacuum pump.
- 3 Degasser chamber empty or connected to air.

#### Suggested actions

Please contact your Agilent service representative.

Please contact your Agilent service representative.

Block unused degasser channels.

## Solvent counter exceeded limit

### Error ID: 22055

The counter for the solvent volume has exceeded the limit, which has been set in the user interface.

Parameter:

- Without Solvent Selection Valve:  
 0 for channel A, 1 for channel B
- With Solvent Selection Valve:  
 2 for channel A1, 3 for channel B1, 4 for channel A2, 5 for channel B2

#### Probable cause

- 1 No solvent present.
- 2 Inappropriate setting.

#### Suggested actions

Refill solvent bottle.

Check solvent counter setting in user interface.

## Waste counter limit exceeded

### Error ID: 22056

The counter for the waste volume has exceeded the limit, which has been set in the user interface.

- Parameter: None

#### Probable cause

- 1 The waste container is full.
- 2 Inappropriate setting for waste counter.

#### Suggested actions

- Empty waste container.
- Reset waste counter.
  - Adjust waste counter limit.

## Flow rate limit exceeded

### Error ID: 22064

The flow rate of the binary pump has exceeded the limit, while the pump runs in pressure controlled mode, e.g. during a pressure test.

- Parameter: None

#### Probable cause

- 1 Leak
- 2 Bottle empty.
- 3 Shutoff valve closed.
- 4 Drift of pressure sensor (unlikely for short tests taking some minutes).

#### Suggested actions

- Check for leaks in the pump and flow path.
- Fill solvent bottle.
- Open shutoff valve.
- Replace pressure sensor.



## Binary pump shutdown during analysis

### Error ID: 22065

The binary pump has been shut down by the control software or control module during an analysis.

- Parameter: 0 for off, 1 for standby.

#### Probable cause

- 1 Pump has been shut down.

#### Suggested actions

Restart pump.

## Reading the pump encoder tag failed

### Error ID: 22402

Reading the pump encoder tag has failed.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Defect connection between encoder and main board.
- 2 Missing or defect tag Defect connection between tag and encoder.

#### Suggested actions

Please contact your Agilent service representative.

Please contact your Agilent service representative.

## Writing the pump encoder tag failed

### Error ID: 22405

Writing the pump encoder tag has failed.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Defect connection between encoder and main board.
- 2 Defect tag Defect connection between tag and encoder.

#### Suggested actions

Please contact your Agilent service representative.

Please contact your Agilent service representative.

## Pump drive blocked or encoder failed

### Error ID: 22406

Pump drive blocked or encoder failed.

- Parameter: None

#### Probable cause

- 1 Blockage of the pump drive Drive encoder failed.

#### Suggested actions

Please contact your Agilent service representative.

## Drive current too low

### Error ID: 22407

The current consumption of the pump drive is too low.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Drive motor defect.
- 2 Wrong/missing connection of pump drive to main board.

#### Suggested actions

Please contact your Agilent service representative.

Please contact your Agilent service representative.

## Drive Encoder failed

### Error ID: 22408

Drive encoder failed during pump drive calibration.

#### Probable cause

- 1 Internal error.

#### Suggested actions

Contact Agilent support.

## Drive current too high

### Error ID: 22409

The current consumption of the pump drive is too high.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Blockage of system before pressure sensor.
- 2 Drive motor defect.

#### Suggested actions

Check for blockage of e.g. outlet valve filter frit, purge valve, heat exchanger.

Please contact your Agilent service representative.

## Drive timeout

### Error ID: 22410

Drive is blocked mechanically, fails during initialization.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Blockage of pump drive Drive motor defect.

#### Suggested actions

Please contact your Agilent service representative.

## Overcurrent of pump drive

### Error ID: 22411

The current consumption of the pump drive is too high.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Blockage of system before pressure sensor.
- 2 Drive motor defect.

#### Suggested actions

Check for blockage of e.g. outlet valve filter frit, purge valve, heat exchanger.

Please contact your Agilent service representative.

## Overcurrent of solvent selection valve (SSV)

**Error ID: 22412**

Overcurrent of solvent selection valve (SSV).

- Parameter: None

**Probable cause**

- 1 Valve defect.

**Suggested actions**

Replace the solvent selection valve.

## Deliver underrun

**Error ID: 22413**

Internal error.

- Parameter: None

**Probable cause**

- 1 Internal error.

**Suggested actions**

Please contact your Agilent service representative.

## Defect connection between main board and pump drive encoder

**Error ID: 22414**

Defect connection between main board and pump drive encoder.

- Parameter: 1 – 4 referring to pump drive

**Probable cause**

- 1 Defect connection between main board and pump drive encoder.
- 2 Defect encoder.

**Suggested actions**

Please contact your Agilent service representative.

Please contact your Agilent service representative.

## Pump drive encoder defect

### Error ID: 22415

Defect pump drive encoder.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Defect encoder.

#### Suggested actions

Please contact your Agilent service representative.

## Purge valve failed

### Error ID: 22417

Lost steps of the purge valve encoder.

- Parameter: None

#### Probable cause

- 1 Purge valve drive mechanically blocked or defect.

#### Suggested actions

- Check installation of purge valve head.
- Please contact your Agilent service representative.

## Reading of purge valve tag failed

### Error ID: 22420

Reading the purge valve tag failed.

- Parameter: None

#### Probable cause

- 1 Reading of purge valve tag failed.
- 2 Purge valve head tag defect or empty.
- 3 Purge valve tag reader is defect.

#### Suggested actions

Check cable connection.  
Replace purge valve head.  
Please contact your Agilent service representative.

## Pump drive encoder rollover

**Error ID: 22424**

Invalid pump drive encoder signals have been detected.

- Parameter: 1 – 4 referring to pump drive

**Probable cause**

- 1 Pump drive encoder is defect.

**Suggested actions**

Please contact your Agilent service representative.

## Drive position limit

**Error ID: 22425**

Internal error.

- Parameter: 1 – 4 referring to pump drive

**Probable cause**

- 1 Internal error.

**Suggested actions**

Please contact your Agilent service representative.

## Insufficient power of drive encoder LED

**Error ID: 22426**

Insufficient power of drive encoder LED.

- Parameter: 1 – 4 referring to pump drive

**Probable cause**

- 1 Pump drive encoder is defect.

**Suggested actions**

Please contact your Agilent service representative.

## Drive encoder error

### Error ID: 22427- 22430

An error has occurred for the pump drive encoder.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Pump drive encoder is defect.

#### Suggested actions

Please contact your Agilent service representative.

## Writing the purge valve tag failed

### Error ID: 22431

Writing the purge valve tag failed.

- Parameter: None

#### Probable cause

- 1 Purge valve head tag defect.
- 2 Purge valve tag reader is defect.

#### Suggested actions

Replace purge valve head.

Please contact your Agilent service representative.

## Current of primary pump drive too high

### Error ID: 22433

The current of the primary pump drive is too high.

- Parameter: 1 or 4 referring to pump drive.

#### Probable cause

- 1 Blockage of flow path between primary pump head and pressure sensor, e.g. of the heat exchanger.

- 2 Primary pump drive is defect.

#### Suggested actions

- Check for blockages in flow path.
- Please contact your Agilent service representative.

Please contact your Agilent service representative.

## Current of secondary pump drive too high

### Error ID: 22434

The current of the secondary pump drive is too high.

- Parameter: 2 or 3 referring to pump drive

#### Probable cause

- 1 Blockage of flow path between secondary pump head and pressure sensor, e.g. of the heat exchanger.

- 2 Secondary pump drive is defect.

#### Suggested actions

- Check for blockages in the flow path.
- Please contact your Agilent service representative.

Please contact your Agilent service representative.



## Unknown purge valve type

### Error ID: 22435

The type information of the purge valve is invalid.

- Parameter: None

#### Probable cause

- 1 Wrong valve head installed.
- 2 Valve head has invalid RFID tag content.

#### Suggested actions

- Check or replace purge valve head.
- Check or replace purge valve head.

## Pump drive encoder error

### Error ID: 22437

The pump drive encoder has generated no signal.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Pump drive encoder is defect.

#### Suggested actions

Please contact your Agilent service representative.

## Pump drive error

### Error ID: 22438, 22439

The pump drive failed during calibration.

- Parameter: 1 – 4 referring to pump drive

#### Probable cause

- 1 Pump drive motor defect or mechanically blocked.

#### Suggested actions

Please contact your Agilent service representative.

## Pump drive stroke blocked

**Error ID: 22441**

The pump drive movement is blocked.

- Parameter: 1 – 4 referring to pump drive

**Probable cause**

- 1 Pump drive motor is mechanically blocked.

**Suggested actions**

Please contact your Agilent service representative.

## Pump drive stop not found

**Error ID: 22442**

The pump drive stop has not been found.

- Parameter: 1 – 4 referring to pump drive

**Probable cause**

- 1 Pump drive spindle is defect.

**Suggested actions**

Please contact your Agilent service representative.

## Pressure sensor calibration wrong or missing

**Error ID: 22443**

Pressure sensor calibration wrong or missing.

- Parameter: None

**Probable cause**

- 1 Pressure sensor calibration wrong or missing.

**Suggested actions**

- Replace pressure sensor.
- Please contact your Agilent service representative.

## Seal wash pump was missing when tried to turn on

### Error ID: 22499

The seal wash pump has not been detected (while being configured or detected before)

#### Probable cause

- 1 Defect cable connection to seal wash pump.
- 2 Defect seal wash pump motor.
- 3 Defective main board.

#### Suggested actions

- Check cable connection.
- Please contact your Agilent service representative.
- Please contact your Agilent service representative.

## **7 Error Information**

### **Pump Error Messages**



## 8 Test Functions and Calibrations

Introduction	102
System Pressure Test	103
System Pressure Test Evaluation	104
Pump Head Leak Test	105
Pump Head Leak Test Evaluation	107

This chapter describes the tests for the module.



## Introduction

All tests are described based on the Agilent Lab Advisor Software B.01.04. Other user interfaces may not provide any test or just a few.

**Table 8** Interfaces and available test functions

Interface	Comment	Available Function
Agilent Instrument Utilities	Maintenance tests available	<ul style="list-style-type: none"> <li>• System pressure test</li> </ul>
Agilent Lab Advisor	All tests are available	<ul style="list-style-type: none"> <li>• System pressure test</li> <li>• Pump head leak test</li> </ul>
Agilent ChemStation	No tests available Adding of pressure to chromatographic signals possible	<ul style="list-style-type: none"> <li>• Pressure</li> <li>• Pressure ripple</li> <li>• Temperature main board</li> </ul>
Agilent Instant Pilot	Some tests are available	<ul style="list-style-type: none"> <li>• System pressure test</li> <li>• Monitoring of values                             <ul style="list-style-type: none"> <li>• Pressure</li> <li>• Pressure ripple</li> <li>• Flow (in case of operating pressure)</li> </ul> </li> </ul>

For details on the use of the interface refer to the interface documentation.

# System Pressure Test

The test determines the leak rate of the system between pump outlet valves and a blank nut. The blank nut can be positioned at different locations in the system before the flow cell, to determine and verify the leak rate of individual modules and components. The test allows for setting the pressure at which the test is performed. The leak rate of high pressure parts is not always a linear function and therefore it is recommended to perform the test at a pressure that corresponds to the normal operating pressure of the system.

**When** In case of a suspected leak. To verify successful execution of maintenance tasks.

Parts required	#	p/n	Description
	1	01080-83202	Blank nut

**Preparations** Solvents must be present in both channels.

- 1 Run the **System Pressure Test** with the recommended user interface (for further information see Online-Help of user interface).




Figure 11 System Pressure Test – Result




Figure 12 System Pressure Test – dynamic pressure input

## System Pressure Test Evaluation

### Test Failed

#### Probable cause

- 1 Damaged blank nut (poorly shaped from over tightening)
- 2 Pump leakages
- 3 Loose or leaky fittings
- 4 Autosampler leakages
- 5 Thermostatted Column Compartment valve leakages

#### Suggested actions

- Before investigating any other possible sources of failure make sure that the blank nut you are using is in a good condition and properly tightened.
- Perform the Pump Head Leak test.
- Tighten the fittings or replace capillaries.
- Perform the Autosampler Leak test.
- Replace the TCC valve rotor seal.

### NOTE

Notice the difference between *error* in the test and a *failed* result! An *error* is caused by an abnormal termination during the operation of the test, whereas a *failed* result indicates that the test result were not within the specified limits.



## Pump Head Leak Test

The test determines the leakage of the individual pump heads, by blocking each pump head separately and pressurizing to a specified level. The flow is delivered by the secondary piston and thereby the parts between the outlet ball valve and the purge valve are tested for leak tightness. The test allows for setting the pressure at which the test is performed. The leak rate of high pressure parts is not always a linear function and therefore it is recommended to perform the test at a pressure that corresponds to the normal operating pressure of the system.

<b>When</b>	Excessive pressure ripple or suspected pump performance problems.
<b>Preparations</b>	Solvents must be present in both channels.

# 8 Test Functions and Calibrations

## Pump Head Leak Test

- 1 Run the **Pump Head Leak Test** with the recommended user interface (for further information see Online-Help of user interface).

Name	Value
Maximum system pressure	1200 bar
Channel A2	Done
Channel B2	Done
Remaining purge time	0 min 0 sec

Figure 13 Pump Head Leak Test – Results

Enter the test pressure

Enter the pressure at which the test will be executed 1200 bar

OK Cancel

Figure 14 Pump Head Leak Test - dynamic pressure input

# Pump Head Leak Test Evaluation

## Test Failed

### Probable cause

- 1 Loose or leaky fittings
- 2 Leaky High Pressure Filter Assembly
- 3 Damaged pump seals or pistons
- 4 Leaky outlet ball valve
- 5 Leaky purge valve

### Suggested actions

- Tighten the fittings or replace capillaries.
- Tighten Filter housing and fitting.
  - NOTE - Fittings should only be tightened in a system that is not pressurized, for better chance of success.
  - Re-run test.

- Exchange the pump head or pump seals.
  - NOTE - A wrongly installed pump seal will reduce the life time of the pump seal considerably.


Exchange pump head or outlet ball valve.

Exchange Purge valve rotor seal.
- NOTE
- Notice the difference between *error* in the test and a *failed* result! An *error* is caused by an abnormal termination during the operation of the test, whereas a *failed* result indicates that the test result were not within the specified limits.
- Agilent 1290 Infinity Binary Pump User Manual

107

## **8 Test Functions and Calibrations**

### **Pump Head Leak Test**



## 9 Maintenance


Introduction to Maintenance	110
Warnings and Cautions	112
Overview of Maintenance	113
Cleaning the Module	114
Installing Fittings and Capillaries	115
Replacing the Shutoff Valve Panel	116
Replacing the Pressure Sensor	118
Replacing the Inlet valve	120
Replacing the Outlet valve	122
Replacing the Solvent Selection Valve (SSV)	125
Changing configuration or replacing the Jet Weaver	127
Replacing the Seal Wash Pump	130
Releasing a Stuck Inlet Valve	131
Replacing the Pump Heads	134
Disassembling the Pump Head	139
Disassembling the Primary Pump Head	140
Disassembling the Secondary Pump Head	143
Assembling the Pump Head	146
Replacing the Purge Valve Head	154
Replacing Parts of the High Pressure Filter Assembly	157
Installing the Valve Rail Kit	159
Replacing the Main Power Fuses	160
Replacing Module Firmware	162
Preparing the Pump Module for Transport	163

This chapter describes the maintenance of the Agilent 1290 Infinity Binary Pump.




## Introduction to Maintenance

Figure 15 on page 110 shows the main user accessible assemblies of the Agilent 1290 Infinity Binary Pump. These parts can be accessed from the front (simple repairs) and don't require to remove the pump from the system stack.



**Figure 15** Maintenance Parts

Figure 16 on page 111 shows the flow connections between these main assemblies.



**Figure 16** Flow Connections

## Warnings and Cautions

### WARNING

#### Toxic, flammable and hazardous solvents, samples and reagents

The handling of solvents, samples and reagents can hold health and safety risks.

- When working with these substances observe appropriate safety procedures (for example by wearing goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the vendor, and follow good laboratory practice.
  - The volume of substances should be reduced to the minimum required for the analysis.
  - Do not operate the instrument in an explosive atmosphere.
- 

### WARNING

#### Electrical shock

Repair work at the module can lead to personal injuries, e.g. shock hazard, when the cover is opened.

- Do not remove the cover of the module.
  - Only certified persons are authorized to carry out repairs inside the module.
- 

### WARNING

#### Personal injury or damage to the product

Agilent is not responsible for any damages caused, in whole or in part, by improper use of the products, unauthorized alterations, adjustments or modifications to the products, failure to comply with procedures in Agilent product user guides, or use of the products in violation of applicable laws, rules or regulations.

- Use your Agilent products only in the manner described in the Agilent product user guides.
- 

### CAUTION

#### Safety standards for external equipment

- If you connect external equipment to the instrument, make sure that you only use accessory units tested and approved according to the safety standards appropriate for the type of external equipment.
-



## Overview of Maintenance

The following pages describe maintenance (simple repairs) of the module that can be carried out without opening the main cover.

## Cleaning the Module

To keep the module case clean, use a soft cloth slightly dampened with water, or a solution of water and mild detergent.

### WARNING

**Liquid dripping into the electronic compartment of your module can cause shock hazard and damage the module**

- Do not use an excessively damp cloth during cleaning.
  - Drain all solvent lines before opening any connections in the flow path.
-

## Installing Fittings and Capillaries

### WARNING

**Solvent can spray under high pressure.**

→ Observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing), when opening flow path.

---

### CAUTION

**Deformation of fittings and seals**

Liquid drops under high pressure up to 1200 bar act like solid parts. Tightening connections under high pressure can deform or destroy fittings and seals.

→ Never tighten flow connections under pressure.

---

### NOTE

The lifetime of a fitting depends on how firmly it has been tightened; firm tightening reduces the lifetime.

If fitting has been overtightened, replace it.

---


- 1 Install fittings and capillaries.
- 2 Tighten fittings and capillaries.

## Replacing the Shutoff Valve Panel


**When** If a shutoff valve is damaged or the panel needs to be removed for other repair procedures.

Parts required	#	p/n	Description
	1	5067-4124	Shutoff valve
	1	G4220-40004	Shutoff valve panel
	2	G4220-60035	Tubing kit 140 mm SSV to shutoff valve or degassing unit (2 tubes)


**Preparations** In order to avoid leaks, remove tubings from the solvent bottles.



**3** Remove the shutoff valve panel by pulling it downwards



**4** After replacing the panel or after completion of other maintenance, re-install the panel and all flow connections.



## Replacing the Pressure Sensor

**When** No or invalid pressure signal

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-2412	Hex key 2.5 mm, 15 cm long, straight handle
	8710-0510	Wrench open 1/4 — 5/16 inch
		Screwdriver

<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	G4220-60001	Pressure sensor 1200 bar


**Preparations** Turn off pump flow, switch off pump

### NOTE


This procedure describes how to replace the pressure sensor.

In case the cable to the sensor shall be replaced as well, please contact your Agilent service representative.


- 1** Remove capillary connections between the pressure sensor and purge valve.




- 2** Remove the screws that fix the pressure sensor to the chassis.




- 3** Carefully pull out the pressure sensor for about 2 cm. Then unscrew the cable from the pressure sensor.




- 4** Connect the new pressure sensor to the pressure sensor connector.



- 5** Fix the pressure sensor to the instrument chassis.



- 6** Connect the capillaries from the valve to the pressure sensor: connect port 3 to the pressure sensor inlet and port 2 to the outlet.



## Replacing the Inlet valve

**When** If Inlet valve is defective.

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
		Wrench, 14 mm
	G4220-20012	Torque wrench 2 – 25 Nm

<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	G4220-60022	Inlet valve (primary pump head)


### CAUTION

Potential damage of inlet valve


Overtightening the valve will damage it, a lower torque will create leaks.

→ Use a torque wrench (10 Nm) for fixing the inlet valve.

**1** Close the shut off valves to avoid solvent leaks.




**2** Unscrew the tubing at the inlet valve.







- 3** With a 14 mm wrench, unscrew the inlet valve and remove it.




- 4** Install inlet valve and tighten it at 10 Nm with a torque wrench (14 mm).



- 5** Attach the inlet tubing at the inlet valve.



- 6** Open the shut off valves and purge the system to remove air.




## Replacing the Outlet valve

**When** If Outlet valve is defective.


Tools required	p/n	Description
		Wrench, 14 mm
	G4220-20012	Torque wrench 2 – 25 Nm
	G4220-20014	2.5 mm Hex Bit

Parts required	p/n	Description
		Outlet valve 1290 Infinity Pump (primary pump head)
	G4220-20020	Internal gold seal for 1290 Infinity outlet valve


**1** Close the shut off valves to avoid solvent leaks.




**2** Lift up the capillary and remove it from the primary pump head.




- 3** A gold seal between outlet valve and heat exchanger capillary is used for a tight connection. The seal can be replaced separately as needed.




- 4** Unscrew the outlet valve with a 14 mm wrench.



- 5** Insert the outlet valve into the pump head. Using a torque wrench, set 10 Nm and close the outlet valve.




- 6** Insert the heat exchanger capillary into the outlet of the outlet valve. Using a torque wrench with a 2.5 mm hex bit, set 3 Nm and close the hex screw at the top of the outlet.



## 9 Maintenance

### Replacing the Outlet valve

- 7 Open the shut off valves and purge the system to remove air.




## Replacing the Solvent Selection Valve (SSV)


**When** In case of problems with the solvent selection valve

Parts required	#	p/n	Description
	1	G4280-60029	Solvent selection valve
	4	G4220-60035	Tubing kit 140 mm SSV to shutoff valve or degassing unit (2 tubes)

- 1** Close shut-off valve. Remove tubing connections between the SSV and the solvent shut-off valves and the SSV and the degassing unit inlets.




- 2** Push down the SSV panel for removing it.




## 9 Maintenance

### Replacing the Solvent Selection Valve (SSV)

- 3** Remove the connector by pushing up the small clip at the bottom of the connector.



- 4** Install a new SSV by inserting the connector and clipping the SSV panel to the module top panel. Then re-install all tubing connections, open shut-off valve and purge valve.



## Changing configuration or replacing the Jet Weaver

**When** For optimizing the pump configuration to mixing performance or low delay volumes/fast gradients, see chapter *Optimizing Performance*.


<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	8710-0510	¼ inch wrench
		3 mm hex key

<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	G4220-60006	Jet Weaver 35 µL/100 µL
	1	G4220-60012	Jet Weaver 380 µL (OPTIONAL)
	1	G4220-87000	Capillary SST Valve to Jet Weaver 300 mm x 0.17 mm I.D.


## 9 Maintenance

### Changing configuration or replacing the Jet Weaver

- 1 Remove capillary connections from the Jet Weaver.



- 2 Remove the hex screws that fix the Jet Weaver to the pump housing.



#### NOTE


The standard Jet Weaver (Jet Weaver 35  $\mu$ L/100  $\mu$ L (G4220-60006)) has a front and a rear side with different internal volumes (35 / 100  $\mu$ L) that are optimized for a low delay volume or best mixing performance. Please refer to recommendations in the *Agilent 1290 Infinity System Manual*.

The optional Jet Weaver (Jet Weaver 380  $\mu$ L (G4220-60012)) is recommended for applications which are challenging with respect to mixing noise (e.g. TFA applications) and has just one side.




## Changing configuration or replacing the Jet Weaver

- 3** Install new Jet Weaver or flip the Jet Weaver for backside.



- 4** Reinstall the capillary connections.



The inlet at the bottom of the Jet Weaver is connected to the central port of the pump valve by a capillary (length 300 mm, 0.17 mm i.d.). The outlet at the top is connected to the autosampler.


## Replacing the Seal Wash Pump

**When** In case of wear of the seal wash pump


<b>Parts required</b>	<b>p/n</b>	<b>Description</b>
	5042-8507	Peristaltic Pump for Seal Wash

**Preparations** Remove the shutoff valve panel ("[Replacing the Shutoff Valve Panel](#)" on page 116) and flow connections from and to the seal wash pump

- 1** For removing the seal wash pump, press the clips (1) and pull the pump to the front (2).



- 2** Insert the pump clips to the holes in the binary pump housing.




## Releasing a Stuck Inlet Valve

**When** If inlet valve is stuck, or if pump is not generating pressure after being turned off for an extended period of time.


**NOTE**

Before the system is turned off for an extended period of time, it should be flushed with at least 10 % IPA to prevent inlet valves from getting stuck.

- 1** Remove the capillary connection from the outlet of the secondary pump head.




- 2** Unscrew the tubing at the inlet valve.




## 9 Maintenance

### Releasing a Stuck Inlet Valve


- 3** Attach a Luer lock syringe with adapter to the tubing and fill it with solvent.




- 4** Reconnect tubing to inlet valve.




- 5** Unscrew tubing at degassing unit and attach the syringe to it.




- 6** Push solvent with syringe until it comes out at the top of the High Pressure Filter Assembly.



- 7** Detach the syringe and reconnect the tubing into the degassing unit.




- 8** Reinstall the capillary connection to the High Pressure Filter Assembly.



- 9** Purge the system to remove air.

## Replacing the Pump Heads

The Agilent 1290 Infinity Binary Pump has two pump assemblies for two solvent channels A and B which both consist of two pump drives and pump heads. The solvent enters each pump through the primary pump head, is transferred to the secondary pump head and leaves the outlet of the secondary pump head, which is connected to the pump valve.



**Figure 17** Pump head assembly (left)

**When** In case of problems with the pump performance.

<b>Tools required</b>	<b>p/n</b>	<b>Description</b>
	G4220-20012	Torque wrench 2 – 25 Nm
	G4220-20013	4 mm hex bit
	G4220-20015	Adapter ¼ in square to hex

<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	G4220-60200	Pump Head Assembly Channel A (left) with Seal Wash Option
	1	G4220-60210	Pump Head Assembly Channel B (right) with Seal Wash Option
	1	G4220-60400	Pump Head Assembly Channel A (left) without Seal Wash Option
	1	G4220-60410	Pump Head Assembly Channel B (right) without Seal Wash Option

**Preparations**

- Switch off pump at the main power switch
- Remove the front cover

### CAUTION

Limitation of life time

The pump head assembly is an exchange part which cannot be reassembled with standard tools. Disassembling the pump head will strongly limit its life time.

→ Do not disassemble the pump head assembly.

### CAUTION

Damage of connections

Disassembling the flow connection between the both pump heads of one pump head assembly (solvent channel) can damage the connection and cause leaks.

→ Do not disconnect the flow connection between the pump heads.

### CAUTION

Damage of internal parts

→ Do not apply a strong force to the screws of the pump head.

→ Use a torque hex key for that purpose.

## 9 Maintenance

### Replacing the Pump Heads

#### CAUTION

Damage of the pump piston

Removing pump heads in a position other than the maintenance position can damage the pump piston.

→ Before switching off the pump, bring it to the maintenance position.


#### NOTE

This procedure describes the replacement of the left pump head assembly. Similarly, the right pump head assembly can be replaced.

One pump head assembly consists of two pump heads, which are both removed at the same time.


**1** In Lab Advisor go to **Tools > Remove/Install Pump Head** and follow instructions given on the screen.

**2** Close the shut-off valve of the respective pump channel.






- 3** Remove the flow connection between the degassing unit and the primary pump head inlet. Remove the capillary connection at the top of the secondary pump head to the pump valve. If installed, remove the tubings of the seal wash option.



**DO NOT REMOVE** the capillary connection between the pump heads marked by the red X.

- 4** Open the 4 screws holding the pump heads.




**NOTE**

Open all screws step by step, not screw by screw.

## 9 Maintenance

### Replacing the Pump Heads


- 5** Remove the complete pump head assembly by holding it with both hands and pulling it to the front.




#### NOTE

Do not further disassemble the pump head.

- 7** Install flow connections.
- Connect the degassing unit outlet to the inlet of the primary pump head and the outlet of the secondary pump head to the inlet of the purge valve.
  - Channel A (left pump head assembly) is connected to port 4, channel B (right pump head assembly) to port 1.



- 6** Install the pump head assembly by fixing the 4 screws using a 4 mm hex key and a torque wrench set to 5 Nm. Install screws step by step, not screw by screw.



- 8** After the complete reassembly of the module, the pump must be purged. For 10 min, apply a flow of 10 ml/min using solvents suitable for your application and a composition 50:50 for channels A and B, for example 50 % acetonitrile, 50 % water.

## Disassembling the Pump Head

**When** If parts inside the pump head need to be replaced

Tools required	p/n	Description
	G4220-20012	Torque wrench 2 – 25 Nm
	G4220-20013	4 mm hex bit
	G4220-20014	2.5 mm Hex Bit
	G4220-20015	Adapter ¼ in square to hex
	8710-0510	Wrench open 1/4 — 5/16 inch


**Preparations** Remove the pump head assembly as described in “Replacing the Pump Heads” on page 134, steps 1-4.

### NOTE

This procedure describes replacements for the pump heads of channel A. Replacement for channel B can be done accordingly. The primary pump head does not have a heat exchanger. Seal wash parts are optional for both pump heads.

- 1 Open the 2.5 mm hex screw at the top of the primary pump head, which fixes the connection capillary of the heat exchanger. Then lift up the capillary and remove it from the primary pump head.


Heat exchanger with connection capillary



- 2 Remove both pump heads from the link plate by pushing the clips at the rear of the plate.

Link plate

Clips



## Disassembling the Primary Pump Head


### CAUTION

Damage of pump piston


The pump piston is made of  $ZrO_2$ -based ceramic, which is a very hard and resistant material, but it is sensitive to shearing forces from the side.

- Do not try to remove the pump piston from the rear.
- Do not use the piston for removing pump seals.


- 1 Remove the outlet valve at the top of the pump head and the inlet valve at the bottom of the pump head.




- 2 For disassembling the pump head, remove the 4 hex screws at the rear of the pump head.



- 3** Remove the front part of the pump head including pump chamber housing with pump seal and seal holder. If the seal wash option is installed, also remove the support ring with wash seal and gasket.



- 4** Remove the piston from the piston housing by pushing it to the rear, then pull it out from the rear.



## 9 Maintenance

### Disassembling the Primary Pump Head

- 5** Check the pump pistons for scratches, grooves and dents when changing the piston seals.


#### NOTE

Damaged pistons cause micro leaks and will decrease the lifetime of the seals.

- 6** Use the soft plastic side of the insert tool for removing the pump seal from the pump chamber housing.

Insert tool

Pump seal




#### NOTE

Do not use the pump piston for that purpose!

- 7** If the wash seal shall be replaced by a new one, use the steel side of the insert tool for removing it.

Wash seal



## Disassembling the Secondary Pump Head


### CAUTION

Damage of pump piston


The pump piston is made of  $\text{ZrO}_2$ -based ceramic, which is a very hard and resistant material, but it is sensitive to shearing forces from the side.

- Do not try to remove the pump piston from the rear.
- Do not use the piston for removing pump seals.

- 1** Remove the filter assembly at the top of the pump head. See [“Replacing Parts of the High Pressure Filter Assembly”](#) on page 157 for replacing the filter frit.




- 2** For disassembling the pump head, remove the 4 hex screws at the rear of the pump head.




## 9 Maintenance

### Disassembling the Secondary Pump Head

- 3** Remove the front part of the pump head including pump chamber housing with pump seal and seal holder. If the seal wash option is installed, also remove the support ring with wash seal and gasket.



- 4** Remove the piston from the piston housing by pushing it to the rear, then pull it out from the rear.




- 5** Check the pump pistons for scratches, grooves and dents when changing the piston seals.

#### NOTE

Damaged pistons cause micro leaks and will decrease the lifetime of the seals.

- 6** Use the soft plastic side of the insert tool for removing the pump seal from the pump chamber housing.




#### NOTE

Do not use the pump piston for that purpose!



- 7 If the wash seal shall be replaced by a new one, use the steel side of the insert tool for removing it.



# Assembling the Pump Head

**When** Before installing the pump head.

Tools required	p/n	Description
		Pump head alignment tool
	G4220-20012	Torque wrench 2 – 25 Nm
	G4220-20013	4 mm hex bit
	G4220-20014	2.5 mm Hex Bit
	G4220-20015	Adapter ¼ in square to hex
	01018-23702	Insert tool

Parts required	p/n	Description
		PE seals (pack of 2)
	0905-1718	Wash Seal PE
	5062-2484	Gasket, seal wash (pack of 6)

See chapter *"Parts"* for details.

**CAUTION**

Limited life time of the pump head  
Inserting the backup seal wrongly may limit the life time of the pump head.  
→ Please note the correct orientation of the backup seal.

**CAUTION**

Damage of the pump piston  
The pump piston is very sensitive to shearing forces from the side.  
→ Use the alignment piston of the pump head alignment tool for the alignment procedure described below.

## CAUTION

Wrong orientation of pins on support ring

Assembling the pump head without paying attention to the correct orientation of the pins on the support ring can lead to leaks or damage of the piston and pump head.

→ Observe pins on the support ring, which help assembling the parts of the pump head in the correct orientation.

## CAUTION

Damage of the pump head assembly


When installing the pump head assembly, the pump drives need to be in maintenance position, where they are retracted. Using the pump drive in default position will damage the pump head assembly.

→ Bring the pump drive to the maintenance position.


## NOTE

This procedure describes how to assemble the secondary pump head using the pump head alignment tool. Assembling the primary pump head can be done accordingly. The secondary pump head has the heat exchanger capillary, which must fit into the openings of the alignment tool, whereas the primary pump head does not have a heat exchanger.

- 1 Use the soft plastic side of the insert tool for inserting the piston seal to the pump chamber housing.




- 2 If the seal wash option is installed, use the soft plastic side of the insert tool for inserting the wash seal into the support ring.




## 9 Maintenance

### Assembling the Pump Head


- 3** Initial versions of the seal holder and support ring allowed removing the backup ring. In any case, avoid doing so. If the backup seal by accident has been removed from the seal holder, insert it in the correct orientation as shown below. Please ensure that the cut end faces the seal holder.




- 4** If the seal wash option is installed, put the gasket into the support ring and insert the seal holder.



- 5** Insert the support ring and pump head ferrules into the piston housing. Observe the pins on the support ring, which help you assembling the pump head correctly.



- 6** Assemble the pump head by putting the pump chamber housing on top of the support ring. Insert and tighten the two pump head screws.




#### NOTE

Do NOT install the inlet and outlet valves (primary pump head) and the high-pressure filter (secondary pump head) at this stage.


## 9 Maintenance

### Assembling the Pump Head

- 7** Loosely close the 4 screws at the rear of the pump head. The screws will be fixed tightly later.




- 8** Insert the alignment piston of the pump head alignment tool.



- 9** Insert the pump head to the pump head alignment tool. There are openings for the seal wash support ring and heat exchanger of the secondary pump head.

Pump head alignment tool (rear)


Pump head (secondary)




- 10** Close the tool by closing the 3 screws at the connection ring.

Heat exchanger capillary


Tool handle



- 11** Using a torque key, set 5 Nm and fix the central alignment screw.




- 12** Using a torque key, set 5 Nm and fix the 4 screws at the rear of the alignment tool. Tighten screws crosswise.




#### NOTE

This procedure will align pump head parts to their correct positions and close the pump head tightly.

- 13** Open the 3 screws which have closed the pump head alignment tool and take out the aligned pump head. In case the pump head sticks inside the alignment tool, you can use the handle and insert it to the rear of the tool for pushing out the pump head.




- 14** For the primary pump head, install the inlet valve and outlet valve using the torque wrench. Set 10 Nm for the inlet valve and 10 Nm for the outlet valve.




## 9 Maintenance

### Assembling the Pump Head


- 15** For the secondary pump head, assemble and install the high pressure filter assembly using the torque wrench (14 mm hex wrench) set to 16 Nm.




- 16** Remove the alignment piston.



- 17** Insert the pump piston.




- 18** Insert both pump heads to the link plate and make sure that the clips snap in that fix the pump heads.







- 19** Insert the heat exchanger capillary into the outlet of the primary pump head. Using a torque key, set 3 Nm and close the hex screw at the top of the outlet.




- 20** Bring the pump drive to the maintenance position using the Lab Advisor user interface, see “Disassembling the Pump Head” on page 139. Both pump drives must be retracted.



- 21** Install the pump head assembly by fixing the 4 screws using a 4 mm hex key and a torque wrench set to 5 Nm. Install screws step by step, not screw by screw.



- 22** Install flow connections as described in “Replacing the Pump Heads” on page 134.



# Replacing the Purge Valve Head

**When** In case of problems with the purge valve

Parts required	#	p/n	Description
	1	5067-4119	Purge valve head
	1	5067-4655	Capillaries to pump head assemblies channel A and B (2x)
	1	G4220-87000	Capillary SST Valve to Jet Weaver 300 mm x 0.17 mm I.D.
	1	5067-4656	Capillaries to pressure sensor (2x)

**Preparations** Remove all capillary connections to the purge valve

## CAUTION

Potential damage of valve head or malfunction of valve

When the pump is switched on, the valve tag is accessed (read/write) and used for correctly positioning the valve.

If the valve head is replaced while the pump is on, invalid information may be written to the valve head making it unusable, or positioning may be wrong resulting in wrong flow connections inside the valve potentially damaging parts.

→ Switch off the pump before working on the purge valve.


## CAUTION

Bias measurement results


The valve drive contains sensitive optical parts. Pollution of these parts can impair the accurate selection of valve ports and therefore bias measurement results.

→ Protect the optical parts from dust and other pollutions.

- 1** Remove all capillary connections. Then unscrew the black union nut and remove the head of the purge valve by pulling it to the front.



- 2** Put the new valve head onto the valve drive such that the lobe fits to the groove. Screw the valve head onto the valve drive using the union nut.



## 9 Maintenance

### Replacing the Purge Valve Head

#### 3 Install all flow connections:

- Port 1 is connected to the outlet of the secondary pump head of channel B
- Port 2 is connected to the outlet of the pressure sensor
- Port 3 is connected to the inlet of the pressure sensor
- Port 4 is connected the outlet of the secondary pump head of channel A
- Ports 5 and 6 are connected to waste capillaries
- The central port is connected to the Jet Weaver inlet


To pressure sensor

From pump channel A

To Jet Weaver

From pump channel B

To waste



# Replacing Parts of the High Pressure Filter Assembly

**When** For removing blockages and leaks in the high pressure filter assembly. The filter frit in the outlet valve should be replaced regularly depending on the system usage. Other parts are covered by the Agilent Preventive Maintenance (PM) Service.

Tools required	p/n	Description
	G4220-20012	Torque wrench 2 – 25 Nm
	G4220-20015	Adapter ¼ in square to hex 14 mm Hex bit

Parts required	#	p/n	Description
	1	01018-22707	PTFE frits (pack of 5)
	1	5067-4728	Seal cap

CAUTION

Leakage or damaged connection

Opening the outlet of the primary pump head may cause leaks or damage the connection between the pump heads.

→ Do not open the outlet of the primary pump head.

NOTE

This procedure describes replacements for channel A (left pump head assembly) and can be applied accordingly to channel B. In both cases, maintenance is done only at the secondary pump head outlet, which hosts the filter frit.


NOTE

When replacing a PTFE frit, consider replacing the seal cap as well in order to prevent leaks.


## 9 Maintenance

### Replacing Parts of the High Pressure Filter Assembly


- 1** Remove the capillary connection from the outlet of the secondary pump head.




- 2** Use a 14 mm hex wrench for opening the filter assembly of the secondary pump head



- 3** Replace the filter frit and seal cap as desired. Please note the correct orientation of the filter frit.



- 4** Re-install the filter assembly using the torque wrench (14 mm hex bit) set to 16 Nm.



# Installing the Valve Rail Kit

**When** This rail is needed for the installation of external valves


**Tools required** **Description**  
Pozidrive screwdriver #1

Parts required	#	p/n	Description
	1	5067-4634	Valve Rail Kit

**NOTE**

The rail can be installed on the left or right side of the pump. This procedure describes the installation on the left side and applies similarly to the right side.

- 1 The valve rail is fixed to the pump cover by 4 screws. The position of the lower screws is marked on the module cover. First tighten these screws, and then tighten the upper screws.



# Replacing the Main Power Fuses

**When** If the main power LED is off while the main power button is pressed (see “Status indicators” on page 71).

**Tools required** **Description**  
Flat head screwdriver

Parts required	#	p/n	Description
	1	2110-1004	Fuse 10 A t

**Preparations** Switch off the instrument and unplug the main power cable.

**WARNING**

**Fire hazard**

**Using wrong fuses can result in fire hazard.**

- For continued protection against fire hazard, replace line fuses only with the same type and ratings.
- Only use the fuses specified for this instrument. The use of other fuses or materials is prohibited.


**NOTE**

There are more fuses inside the instrument.


If replacing the main power fuse does not resolve the issue, please contact your Agilent service representative.



- 1** Use a screwdriver for removing the main fuse carrier from the compartment next to the main power plug. Remove the fuse from the carrier.



- 2** Install the new fuse 10A to the carrier and insert the carrier to the fuse compartment.



## Replacing Module Firmware

<b>When</b>	<p>The installation of newer firmware might be necessary</p> <ul style="list-style-type: none"><li>• if a newer version solves problems of older versions or</li><li>• to keep all systems on the same (validated) revision.</li></ul> <p>The installation of older firmware might be necessary</p> <ul style="list-style-type: none"><li>• to keep all systems on the same (validated) revision or</li><li>• if a new module with newer firmware is added to a system or</li><li>• if third party control software requires a special version.</li></ul>
-------------	---

<b>Tools required</b>	<b>Description</b>
	LAN/RS-232 Firmware Update Tool
OR	Agilent Lab Advisor software
OR	Instant Pilot G4208A (only if supported by module)

<b>Parts required</b>	<b>#</b>	<b>Description</b>
	1	Firmware, tools and documentation from Agilent web site

<b>Preparations</b>	Read update documentation provided with the Firmware Update Tool.
---------------------	---

**To upgrade/downgrade the module's firmware carry out the following steps:**

- 1 Download the required module firmware, the latest LAN/RS-232 FW Update Tool and the documentation from the Agilent web.
  - [http://www.chem.agilent.com/scripts/cag\\_firmware.asp](http://www.chem.agilent.com/scripts/cag_firmware.asp).
- 2 For loading the firmware into the module follow the instructions in the documentation.

### *Module Specific Information*

There is no specific information for this module.

## Preparing the Pump Module for Transport

**When** If the module shall be transported or shipped.

Parts required	#	p/n	Description
	1	9301-0411	Syringe; Plastic
	1	9301-1337	Syringe adapter
	1	G4220-44000	Protective Foam


**Preparations** Flush both solvent channels with isopropanol.

### CAUTION

Mechanical damage

- For shipping the module, insert the Protective Foam to protect the module from mechanical damage.
- Be careful not to damage tubing or capillary connections while inserting the module in the Protective Foam.

- 1 Remove solvent inlets from solvent reservoirs. Disconnect the solvent tubing from the inlet of primary pump heads for both solvent channels. Use a syringe for removing liquid from the solvent tubings between solvent reservoir, shutoff valve panel, solvent selection valve, degassing unit and pump inlets. Switch the solvent selection valve if applicable.




Inlet primary pump head


## 9 Maintenance

### Preparing the Pump Module for Transport


- 2 Remove tubing and capillary connections to other modules and the solvent cabinet. Remove tubing plugs.



- 3 Remove the shutoff valve panel by pulling it downwards



- 4 You may keep internal tubing and capillary connections.




- 5 Remove cable connections to other modules. Remove the module from the stack.


## 9 Maintenance

### Preparing the Pump Module for Transport

- 6 Carefully insert the Protective Foam to the front part of the instrument. Do not damage any tubing or capillary connections.



- 7 Close the front cover.



- 8 For transport or shipment, put the module and accessory kit to the original shipment box.



## 10 Parts and Materials for Maintenance

Overview of Maintenance Parts [168](#)

Capillaries [169](#)

Solvent Cabinet Kit [170](#)

Seal Wash Option [171](#)

Pump Head Assembly Parts [172](#)

Primary Pump Head Parts [174](#)

Secondary Pump Head Parts [178](#)

Purge Valve [182](#)

Cover Parts [183](#)

Leak Parts [184](#)

Fuses [184](#)

Accessory Kit [185](#)

Others [186](#)

This chapter provides information on parts for maintenance.



## Overview of Maintenance Parts

p/n	Description
G4280-60029	Solvent selection valve
G4220-60006	Jet Weaver 35 µL/100 µL
G4220-60012	Jet Weaver 380 µL (OPTIONAL)
G4220-60001	Pressure sensor 1200 bar
G4220-60200	Pump Head Assembly Channel A (left) with Seal Wash Option
G4220-60210	Pump Head Assembly Channel B (right) with Seal Wash Option
G4220-60400	Pump Head Assembly Channel A (left) without Seal Wash Option
G4220-60410	Pump Head Assembly Channel B (right) without Seal Wash Option




## Capillaries

p/n	Description
5067-4655	Capillaries to pump head assemblies channel A and B (2x)
G4220-87000	Capillary SST Valve to Jet Weaver 300 mm x 0.17 mm I.D.
5067-4656	Capillaries to pressure sensor (2x)

## Solvent Cabinet Kit

p/n	Description
5067-4644	Solvent Cabinet Kit 1290 Infinity Pump includes the following parts:
5065-9981	Solvent cabinet 1200 Infinity, including all plastic parts
9301-1420 (3x)	Solvent bottle, transparent
9301-1450	Solvent bottle, amber
5067-4124 (4x)	Shutoff valve
G4220-60007 (4x)	Bottle Head Assembly
G4220-60035 (4x)	Tubing kit 140 mm SSV to shutoff valve or degassing unit (2 tubes)
G4220-40004	Shutoff valve panel
5042-9967	Tubing clip (set of 5 clips)


## Seal Wash Option



**Figure 18** Seal wash pump

p/n	Description
5042-8507	Seal wash pump cartridge
5065-9978	Tubing, 1 mm i.d., 3 mm o.d., silicone, 5 m, for seal wash option

## Pump Head Assembly Parts




**Figure 19** Pump head assembly parts

<b>Item</b>	<b>p/n</b>	<b>Description</b>
1	G4280-60026	High Pressure Filter Assembly (secondary pump head)
2	G4220-60028	Outlet valve 1290 Infinity Pump (primary pump head)
3	G4220-60022	Inlet valve (primary pump head)

## Primary Pump Head Parts


### Primary Pump Head With Seal Wash (Binary Pump)



**Figure 20** Primary pump head (Binary Pump) with seal wash

<b>Item</b>	<b>p/n</b>	<b>Description</b>
1	0515-1218	Screw M5, 40 mm long
2		Pump Chamber Housing (order pump head)
3	G4220-60022	Inlet valve (primary pump head)
4	G4220-60028	Outlet valve 1290 Infinity Pump (primary pump head)
5	G4220-20020	Internal gold seal for 1290 Infinity outlet valve
6	0905-1420	PE seals (pack of 2)
7	G4220-60016	Seal holder including backup ring
8	5062-2484	Gasket, seal wash (pack of 6)
9	0905-1718	Wash Seal PE
10	G4220-63010	Support Ring (Seal Wash)
11		Pump Head Ferrules (order pump head)
12		Piston Housing (order pump head)
13	5067-5678	Piston 1290 Infinity Pumps, ceramic

## Primary Pump Head Without Seal Wash (Binary Pump)




**Figure 21** Primary pump head (Binary Pump) without seal wash



<b>Item</b>	<b>p/n</b>	<b>Description</b>
1	0515-1218	Screw M5, 40 mm long
2		Pump Chamber Housing (order pump head)
3	G4220-60022	Inlet valve (primary pump head)
4	G4220-60028	Outlet valve 1290 Infinity Pump (primary pump head)
5	G4220-20020	Internal gold seal for 1290 Infinity outlet valve
6	0905-1420	PE seals (pack of 2)
7	G4220-60015	Support ring including backup ring
8		Pump Head Ferrules (order pump head)
9		Piston Housing (order pump head)
10	5067-5678	Piston 1290 Infinity Pumps, ceramic

## Secondary Pump Head Parts


### Secondary Pump Head With Seal Wash (Binary Pump)



**Figure 22** Secondary pump head (Binary Pump) with seal wash

<b>Item</b>	<b>p/n</b>	<b>Description</b>
1	0515-1218	Screw M5, 40 mm long
2	G4220-20003	Pump Head Screw
3		Pump Head Front Plate (order pump head)
4	G4220-81013	Heat Exchanger (secondary pump head only) Channel A
5	0515-5237	Heat Exchanger Screw
6	G4220-20001	Spacer Fitting
7		Pump Chamber Housing (order pump head)
8	0905-1420	PE seals (pack of 2)
9	G4220-60016	Seal holder including backup ring
10	5062-2484	Gasket, seal wash (pack of 6)
11	0905-1718	Wash Seal PE
12	G4220-63010	Support Ring (Seal Wash)
13	(2x)	Pump Head Ferrules (order pump head)
14		Piston Housing (order pump head)
15	5067-5678	Piston 1290 Infinity Pumps, ceramic
16	G4280-60026	High Pressure Filter Assembly (secondary pump head)
17	01018-22707	PTFE frits (pack of 5)
18	5067-4728	Seal cap


## Secondary Pump Head Without Seal Wash (Binary Pump)



**Figure 23** Secondary pump head (Binary Pump) without seal wash

<b>Item</b>	<b>p/n</b>	<b>Description</b>
1	0515-1218	Screw M5, 40 mm long
2	G4220-20003	Pump Head Screw
3		Pump Head Front Plate (order pump head)
4	G4220-81013	Heat Exchanger (secondary pump head only) Channel A
5	0515-5237	Heat Exchanger Screw
6	G4220-20001	Spacer Fitting
7		Pump Chamber Housing (order pump head)
8	0905-1420	PE seals (pack of 2)
9	G4220-60015	Support ring including backup ring
10	(2x)	Pump Head Ferrules (order pump head)
11		Piston Housing (order pump head)
12	5067-5678	Piston 1290 Infinity Pumps, ceramic
13	G4280-60026	High Pressure Filter Assembly (secondary pump head)
14	01018-22707	PTFE frits (pack of 5)
15	5067-4728	Seal cap

# Purge Valve



**Figure 24**   Purge valve parts

Item	p/n	Description
	5067-4119	Purge valve head
1	1535-4857	Stator screws
2	5068-0004	Purge Valve Stator
3	5068-0005	Purge Valve Rotor Seal, 1200 bar
4	1535-4045	Bearing ring

Cover Parts




Figure 25    Cover Parts

Item	p/n	Description
1	5067-4613	Cabinet Kit (Side Covers left/right, Top Cover, Tubing Plug, Base Cover and Leak Seal)
2	5042-9949	Tubing Plug, Plastic
	5042-9972	Tubing grommet
3	5067-4612	Front Cover 1290 Infinity Binary Pump
4	5067-4634	Valve Rail Kit

# Leak Parts



**Figure 26**   Leak funnel

p/n	Description
5041-8388	Leak funnel

# Fuses

p/n	Description
2110-1004 (2x)	Fuse 10 A t




## Accessory Kit

p/n	Description
G4220-68705	Accessory kit
5042-9974	Tubing Flex (1.5 m)
8710-0510	Wrench open 1/4 — 5/16 inch
8710-1924	Wrench open 14 mm
5023-0240	Hex driver, 1/4", slitted
8710-2392	Hex key 4 mm15 cm long T-handle
9301-0411	Syringe; Plastic
9301-1337	Syringe adapter
5067-4657	SST Capillary 300 x 0.17 mm, Pump to Autosampler
5067-4658	SST Capillary 450 x 0.17 mm, Pump to Thermostatted Autosampler
0100-1710	Mounting Tool for Tubing Connections
G4220-67000	Waste tubes
8710-1534	Wrench, 4 mm both ends, open end
5065-9978	Tubing, 1 mm i.d., 3 mm o.d., silicone, 5 m, for seal wash option
5181-1519	CAN cable, Agilent module to module, 1 m
01018-22707	PTFE frits (pack of 5)
01018-23702	Insert tool
5042-9972	Tubing grommet


## Others

Item	p/n	Description
	5067-4699	1290 Infinity pump service kit
1	G4220-20012	Torque wrench 2 – 25 Nm
2	G4220-20013	4 mm hex bit
3	G4220-20014	2.5 mm Hex Bit
4	G4220-20015	Adapter ¼ in square to hex
	G4220-44000	Protective Foam
	5023-0285	Replacement kit for 1290 Infinity pump head alignment tool (piston/handle)
	G4203-68708	HPLC System Tool Kit

The 1290 Infinity pump service kit (5067-4699) includes pump head alignment tool and items 1 – 4 .



**Figure 27** Pump head alignment tool



**Figure 28** Replacement kit for 1290 Infinity pump head alignment tool (piston/handle)



**Figure 29** HPLC System Tool Kit

**10 Parts and Materials for Maintenance**  
**Others**



## 11 Identifying Cables

Cable Overview	190
Analog cables	192
Remote Cables	194
BCD Cables	197
CAN/LAN Cable	199
RS-232 Cable Kit	200
Agilent 1200 Module to Printer	201

This chapter summarizes information on all cables.



# Cable Overview

**NOTE**

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

**Analog cables**

p/n	Description
35900-60750	Agilent module to 3394/6 integrators
35900-60750	Agilent 35900A A/D converter
01046-60105	Analog cable (BNC to general purpose, spade lugs)

**Remote cables**

p/n	Description
03394-60600	Agilent module to 3396A Series I integrators  3396 Series II / 3395A integrator, see details in section <a href="#">“Remote Cables”</a> on page 194
03396-61010	Agilent module to 3396 Series III / 3395B integrators
5061-3378	Remote Cable
01046-60201	Agilent module to general purpose

**BCD cables**

p/n	Description
03396-60560	Agilent module to 3396 integrators
G1351-81600	Agilent module to general purpose

### CAN cables

p/n	Description
5181-1516	CAN cable, Agilent module to module, 0.5 m
5181-1519	CAN cable, Agilent module to module, 1 m


### LAN cables

p/n	Description
5023-0203	Cross-over network cable, shielded, 3 m (for point to point connection)
5023-0202	Twisted pair network cable, shielded, 7 m (for point to point connection)

### RS-232 cables


p/n	Description
G1530-60600	RS-232 cable, 2 m
RS232-61601	RS-232 cable, 2.5 m Instrument to PC, 9-to-9 pin (female). This cable has special pin-out, and is not compatible with connecting printers and plotters. It's also called "Null Modem Cable" with full handshaking where the wiring is made between pins 1-1, 2-3, 3-2, 4-6, 5-5, 6-4, 7-8, 8-7, 9-9.
5181-1561	RS-232 cable, 8 m

# Analog cables




One end of these cables provides a BNC connector to be connected to Agilent modules. The other end depends on the instrument to which connection is being made.

## Agilent Module to 3394/6 Integrators


p/n 35900-60750	Pin 3394/6	Pin Agilent module	Signal Name
	1		Not connected
	2	Shield	Analog -
	3	Center	Analog +

## Agilent Module to BNC Connector

p/n 8120-1840	Pin BNC	Pin Agilent module	Signal Name
	Shield	Shield	Analog -
	Center	Center	Analog +



Agilent Module to General Purpose


p/n 01046-60105	Pin	Pin Agilent module	Signal Name
	1		Not connected
	2	Black	Analog -
	3	Red	Analog +

# Remote Cables



One end of these cables provides a Agilent Technologies APG (Analytical Products Group) remote connector to be connected to Agilent modules. The other end depends on the instrument to be connected to.


## Agilent Module to 3396A Integrators

p/n 03394-60600	Pin 3396A	Pin Agilent module	Signal Name	Active (TTL)
	9	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	3	3 - Gray	Start	Low
	NC	4 - Blue	Shut down	Low
	NC	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	5,14	7 - Red	Ready	High
	1	8 - Green	Stop	Low
	NC	9 - Black	Start request	Low
	13, 15		Not connected	


## Agilent Module to 3396 Series II / 3395A Integrators

Use the cable Agilent module to 3396A Series I integrators (03394-60600) and cut pin #5 on the integrator side. Otherwise the integrator prints START; not ready.


Agilent Module to 3396 Series III / 3395B Integrators

p/n 03396-61010	Pin 33XX	Pin Agilent module	Signal Name	Active (TTL)
	9	1 - White	Digital ground	
	NC	2 - Brown	Prepare run	Low
	3	3 - Gray	Start	Low
	NC	4 - Blue	Shut down	Low
	NC	5 - Pink	Not connected	
	NC	6 - Yellow	Power on	High
	14	7 - Red	Ready	High
	4	8 - Green	Stop	Low
	NC	9 - Black	Start request	Low
	13, 15		Not connected	

### Agilent Module to Agilent 35900 A/D Converters

p/n 5061-3378	Pin 35900 A/D	Pin Agilent module	Signal Name	Active (TTL)
	1 - White	1 - White	Digital ground	
	2 - Brown	2 - Brown	Prepare run	Low
	3 - Gray	3 - Gray	Start	Low
	4 - Blue	4 - Blue	Shut down	Low
	5 - Pink	5 - Pink	Not connected	
	6 - Yellow	6 - Yellow	Power on	High
	7 - Red	7 - Red	Ready	High
	8 - Green	8 - Green	Stop	Low
	9 - Black	9 - Black	Start request	Low

### Agilent Module to General Purpose

p/n 01046-60201	Wire Color	Pin Agilent module	Signal Name	Active (TTL)
	White	1	Digital ground	
	Brown	2	Prepare run	Low
	Gray	3	Start	Low
	Blue	4	Shut down	Low
	Pink	5	Not connected	
	Yellow	6	Power on	High
	Red	7	Ready	High
	Green	8	Stop	Low
	Black	9	Start request	Low

# BCD Cables




One end of these cables provides a 15-pin BCD connector to be connected to the Agilent modules. The other end depends on the instrument to be connected to


## Agilent Module to General Purpose

p/n G1351-81600	Wire Color	Pin Agilent module	Signal Name	BCD Digit
	Green	1	BCD 5	20
	Violet	2	BCD 7	80
	Blue	3	BCD 6	40
	Yellow	4	BCD 4	10
	Black	5	BCD 0	1
	Orange	6	BCD 3	8
	Red	7	BCD 2	4
	Brown	8	BCD 1	2
	Gray	9	Digital ground	Gray
	Gray/pink	10	BCD 11	800
	Red/blue	11	BCD 10	400
	White/green	12	BCD 9	200
	Brown/green	13	BCD 8	100
	not connected	14		
	not connected	15	+ 5 V	Low

**Agilent Module to 3396 Integrators**

p/n 03396-60560	Pin 3396	Pin Agilent module	Signal Name	BCD Digit
	1	1	BCD 5	20
	2	2	BCD 7	80
	3	3	BCD 6	40
	4	4	BCD 4	10
	5	5	BCD0	1
	6	6	BCD 3	8
	7	7	BCD 2	4
	8	8	BCD 1	2
	9	9	Digital ground	
	NC	15	+ 5 V	Low

# CAN/LAN Cable



Both ends of this cable provide a modular plug to be connected to Agilent modules CAN or LAN connectors.

## CAN Cables

p/n	Description
5181-1516	CAN cable, Agilent module to module, 0.5 m
5181-1519	CAN cable, Agilent module to module, 1 m

## LAN Cables

p/n	Description
5023-0203	Cross-over network cable, shielded, 3 m (for point to point connection)
5023-0202	Twisted pair network cable, shielded, 7 m (for point to point connection)

# RS-232 Cable Kit

p/n	Description
G1530-60600	RS-232 cable, 2 m
RS232-61601	RS-232 cable, 2.5 m Instrument to PC, 9-to-9 pin (female). This cable has special pin-out, and is not compatible with connecting printers and plotters. It's also called "Null Modem Cable" with full handshaking where the wiring is made between pins 1-1, 2-3, 3-2, 4-6, 5-5, 6-4, 7-8, 8-7, 9-9.
5181-1561	RS-232 cable, 8 m



## Agilent 1200 Module to Printer

p/n	Description
5181-1529	Cable Printer Serial & Parallel, is a SUB-D 9 pin female vs. Centronics connector on the other end (NOT FOR FW UPDATE). For use with G1323 Control Module.

## **11 Identifying Cables**

Agilent 1200 Module to Printer



## 12 Hardware Information

Firmware Description	204
Electrical Connections	207
Rear View of the Module	208
Interfaces	209
Overview Interfaces	212
Setting the 8-bit Configuration Switch	216
Special Settings	218
Instrument Layout	220
Early Maintenance Feedback	221

This chapter describes the pump in more detail on hardware and electronics.



## Firmware Description

The firmware of the instrument consists of two independent sections:

- a non-instrument specific section, called *resident system*
- an instrument specific section, called *main system*

### Resident System

This resident section of the firmware is identical for all Agilent 1100/1200/1220/1260/1290 series modules. Its properties are:

- the complete communication capabilities (CAN, LAN and RS-232C)
- memory management
- ability to update the firmware of the 'main system'

### Main System

Its properties are:

- the complete communication capabilities (CAN, LAN and RS-232C)
- memory management
- ability to update the firmware of the 'resident system'

In addition the main system comprises the instrument functions that are divided into common functions like

- run synchronization through APG remote,
- error handling,
- diagnostic functions,
- or module specific functions like
  - internal events such as lamp control, filter movements,
  - raw data collection and conversion to absorbance.

## Firmware Updates

Firmware updates can be done using your user interface:

- PC and Firmware Update Tool with local files on the hard disk
- Instant Pilot (G4208A) with files from a USB Flash Disk
- Agilent Lab Advisor software B.01.03 and above

The file naming conventions are:

PPPP\_RVVV\_XXX.dlb, where

PPPP is the product number, for example, 1315AB for the G1315A/B DAD,

R the firmware revision, for example, A for G1315B or B for the G1315C DAD,

VVV is the revision number, for example 102 is revision 1.02,


XXX is the build number of the firmware.

For instructions on firmware updates refer to section *Replacing Firmware* in chapter "Maintenance" or use the documentation provided with the *Firmware Update Tools*.

### NOTE

Update of main system can be done in the resident system only. Update of the resident system can be done in the main system only.

Main and resident firmware must be from the same set.



**Figure 30** Firmware Update Mechanism

#### NOTE

Some modules are limited in downgrading due to their main board version or their initial firmware revision. For example, a G1315C DAD SL cannot be downgraded below firmware revision B.01.02 or to a A.xx.xx.

Some modules can be re-branded (e.g. G1314C to G1314B) to allow operation in specific control software environments. In this case the feature set of the target type are use and the feature set of the original are lost. After re-branding (e.g. from G1314B to G1314C), the original feature set is available again.

All these specific informations are described in the documentation provided with the firmware update tools.

---

The firmware update tools, firmware and documentation are available from the Agilent web.

- <http://www.chem.agilent.com/EN-US/SUPPORT/DOWNLOADS/FIRMWARE/Pages/LC.aspx>

## Electrical Connections


- The CAN bus is a serial bus with high speed data transfer. The two connectors for the CAN bus are used for internal module data transfer and synchronization.
- The REMOTE connector may be used in combination with other analytical instruments from Agilent Technologies if you want to use features such as start, stop, common shut down, prepare, and so on.
- With the appropriate software, the RS-232C connector may be used to control the module from a computer through a RS-232C connection. This connector is activated and can be configured with the configuration switch.
- The power input socket accepts a line voltage of 100 – 240 VAC  $\pm$  10 % with a line frequency of 50 or 60 Hz. Maximum power consumption varies by module. There is no voltage selector on your module because the power supply has wide-ranging capability. There are no externally accessible fuses, because automatic electronic fuses are implemented in the power supply.

### NOTE

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

---

## Rear View of the Module



**Figure 31** Rear of Binary Pump



# Interfaces

The Agilent 1200 Infinity Series modules provide the following interfaces:

**Table 9** Agilent 1200 Infinity Series Interfaces

Module	CAN	LAN/BCD (optional)	LAN (on-board)	RS-232	Analog	APG Remote	Special
<b>Pumps</b>							
G1310B Iso Pump G1311B Quat Pump G1311C Quat Pump VL G1312B Bin Pump G1312C Bin Pump VL 1376A Cap Pump G2226A Nano Pump G5611A Bio-inert Quat Pump	2	Yes	No	Yes	1	Yes	
G4220A/B Bin Pump G4204A Quat Pump	2	No	Yes	Yes	No	Yes	CAN-DC- OUT for CAN slaves
G1361A Prep Pump	2	Yes	No	Yes	No	Yes	CAN-DC- OUT for CAN slaves
<b>Samplers</b>							
G1329B ALS G2260A Prep ALS	2	Yes	No	Yes	No	Yes	THERMOSTAT for G1330B
G1364B FC-PS G1364C FC-AS G1364D FC- $\mu$ S G1367E HiP ALS G1377A HiP micro ALS G2258A DL ALS G5664A Bio-inert FC-AS G5667A Bio-inert Autosampler	2	Yes	No	Yes	No	Yes	THERMOSTAT for G1330B CAN-DC- OUT for CAN slaves
G4226A ALS	2	Yes	No	Yes	No	Yes	

## 12 Hardware Information

### Interfaces

**Table 9** Agilent 1200 Infinity Series Interfaces

Module	CAN	LAN/BCD (optional)	LAN (on-board)	RS-232	Analog	APG Remote	Special
<b>Detectors</b>							
G1314B VWD VL G1314C VWD VL+	2	Yes	No	Yes	1	Yes	
G1314E/F VWD	2	No	Yes	Yes	1	Yes	
G4212A/B DAD	2	No	Yes	Yes	1	Yes	
G1315C DAD VL+ G1365C MWD G1315D DAD VL G1365D MWD VL	2	No	Yes	Yes	2	Yes	
G1321B FLD G1362A RID	2	Yes	No	Yes	1	Yes	
G4280A ELSD	No	No	No	Yes	Yes	Yes	EXT Contact AUTOZERO
<b>Others</b>							
G1170A Valve Drive	2	No	No	No	No	No	Requires a HOST module with on-board LAN (e.g. G4212A or G4220A with minimum firmware B.06.40 or C.06.40) or with additional G1369C LAN Card
G1316A/C TCC	2	No	No	Yes	No	Yes	
G1322A DEG	No	No	No	No	No	Yes	AUX
G1379B DEG	No	No	No	Yes	No	Yes	
G4225A DEG	No	No	No	Yes	No	Yes	
G4227A Flex Cube	2	No	No	No	No	No	
G4240A CHIP CUBE	2	Yes	No	Yes	No	Yes	CAN-DC- OUT for CAN slaves THERMOSTAT for G1330A/B (NOT USED)

**NOTE**

The detector (DAD/MWD/FLD/VWD/RID) is the preferred access point for control via LAN. The inter-module communication is done via CAN.

---

- CAN connectors as interface to other modules
- LAN connector as interface to the control software
- RS-232C as interface to a computer
- REMOTE connector as interface to other Agilent products
- Analog output connector(s) for signal output

## Overview Interfaces

### CAN

The CAN is inter-module communication interface. It is a 2-wire serial bus system supporting high speed data communication and real-time requirement.

### LAN

The modules have either an interface slot for an LAN card (e.g. Agilent G1369B/C LAN Interface) or they have an on-board LAN interface (e.g. detectors G1315C/D DAD and G1365C/D MWD). This interface allows the control of the module/system via a PC with the appropriate control software.

#### NOTE

If an Agilent detector (DAD/MWD/FLD/VWD/RID) is in the system, the LAN should be connected to the DAD/MWD/FLD/VWD/RID (due to higher data load). If no Agilent detector is part of the system, the LAN interface should be installed in the pump or autosampler.

### RS-232C (Serial)

The RS-232C connector is used to control the module from a computer through RS-232C connection, using the appropriate software. This connector can be configured with the configuration switch module at the rear of the module. Refer to *Communication Settings for RS-232C*.

#### NOTE

There is no configuration possible on main boards with on-board LAN. These are pre-configured for

- 19200 baud,
- 8 data bit with no parity and
- one start bit and one stop bit are always used (not selectable).

The RS-232C is designed as DCE (data communication equipment) with a 9-pin male SUB-D type connector. The pins are defined as:

Table 10 RS-232C Connection Table

Pin	Direction	Function
1	In	DCD
2	In	RxD
3	Out	TxD
4	Out	DTR
5		Ground
6	In	DSR
7	Out	RTS
8	In	CTS
9	In	RI




Figure 32 RS-232 Cable

Analog Signal Output

The analog signal output can be distributed to a recording device. For details refer to the description of the module’s main board.

#### APG Remote

The APG Remote connector may be used in combination with other analytical instruments from Agilent Technologies if you want to use features as common shut down, prepare, and so on.

Remote control allows easy connection between single instruments or systems to ensure coordinated analysis with simple coupling requirements.

The subminiature D connector is used. The module provides one remote connector which is inputs/outputs (wired- or technique).

To provide maximum safety within a distributed analysis system, one line is dedicated to **SHUT DOWN** the system's critical parts in case any module detects a serious problem. To detect whether all participating modules are switched on or properly powered, one line is defined to summarize the **POWER ON** state of all connected modules. Control of analysis is maintained by signal readiness **READY** for next analysis, followed by **START** of run and optional **STOP** of run triggered on the respective lines. In addition **PREPARE** and **START REQUEST** may be issued. The signal levels are defined as:

- standard TTL levels (0 V is logic true, + 5.0 V is false),
- fan-out is 10 ,
- input load is 2.2 kOhm against + 5.0 V, and
- output are open collector type, inputs/outputs (wired- or technique).

#### NOTE

All common TTL circuits operate with a 5 V power supply. A TTL signal is defined as "low" or L when between 0 V and 0.8 V and "high" or H when between 2.0 V and 5.0 V (with respect to the ground terminal).

**Table 11** Remote Signal Distribution

Pin	Signal	Description
1	DGND	Digital ground
2	PREPARE	(L) Request to prepare for analysis (for example, calibration, detector lamp on). Receiver is any module performing pre-analysis activities.
3	START	(L) Request to start run / timetable. Receiver is any module performing run-time controlled activities.
4	SHUT DOWN	(L) System has serious problem (for example, leak: stops pump). Receiver is any module capable to reduce safety risk.
5		Not used
6	POWER ON	(H) All modules connected to system are switched on. Receiver is any module relying on operation of others.
7	READY	(H) System is ready for next analysis. Receiver is any sequence controller.
8	STOP	(L) Request to reach system ready state as soon as possible (for example, stop run, abort or finish and stop injection). Receiver is any module performing run-time controlled activities.
9	START REQUEST	(L) Request to start injection cycle (for example, by start key on any module). Receiver is the autosampler.

## Special Interfaces

The module includes a DC-Out (24 VDC) power line that is intended to be used with certain modules that operate as CAN slaves, for example external valves. The line has a limited output of 0.5 A (1.7 A as of August 2011) and is self resetting.

## Setting the 8-bit Configuration Switch


The 8-bit configuration switch is located at the rear of the module. Switch settings provide configuration parameters for LAN, serial communication protocol and instrument specific initialization procedures.

All modules with on-board LAN, e.g. G1315/65C/D, G1314D/E/F, G4212A/B, G4220A/B:

- Default is ALL switches DOWN (best settings).
  - Bootp mode for LAN and
  - 19200 baud, 8 data bit / 1 stop bit with no parity for RS-232
- For specific LAN modes switches 3-8 must be set as required.
- For boot/test modes switches 1+2 must be UP plus required mode.

### NOTE

For normal operation use the default (best) settings.



**Figure 33** Location of Configuration Switch (example shows a G4212A DAD)



**NOTE**

To perform any LAN configuration, SW1 and SW2 must be set to OFF. For details on the LAN settings/configuration refer to chapter LAN Configuration.

**Table 12** 8-bit Configuration Switch (with on-board LAN)

	Mode		Function					
	SW 1	SW 2	SW 3	SW 4	SW 5	SW 6	SW 7	SW 8
LAN	0	0	Link Configuration			Init Mode Selection		
Auto-negotiation			0	x	x	x	x	x
10 MBit, half-duplex			1	0	0	x	x	x
10 MBit, full-duplex			1	0	1	x	x	x
100 MBit, half-duplex			1	1	0	x	x	x
100 MBit, full-duplex			1	1	1	x	x	x
Bootp			x	x	x	0	0	0
Bootp & Store			x	x	x	0	0	1
Using Stored			x	x	x	0	1	0
DHCP			x	x	x	1	0	0
Using Default			x	x	x	0	1	1
TEST	1	1	System					NVRAM
Boot Resident System			1					x
Revert to Default Data (Coldstart)			x	x	x			1

**Legend:**

0 (switch down), 1 (switch up), x (any position)

**NOTE**

When selecting the mode TEST, the LAN settings are: Auto-Negotiation & Using Stored.

**NOTE**

For explanation of "Boot Resident System" and "Revert to Default Data (Coldstart)" refer to "Special Settings" on page 218.

## Special Settings

The special settings are required for specific actions (normally in a service case).

**NOTE** The tables include both settings for modules – with on-board LAN and without on-board LAN. They are identified as LAN and no LAN.

### Boot-Resident

Firmware update procedures may require this mode in case of firmware loading errors (main firmware part).

If you use the following switch settings and power the instrument up again, the instrument firmware stays in the resident mode. It is not operable as a module. It only uses basic functions of the operating system for example, for communication. In this mode the main firmware can be loaded (using update utilities).

**Table 13** Boot Resident Settings (On-board LAN)

Mode Select	SW1	SW2	SW3	SW4	SW5	SW6	SW7	SW8
TEST/BOOT	1	1	1	0	0	0	0	0

### Forced Cold Start

A forced cold start can be used to bring the module into a defined mode with default parameter settings.

CAUTION

Loss of data

Forced cold start erases all methods and data stored in the non-volatile memory. Exceptions are calibration settings, diagnosis and repair log books which will not be erased.

→ Save your methods and data before executing a forced cold start.

If you use the following switch settings and power the instrument up again, a forced cold start has been completed.

**Table 14** Forced Cold Start Settings (On-board LAN)

Mode Select	SW1	SW2	SW3	SW4	SW5	SW6	SW7	SW8
TEST/BOOT	1	1	0	0	0	0	0	1

## Instrument Layout

The industrial design of the module incorporates several innovative features. It uses Agilent's E-PAC concept for the packaging of electronics and mechanical assemblies. This concept is based upon the use of expanded polypropylene (EPP) layers of foam plastic spacers in which the mechanical and electronic boards components of the module are placed. This pack is then housed in a metal inner cabinet which is enclosed by a plastic external cabinet. The advantages of this packaging technology are:

- virtual elimination of fixing screws, bolts or ties, reducing the number of components and increasing the speed of assembly/disassembly,
- the plastic layers have air channels molded into them so that cooling air can be guided exactly to the required locations,
- the plastic layers help cushion the electronic and mechanical parts from physical shock, and
- the metal inner cabinet shields the internal electronics from electromagnetic interference and also helps to reduce or eliminate radio frequency emissions from the instrument itself.

## Early Maintenance Feedback

Maintenance requires the exchange of components which are subject to wear or stress. Ideally, the frequency at which components are exchanged should be based on the intensity of usage of the module and the analytical conditions, and not on a predefined time interval. The early maintenance feedback (**EMF**) feature monitors the usage of specific components in the instrument, and provides feedback when the user-selectable limits have been exceeded. The visual feedback in the user interface provides an indication that maintenance procedures should be scheduled.

### EMF Counters

**EMF counters** increment with use and can be assigned a maximum limit which provides visual feedback in the user interface when the limit is exceeded. Some counters can be reset to zero after the required maintenance procedure.

### Using the EMF Counters

The user-settable **EMF** limits for the **EMF Counters** enable the early maintenance feedback to be adapted to specific user requirements. The useful maintenance cycle is dependent on the requirements for use. Therefore, the definition of the maximum limits need to be determined based on the specific operating conditions of the instrument.

### Setting the EMF Limits

The setting of the **EMF** limits must be optimized over one or two maintenance cycles. Initially the default **EMF** limits should be set. When instrument performance indicates maintenance is necessary, take note of the values displayed by the **EMF counters**. Enter these values (or values slightly less than the displayed values) as **EMF** limits, and then reset the **EMF counters** to zero. The next time the **EMF counters** exceed the new **EMF** limits, the **EMF** flag will be displayed, providing a reminder that maintenance needs to be scheduled.

## **12 Hardware Information**

### **Early Maintenance Feedback**



## 13 LAN Configuration

What You Have to Do First	224
TCP/IP parameter configuration	226
Configuration Switch	227
Initialization mode selection	228
Dynamic Host Configuration Protocol (DHCP)	232
General Information (DHCP)	232
Setup (DHCP)	234
Link configuration selection	236
Automatic configuration with Bootp	237
About Agilent BootP Service	237
How BootP Service Works	238
Situation: Cannot Establish LAN Communication	238
Installation of BootP Service	239
Two Methods to Determine the MAC Address	241
Assigning IP Addresses Using the Agilent BootP Service	242
Changing the IP Address of an Instrument Using the Agilent BootP Service	245
Manual Configuration	247
With Telnet	248
With the Instant Pilot (G4208A)	252
PC and User Interface Software Setup	253
PC Setup for Local Configuration	253
User Interface Software Setup	254

This chapter provides information on connecting the detector to the Agilent ChemStation PC.



# What You Have to Do First

The module has an on-board LAN communication interface.

- 1 Note the MAC (Media Access Control) address for further reference. The MAC or hardware address of the LAN interfaces is a world wide unique identifier. No other network device will have the same hardware address. The MAC address can be found on a label at the rear of the module (see [Figure 35](#) on page 225).




Part number of the pump main board  
Revision Code, Vendor, Year and Week of assembly  
MAC address  
Country of Origin

**Figure 34** MAC-Label



- 2 Connect the instrument's LAN interface (see [Figure 35](#) on page 225) to
  - the PC network card using a crossover network cable (point-to-point) or
  - a hub or switch using a standard LAN cable.



**Figure 35** Location of LAN interfaces and MAC label

## TCP/IP parameter configuration

To operate properly in a network environment, the LAN interface must be configured with valid TCP/IP network parameters. These parameters are:

- IP address
- Subnet Mask
- Default Gateway


The TCP/IP parameters can be configured by the following methods:

- by automatically requesting the parameters from a network-based BOOTP Server (using the so-called Bootstrap Protocol)
- by automatically requesting the parameters from a network-based DHCP Server (using the so-called Dynamic Host Configuration Protocol). This mode requires a LAN-onboard Module or a G1369C LAN Interface card, see [“Setup \(DHCP\)”](#) on page 234
- by manually setting the parameters using Telnet
- by manually setting the parameters using the Instant Pilot (G4208A)

The LAN interface differentiates between several initialization modes. The initialization mode (short form ‘init mode’) defines how to determine the active TCP/IP parameters after power-on. The parameters may be derived from a Bootp cycle, non-volatile memory or initialized with known default values. The initialization mode is selected by the configuration switch, see [Table 16](#) on page 228.

## Configuration Switch

The configuration switch can be accessed at the rear of the module.



**Figure 36** Location of Configuration Switch

The module is shipped with all switches set to OFF, as shown above.

### NOTE

To perform any LAN configuration, SW1 and SW2 must be set to OFF.


**Table 15** Factory Default Settings

Initialization ('Init') Mode	Bootp, all switches down. For details see <a href="#">"Initialization mode selection"</a> on page 228
Link Configuration	speed and duplex mode determined by auto-negotiation, for details see <a href="#">"Link configuration selection"</a> on page 236

# Initialization mode selection

The following initialization (init) modes are selectable:

Table 16 Initialization Mode Switches

	SW 6	SW 7	SW 8	Init Mode
	OFF	OFF	OFF	Bootp
	OFF	OFF	ON	Bootp & Store
	OFF	ON	OFF	Using Stored
	OFF	ON	ON	Using Default
	ON	OFF	OFF	DHCP <sup>1</sup>

<sup>1</sup> Requires firmware B.06.40 or above. Modules without LAN on board, see G1369C LAN Interface Card

## Bootp

When the initialization mode **Bootp** is selected, the module tries to download the parameters from a **Bootp** Server. The parameters obtained become the active parameters immediately. They are not stored to the non-volatile memory of the module. Therefore, the parameters are lost with the next power cycle of the module.





Figure 37 Bootp (Principle)

## Bootp & Store

When **Bootp & Store** is selected, the parameters obtained from a **Bootp** Server become the active parameters immediately. In addition, they are stored to the non-volatile memory of the module. Thus, after a power cycle they are still available. This enables a kind of bootp once configuration of the module.

*Example:* The user may not want to have a **Bootp** Server be active in his network all the time. But on the other side, he may not have any other configuration method than **Bootp**. In this case he starts the **Bootp** Server temporarily, powers on the module using the initialization mode **Bootp & Store**, waits for the **Bootp** cycle to be completed, closes the **Bootp** Server and powers off the module. Then he selects the initialization mode Using Stored and powers on the module again. From now on, he is able to establish the TCP/IP connection to the module with the parameters obtained in that single **Bootp** cycle.




**Figure 38** Bootp & Store (Principle)

### NOTE

Use the initialization mode **Bootp & Store** carefully, because writing to the non-volatile memory takes time. Therefore, when the module shall obtain its parameters from a **Bootp** Server every time it is powered on, the recommended initialization mode is **Bootp**!

#### Using Stored


When initialization mode **Using Stored** is selected, the parameters are taken from the non-volatile memory of the module. The TCP/IP connection will be established using these parameters. The parameters were configured previously by one of the described methods.



**Figure 39** Using Stored (Principle)

## Using Default

When **Using Default** is selected, the factory default parameters are taken instead. These parameters enable a TCP/IP connection to the LAN interface without further configuration, see [Table 17](#) on page 231.



**Figure 40** Using Default (Principle)

### NOTE

Using the default address in your local area network may result in network problems. Take care and change it to a valid address immediately.

**Table 17** Using Default Parameters

IP address:	192.168.254.11
Subnet Mask:	255.255.255.0
Default Gateway	not specified

Since the default IP address is a so-called local address, it will not be routed by any network device. Thus, the PC and the module must reside in the same subnet.

The user may open a Telnet session using the default IP address and change the parameters stored in the non-volatile memory of the module. He may then close the session, select the initialization mode Using Stored, power-on again and establish the TCP/IP connection using the new parameters.

When the module is wired to the PC directly (e.g. using a cross-over cable or a local hub), separated from the local area network, the user may simply keep the default parameters to establish the TCP/IP connection.

### NOTE

In the **Using Default** mode, the parameters stored in the memory of the module are not cleared automatically. If not changed by the user, they are still available, when switching back to the mode Using Stored.

## Dynamic Host Configuration Protocol (DHCP)


### General Information (DHCP)

The Dynamic Host Configuration Protocol (DHCP) is an auto configuration protocol used on IP networks. The DHCP functionality is available on all Agilent HPLC modules with on-board LAN Interface and “B”-firmware (B.06.40 or above).

- G1314D/E/F VWD
- G1315C/D DAD
- G1365C/D MWD
- G4212A/B DAD
- G4220A/B Binary Pump
- G1369C LAN Interface Card
- 1120/1220 LC System

When the initialization mode “DHCP” is selected, the card tries to download the parameters from a DHCP Server. The parameters obtained become the active parameters immediately. They are not stored to the non-volatile memory of the card.

Besides requesting the network parameters, the card also submits its hostname to the DHCP Server. The hostname equals the MAC address of the card, e.g. *0030d3177321*. It is the DHCP server's responsibility to forward the hostname/address information to the Domain Name Server. The card does not offer any services for hostname resolution (e.g. NetBIOS).



**Figure 41** DHCP (Principle)



**NOTE**

- 1 It may take some time until the DHCP server has updated the DNS server with the hostname information.
  - 2 It may be necessary to fully qualify the hostname with the DNS suffix, e.g. *0030d3177321.country.company.com*.
  - 3 The DHCP server may reject the hostname proposed by the card and assign a name following local naming conventions.
-

# 13 LAN Configuration


## Dynamic Host Configuration Protocol (DHCP)

### Setup (DHCP)

**Software required** The modules in the stack must have at least firmware from set A.06.34 and the above mentioned modules B.06.40 or above (must from the same firmware set).

- 1 Note the MAC address of the LAN interface (provided with G1369C LAN Interface Card or Main Board). This MAC address is on a label on the card or at the rear of the main board, e.g. *0030d3177321*.

On the Instant Pilot the MAC address can be found under **Details** in the LAN section.



**Figure 42** LAN Setting on Instant Pilot

- 2 Set the Configuration Switch to DHCP either on the G1369C LAN Interface Card or the main board of above mentioned modules.

**Table 18** G1369C LAN Interface Card (configuration switch on the card)

SW 4	SW 5	SW 6	SW 7	SW 8	Initialization Mode
ON	OFF	OFF	OFF	OFF	DHCP

**Table 19** LC Modules inclusive 1120/1220 (configuration switch at rear of the instrument)

SW 6	SW 7	SW 8	Initialization Mode
ON	OFF	OFF	DHCP

- 3** Turn on the module that hosts the LAN interface.
- 4** Configure your Control Software (e.g. Agilent ChemStation, LabAdvisor, Firmware Update Tool) and use MAC address as host name, e.g. *0030d3177321*.

The LC system should become visible in the control software (see Note in section [“General Information \(DHCP\)”](#) on page 232).

## Link configuration selection


The LAN interface supports 10 or 100 Mbps operation in full- or half-duplex modes. In most cases, full-duplex is supported when the connecting network device - such as a network switch or hub - supports IEEE 802.3u auto-negotiation specifications.

When connecting to network devices that do not support auto-negotiation, the LAN interface will configure itself for 10- or 100-Mbps half-duplex operation.

For example, when connected to a non-negotiating 10-Mbps hub, the LAN interface will be automatically set to operate at 10-Mbps half-duplex.

If the module is not able to connect to the network through auto-negotiation, you can manually set the link operating mode using link configuration switches on the module.

**Table 20** Link Configuration Switches

	SW 3	SW 4	SW 5	Link Configuration
	OFF	-	-	speed and duplex mode determined by auto-negotiation
	ON	OFF	OFF	manually set to 10 Mbps, half-duplex
	ON	OFF	ON	manually set to 10 Mbps, full-duplex
	ON	ON	OFF	manually set to 100 Mbps, half-duplex
	ON	ON	ON	manually set to 100 Mbps, full-duplex

## Automatic configuration with Bootp

### NOTE

All examples shown in this chapter will not work in your environment. You need your own IP-, Subnet-Mask- and Gateway addresses.

---

### NOTE

Assure that the detector configuration switch is set properly. The setting should be either **BootP** or **BootP & Store**, see [Table 16](#) on page 228.

---

### NOTE

Assure that the detector connected to the network is powered off.

---

### NOTE

If the Agilent BootP Service program is not already installed on your PC, then install it from your Agilent ChemStation DVD, located in folder **BootP**.

---

## About Agilent BootP Service

The Agilent BootP Service is used to assign the LAN Interface with an IP address.

The Agilent BootP Service is provided on the ChemStation DVD. The Agilent BootP Service is installed on a server or PC on the LAN to provide central administration of IP addresses for Agilent instruments on a LAN. The BootP service must be running TCP/IP network protocol and cannot run a DHCP server.

## How BootP Service Works

When an instrument is powered on, an LAN Interface in the instrument broadcasts a request for an IP address or host name and provides its hardware MAC address as an identifier. The Agilent BootP Service answers this request and passes a previously defined IP address and host name associated with the hardware MAC address to the requesting instrument.

The instrument receives its IP address and host name and maintains the IP address as long as it is powered on. Powering down the instrument causes it to lose its IP address, so the Agilent BootP Service must be running every time the instrument powers up. If the Agilent BootP Service runs in the background, the instrument will receive its IP address on power-up.

The Agilent LAN Interface can be set to store the IP address and will not lose the IP address if power cycled.

## Situation: Cannot Establish LAN Communication

If a LAN communication with BootP service cannot be established, check the following on the PC:

- Is the BootP service started? During installation of BootP, the service is not started automatically.
- Does the Firewall block the BootP service? Add the BootP service as an exception.
- Is the LAN Interface using the BootP-mode instead of "Using Stored" or "Using Default" modes?

## Installation of BootP Service

Before installing and configuring the Agilent BootP Service, be sure to have the IP addresses of the computer and instruments on hand.

- 1 Log on as Administrator or other user with Administrator privileges.
- 2 Close all Windows programs.
- 3 Insert the Agilent ChemStation software DVD into the drive. If the setup program starts automatically, click **Cancel** to stop it.
- 4 Open Windows Explorer.
- 5 Go to the BootP directory on the Agilent ChemStation DVD and double-click **BootPPackage.msi**.
- 6 If necessary, click the **Agilent BootP Service...** icon in the task bar.
- 7 The **Welcome** screen of the **Agilent BootP Service Setup Wizard** appears. Click **Next**.
- 8 The **End-User License Agreement** screen appears. Read the terms, indicate acceptance, then click **Next**.
- 9 The **Destination Folder** selection screen appears. Install BootP to the default folder or click **Browse** to choose another location. Click **Next**.

The default location for installation is:


C:\Program Files\Agilent\BootPService\

- 10 Click **Install** to begin installation.

## 13 LAN Configuration

### Automatic configuration with Bootp

11 Files load; when finished, the **BootP Settings** screen appears.



**Figure 43** BootP Settings screen

12 In the **Default Settings** part of the screen, if known, you can enter the subnet mask and gateway.

Defaults can be used:

- The default subnet mask is 255.255.255.0
- The default gateway is 192.168.254.11

13 On the **BootP Settings** screen, click **OK**. The **Agilent BootP Service Setup** screen indicates completion.

14 Click **Finish** to exit the **Agilent BootP Service Setup** screen.

15 Remove the DVD from the drive.

This completes installation.

16 Start BootP Service in the Windows® services: On the Windows® desktop click right on **Computer** icon, select **Manage > Services and Applications > Services**. Select the **Agilent BootP Service** and click **Start**.




## Two Methods to Determine the MAC Address

### Enabling logging to discover the MAC address using BootP

If you want to see the MAC address, select the **Do you want to log BootP requests?** check box.

- 1 Open BootP Settings from **Start > All Programs > Agilent BootP Service > EditBootPSettings**.
- 2 In **BootP Settings...** check **Do you want to log BootP requests?** to enable logging.



**Figure 44** Enable BootP logging

The log file is located in

C:\Documents and Settings\All Users\Application Data\Agilent\BootP\LogFile

It contains a MAC address entry for each device that requests configuration information from BootP.

- 3 Click **OK** to save the values or **Cancel** to discard them. The editing ends.
- 4 After each modification of the BootP settings (i.e. **EditBootPSettings**) a stop or start of the BootP service is required for the BootP service to accept changes. See [“Stopping the Agilent BootP Service”](#) on page 245 or [“Restarting the Agilent BootP Service”](#) on page 246.
- 5 Uncheck the **Do you want to log BootP requests?** box after configuring instruments; otherwise, the log file will quickly fill up disk space.

### Determining the MAC address directly from the LAN Interface card label

- 1 Turn off the instrument.
- 2 Read the MAC address from the label and record it.  
The MAC address is printed on a label on the rear of the module.  
See [Figure 34](#) on page 224 and [Figure 35](#) on page 225.
- 3 Turn on the instrument.

## Assigning IP Addresses Using the Agilent BootP Service

The Agilent BootP Service assigns the Hardware MAC address of the instrument to an IP address.

### Determining the MAC address of the instrument using BootP Service

- 1 Power cycle the Instrument.
- 2 After the instrument completes self-test, open the log file of the BootP Service using Notepad.
  - The default location for the logfile is C:\Documents and Settings\All Users\Application Data\Agilent\BootP\LogFile.
  - The logfile will not be updated if it is open.

The contents will be similar to the following:

02/25/10 15:30:49 PM

Status: BootP Request received at outermost layer

Status: BootP Request received from hardware address: 0010835675AC

Error: Hardware address not found in BootPTAB: 0010835675AC

Status: BootP Request finished processing at outermost layer


- 3 Record the hardware (MAC) address (for example, 0010835675AC).
- 4 The Error means the MAC address has not been assigned an IP address and the Tab File does not have this entry. The MAC address is saved to the Tab File when an IP address is assigned.
- 5 Close the log file before turning on another instrument.
- 6 Uncheck the **Do you want to log BootP requests?** box after configuring instruments to avoid having the logfile use up excessive disk space.

## Adding each instrument to the network using BootP

- 1 Follow **Start > All Programs > Agilent BootP Service** and select **Edit BootP Settings**. The BootP Settings screen appears.
- 2 Uncheck the **Do you want to log BootP requests?** once all instruments have been added.

The **Do you want to log BootP requests?** box must be unchecked when you have finished configuring instruments; otherwise, the log file will quickly fill up disk space.

- 3 Click **Edit BootP Addresses...** The **Edit BootP Addresses** screen appears.
- 4 Click **Add...** The **Add BootP Entry** screen appears.



**Figure 45** Enable BootP logging

- 5 Make these entries for the instrument:
  - MAC address
  - Host name, Enter a Hostname of your choice.  
The Host Name must begin with "alpha" characters (i.e. LC1260)
  - IP address
  - Comment (optional)

## 13 LAN Configuration

### Automatic configuration with Bootp

- Subnet mask
- Gateway address (optional)

The configuration information entered is saved in the Tab File.

**6** Click **OK**.

**7** Leave **Edit BootP Addresses** by pressing **Close**.

**8** Exit **BootP Settings** by pressing **OK**.

**9** After each modification of the BootP settings (i.e. EditBootPSettings) a stop or start of the BootP service is required for the BootP service to accept changes. See [“Stopping the Agilent BootP Service”](#) on page 245 or [“Restarting the Agilent BootP Service”](#) on page 246.

**10** Power cycle the Instrument.

OR

If you changed the IP address, power cycle the instrument for the changes to take effect.

**11** Use the PING utility to verify connectivity by opening a command window and typing:

Ping 192.168.254.11 for example.

The Tab File is located at

C:\Documents and Settings\All Users\Application Data\Agilent\BootP\TabFile

## Changing the IP Address of an Instrument Using the Agilent BootP Service

Agilent BootP Service starts automatically when your PC reboots. To change Agilent BootP Service settings, you must stop the service, make the changes, and then restart the service.

### Stopping the Agilent BootP Service

- 1 From the Windows control panel, select **Administrative Tools > Services**. The **Services** screen appears.




Figure 46 Windows Services screen


- 2 Right-click **Agilent BootP Service**.
- 3 Select **Stop**.
- 4 Close the **Services and Administrative Tools** screen.

## 13 LAN Configuration

### Automatic configuration with Bootp

#### Editing the IP address and other parameters in EditBootPSettings

- 1 Select **Start > All Programs > Agilent BootP Service** and select **Edit BootP Settings**. The **BootP Settings** screen appears.
- 2 When the **BootP Settings** screen is first opened, it shows the default settings from installation.
- 3 Press **Edit BootP Addresses...** to edit the Tab File.



**Figure 47** Edit BootP Addresses screen


- 4 In the **Edit BootP Addresses...** screen press **Add...** to create a new entry or select an existing line from the table and press **Modify...** or **Delete** to change the IP address, comment, subnet mask, for example, in the Tab File.  
If you change the IP address, it will be necessary to power cycle the instrument for the changes to take effect.
- 5 Leave **Edit BootP Addresses...** by pressing **Close**.
- 6 Exit BootP Settings by pressing OK.

#### Restarting the Agilent BootP Service

- 1 In the Windows control panel, select **Administrative Tools > Services**. The **Services** screen appears, see [Figure 46](#) on page 245.
- 2 Right-click **Agilent BootP Service** and select **Start**.
- 3 Close the **Services** and **Administrative Tools** screens.

## Manual Configuration

Manual configuration only alters the set of parameters stored in the non-volatile memory of the module. It never affects the currently active parameters. Therefore, manual configuration can be done at any time. A power cycle is mandatory to make the stored parameters become the active parameters, given that the initialization mode selection switches are allowing it.




**Figure 48** Manual Configuration (Principle)

## With Telnet

Whenever a TCP/IP connection to the module is possible (TCP/IP parameters set by any method), the parameters may be altered by opening a Telnet session.


- 1 Open the system (DOS) prompt window by clicking on Windows **START** button and select **"Run..."**. Type **"cmd"** and press OK.
- 2 Type the following at the system (DOS) prompt:
  - **c:\>telnet <IP address>** or
  - **c:\>telnet <host name>**



**Figure 49** Telnet - Starting a session

where <IP address> may be the assigned address from a Bootp cycle, a configuration session with the Handheld Controller, or the default IP address (see ["Configuration Switch"](#) on page 227).

When the connection was established successfully, the module responds with the following:



**Figure 50** A connection to the module is made



- 3
- Type  
? and press enter to see the available commands.

```
C:\WINDOWS\system32\cmd.exe - telnet 134.40.27.95
Agilent Technologies G1315C PP00000024
>?
command syntax      description
-----
?                    display help info
/                    display current LAN settings
ip <x.x.x.x>          set IP Address
sm <x.x.x.x>          set Subnet Mask
gw <x.x.x.x>          set Default Gateway
exit                exit shell
>
```

Figure 51 Telnet Commands

Table 21 Telnet Commands

Value	Description
?	displays syntax and descriptions of commands
/	displays current LAN settings
ip <x.x.x.x>	sets new ip address
sm <x.x.x.x>	sets new subnet mask
gw <x.x.x.x>	sets new default gateway
exit	exits shell and saves all changes

- 4
- To change a parameter follows the style:
  - parameter value, for example:  
**ip 134.40.27.230**Then press [Enter], where parameter refers to the configuration parameter you are defining, and value refers to the definitions you are assigning to that parameter. Each parameter entry is followed by a carriage return.

## 13 LAN Configuration

### Manual Configuration

- 5 Use the “/” and press Enter to list the current settings.

```
C:\WINDOWS\system32\cmd.exe - telnet 134.40.27.95
>/
LAN Status Page
-----
MAC Address   : 0030D30A0838
Init Mode     : Using Stored
TCP/IP Properties
- active -
IP Address    : 134.40.27.95
Subnet Mask   : 255.255.248.0
Def. Gateway  : 134.40.24.1
-----
TCP/IP Status : Ready
-----
Controllers   : no connections
>_
```

information about the LAN interface  
MAC address, initialization mode  
Initialization mode is Using Stored  
active TCP/IP settings  
TCP/IP status - here ready  
connected to PC with controller software (e.g. Agilent  
ChemStation), here not connected

Figure 52 Telnet - Current settings in "Using Stored" mode

- 6 Change the IP address (in this example 134.40.27.99) and type “/” to list current settings.


```
C:\WINDOWS\system32\cmd.exe - telnet 134.40.27.95
>ip 134.40.27.99
>/
LAN Status Page
-----
MAC Address   : 0030D30A0838
Init Mode     : Using Stored
TCP/IP Properties
- active -
IP Address    : 134.40.27.95
Subnet Mask   : 255.255.248.0
Def. Gateway  : 134.40.24.1
- stored -
IP Address    : 134.40.27.99
Subnet Mask   : 255.255.248.0
Def. Gateway  : 134.40.24.1
-----
TCP/IP Status : Ready
-----
Controllers   : no connections
>_
```

change of IP setting to  
Initialization mode is Using Stored  
active TCP/IP settings  
stored TCP/IP settings in non-volatile memory

connected to PC with controller software (e.g. Agilent  
ChemStation), here not connected

Figure 53 Telnet - Change IP settings

- 7 When you have finished typing the configuration parameters, type **exit** and press **Enter** to exit with storing parameters.



```

c:\WINDOWS\system32\cmd.exe
Agilent Technologies G4290A DE00000000
>exit

Connection to host lost.
C:\>_
  
```

**Figure 54** Closing the Telnet Session


#### NOTE

If the Initialization Mode Switch is changed now to “Using Stored” mode, the instrument will take the stored settings when the module is re-booted. In the example above it would be 134.40.27.99.

## With the Instant Pilot (G4208A)

To configure the TCP/IP parameters before connecting the module to the network, the Instant Pilot (G4208A) can be used.

- 1 From the Welcome screen press the **More** button.
- 2 Select **Configure**.
- 3 Press the **DAD** button.
- 4 Scroll down to the LAN settings.



**Figure 55** Instant Pilot - LAN Configuration


- 5 Press the **Edit** button (only visible if not in Edit mode), perform the required changes and press the **Done** button.
- 6 Leave the screen by clicking **Exit**.

## PC and User Interface Software Setup Setup

### PC Setup for Local Configuration

This procedure describes the change of the TCP/IP settings on your PC to match the module's default parameters in a local configuration (see also [“Initialization mode selection”](#) on page 228).

FIXED IP address



**Figure 56** Changing the TCP/IP settings of the PC

## 13 LAN Configuration

### PC and User Interface Software Setup Setup

#### User Interface Software Setup

Install you user interface software according the provided *User Interface Software Setup Guide*.



## 14 Appendix

General Safety Information [256](#)

The Waste Electrical and Electronic Equipment (WEEE) Directive  
(2002-96-EC) [259](#)

Radio Interference [260](#)

Sound Emission [261](#)

Agilent Technologies on Internet [262](#)






This chapter provides addition information on safety, legal and web.



# General Safety Information

## Safety Symbols

**Table 22** Safety Symbols

Symbol	Description
	The apparatus is marked with this symbol when the user should refer to the instruction manual in order to protect risk of harm to the operator and to protect the apparatus against damage.
	Indicates dangerous voltages.
	Indicates a protected ground terminal.
	Indicates eye damage may result from directly viewing the light produced by the deuterium lamp used in this product.
	The apparatus is marked with this symbol when hot surfaces are available and the user should not touch it when heated up.

**WARNING**

**A WARNING**

**alerts you to situations that could cause physical injury or death.**

- Do not proceed beyond a warning until you have fully understood and met the indicated conditions.

**CAUTION**

**A CAUTION**

**alerts you to situations that could cause loss of data, or damage of equipment.**

- Do not proceed beyond a caution until you have fully understood and met the indicated conditions.



## General Safety Information

The following general safety precautions must be observed during all phases of operation, service, and repair of this instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended use of the instrument. Agilent Technologies assumes no liability for the customer's failure to comply with these requirements.

### WARNING

**Ensure the proper usage of the equipment.**

**The protection provided by the equipment may be impaired.**

→ The operator of this instrument is advised to use the equipment in a manner as specified in this manual.

---

## Safety Standards

This is a Safety Class I instrument (provided with terminal for protective earthing) and has been manufactured and tested according to international safety standards.

## Operation

Before applying power, comply with the installation section. Additionally the following must be observed.

Do not remove instrument covers when operating. Before the instrument is switched on, all protective earth terminals, extension cords, auto-transformers, and devices connected to it must be connected to a protective earth via a ground socket. Any interruption of the protective earth grounding will cause a potential shock hazard that could result in serious personal injury. Whenever it is likely that the protection has been impaired, the instrument must be made inoperative and be secured against any intended operation.

Make sure that only fuses with the required rated current and of the specified type (normal blow, time delay, and so on) are used for replacement. The use of repaired fuses and the short-circuiting of fuse holders must be avoided.

Some adjustments described in the manual, are made with power supplied to the instrument, and protective covers removed. Energy available at many points may, if contacted, result in personal injury.

Any adjustment, maintenance, and repair of the opened instrument under voltage should be avoided whenever possible. When inevitable, this has to be carried out by a skilled person who is aware of the hazard involved. Do not attempt internal service or adjustment unless another person, capable of rendering first aid and resuscitation, is present. Do not replace components with power cable connected.

Do not operate the instrument in the presence of flammable gases or fumes. Operation of any electrical instrument in such an environment constitutes a definite safety hazard.

Do not install substitute parts or make any unauthorized modification to the instrument.

Capacitors inside the instrument may still be charged, even though the instrument has been disconnected from its source of supply. Dangerous voltages, capable of causing serious personal injury, are present in this instrument. Use extreme caution when handling, testing and adjusting.

When working with solvents, observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet by the solvent vendor, especially when toxic or hazardous solvents are used.

## The Waste Electrical and Electronic Equipment (WEEE) Directive (2002-96-EC)

### Abstract

The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC), adopted by EU Commission on 13 February 2003, is introducing producer responsibility on all Electric and Electronic appliances from 13 August 2005.

#### NOTE



This product complies with the WEEE Directive (2002/96/EC) marking requirements. The affixed label indicates that you must not discard this electrical/electronic product in domestic household waste.

Product Category: With reference to the equipment types in the WEEE Directive Annex I, this product is classed as a "Monitoring and Control instrumentation" product.

*Do not dispose off in domestic household waste*

To return unwanted products, contact your local Agilent office, or see [www.agilent.com](http://www.agilent.com) for more information.

## Radio Interference

Never use cables other than the ones supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

### Test and Measurement

If test and measurement equipment is operated with equipment unscreened cables and/or used for measurements on open set-ups, the user has to assure that under operating conditions the radio interference limits are still met within the premises.

## Sound Emission

### Manufacturer's Declaration

This statement is provided to comply with the requirements of the German Sound Emission Directive of 18 January 1991.

This product has a sound pressure emission (at the operator position) < 70 dB.

- Sound Pressure  $L_p < 70$  dB (A)
- At Operator Position
- Normal Operation
- According to ISO 7779:1988/EN 27779/1991 (Type Test)

## **Agilent Technologies on Internet**

For the latest information on products and services visit our worldwide web site on the Internet at:

<http://www.agilent.com>

# Index

## 8

8-bit configuration switch  
on-board LAN 216

## A

Agilent Diagnostic software 74  
Agilent Lab Advisor software 74  
Agilent Lab Advisor 74  
Agilent  
on internet 262  
User Interface Software Setup 253  
algae 42  
ambient non-operating temperature 19  
ambient operating temperature 19  
analog signal 213  
analog  
cable 192  
apg remote 214  
automatic configuration with Bootp 237

## B

BCD  
cable 197  
bench space 18  
binary pump shutdown during  
analysis 89  
BootP service  
installation 239  
restart 246  
settings 246  
stop 245  
Bootp  
& Store 229  
automatic configuration 237

initialization modes 228  
using default 231  
using stored 230  
buffer application 42

## C

cable  
analog 192  
BCD 197  
CAN 199  
LAN 199  
remote 194  
RS-232 200  
cables  
analog 190  
BCD 190  
CAN 191  
LAN 191  
overview 190  
remote 190  
RS-232 191  
CAN  
cable 199  
capillaries 115  
cleaning 114  
compensation sensor open 81  
compensation sensor short 82  
condensation 18  
configuration switches 227  
configuration  
low flow rates 42  
one stack 26  
two stack front 29  
two stack rear 30  
two stack 29

current of primary pump drive too  
high 96  
current of secondary pump drive too  
high 96

## D

defect connection between main board and  
pump drive encoder 92  
defect on arrival 24  
degasser's pressure limit violation 87  
delay volume  
description 62  
deliver underrun 92  
delivery checklist 25  
DHCP  
general information 232  
setup 234  
Diagnostic software 74  
dimensions 19  
drive current too high 91  
drive current too low 90  
drive encoder error 95  
drive Encoder failed 90  
drive position limit 94  
drive timeout 91

## E

electrical connections  
descriptions of 207  
EMF  
early maintenance feedback 221  
error messages  
binary pump shutdown during  
analysis 89

- compensation sensor open 81
- compensation sensor short 82
- cover violation 84
- current of primary pump drive too high 96
- current of secondary pump drive too high 96
- defect connection between main board and pump drive encoder 92
- degasser's pressure limit violation 87
- deliver underrun 92
- drive current too high 91
- drive current too low 90
- drive encoder error 95
- drive Encoder failed 90
- drive position limit 94
- drive timeout 91
- fan failed 82
- flow rate limit exceeded 88
- ignition without cover 83, 83
- insufficient power of drive encoder LED 94
- leak sensor open 81
- leak sensor short 80
- leak 83
- lost CAN partner 80
- overcurrent of pump drive 91
- overcurrent of solvent selection valve (SSV) 92
- pressure below lower limit 86
- pressure exceeded upper pressure limit 85
- pressure sensor calibration wrong or missing 98
- pump drive blocked or encoder failed 90
- pump drive encoder defect 93
- pump drive encoder error 97
- pump drive encoder rollover 94
- pump drive error 97
- pump drive stop not found 98

- pump drive stroke blocked 98
- pump error 85
- pump head leak test 107
- purge valve failed 93
- reading of purge valve tag failed 93
- reading the pump encoder tag failed 89
- remote timeout 79
- seal wash pump was missing when tried to turn on 99
- shutdown 79
- solvent counter exceeded limit 87
- system pressure test 104
- target pressure not reached for binary pump degasser 86
- timeout 78
- unknown purge valve type 97
- waste counter limit exceeded 88
- writing the pump encoder tag failed 89
- writing the purge valve tag failed 95
- extra-column volume 62

## F

- fan failed 82
- filter
  - solvent inlet 42
- firmware
  - description 204
  - main system 204
  - resident system 204
  - update tool 205
  - updates 162, 205
  - upgrade/downgrade 162
- fittings 115
- flow connections
  - to pump 35
- flow rate limit exceeded 88
- flow rate
  - minimum 42

- frequency range 19

## G

- general error messages 78

## H

- humidity 19

## I

- initialization mode selection 228
- installation
  - bench space 18
  - power considerations 16
  - seal wash option 39
  - site requirements 15
- instrument layout 220
- insufficient power of drive encoder LED 94
- interfaces 209
- internet 262

## J

- Jet Weaver 127
  - remove capillary connections 128

## L

- Lab Advisor Software
  - tests 102
- LAN
  - automatic configuration with Bootp 237
  - Bootp & Store 229
  - Bootp 228
  - cable 199
  - configuration switch 227
  - configuration 223
  - initialization mode selection 228
  - link configuration selection 236



- manual configuration with telnet 248
- manual configuration 247
- PC and User Interface Software Setup 253
- TCP/IP parameter configuration 226
  - using default 231
  - using stored 230
- what you have to do first 224
- leak sensor open 81
- leak sensor short 80
- leak 83
- line frequency 19
- line voltage 19
- link configuration selection 236
- lost CAN partner 80

## M

- MAC address
  - determine 242
- MAC
  - address 224
- maintenance
  - feedback 221
  - introduction 110
  - replacing firmware 162
- manual configuration
  - of LAN 247
- message
  - cover violation 84
  - ignition without cover 83, 83
  - remote timeout 79
- minimum flow rate 42

## N

- non-operating altitude 19
- non-operating temperature 19

## O

- operating Altitude 19
- operating temperature 19
- optimization
  - achieving higher resolution 65
  - stack configuration 26
- outlet valve 122
- overcurrent of pump drive 91
- overcurrent of solvent selection valve (SSV) 92

## P

- packaging
  - damaged 24
- parts
  - damaged 25
  - missing 25
- PC and User Interface Software Setup 253
- performance specifications 20
  - Optimization 61
- physical specifications 19
- power considerations 16
- power consumption 19
- power cords 17
- power supply indicator 71
- power switch 32
- pressure below lower limit 86
- pressure exceeded upper pressure limit 85
- pressure sensor calibration wrong or missing 98
- pressure sensor
  - replace 118
- pump drive blocked or encoder failed 90
- pump drive encoder defect 93
- pump drive encoder error 97
- pump drive encoder rollover 94

- pump drive error 97
- pump drive stop not found 98
- pump drive stroke blocked 98
- pump error messages 85
- Pump head leak test
  - error message 107
- pump
  - features 10
  - flow connections 35
  - overview 11
  - preparing 42
  - principle 12
- purge valve failed 93

## R

- reading of purge valve tag failed 93
- reading the pump encoder tag failed 89
- remote
  - cable 194
- repairs
  - replacing firmware 162
- resolution
  - Optimization 65
- RS-232C
  - cable 200

## S

- safety class I 257
- safety
  - general information 257
  - standards 19
  - symbols 256
- Seal wash option
  - installation 39
- seal wash pump was missing when tried to turn on 99
- seal wash pump
  - replacing 130
- shutdown 79

## Index

- shutoff valve panel
  - replacing 116
- shutoff valves
  - replacing 116
- site requirements 15
  - power cords 17
- solvent cabinet 37, 42
- solvent counter exceeded limit 87
- solvent inlet filter 42
- special interfaces 215
- special settings
  - boot-resident 218
  - forced cold start 219
- specification
  - physical 19
- specifications 15
  - performance 20
- status indicator 72
- System pressure test
  - error message 104
- system setup and installation
  - optimizing stack configuration 26

## T

- target pressure not reached for binary pump degasser 86
- TCP/IP parameter configuration 226
- telnet
  - configuration 248
- temperature sensor 83
- test functions 70
- tests 102
- timeout 78
- troubleshooting
  - error messages 77, 70
  - status indicators 70, 71

## U

- unknown purge valve type 97

- unpacking 24
- user interfaces 73

## V

- vacuum degassing unit 42
- voltage range 19

## W

- waste counter limit exceeded 88
- weight 19
- writing the pump encoder tag failed 89
- writing the purge valve tag failed 95



## **In This Book**

This manual contains technical reference information about the Agilent 1290 Infinity Binary Pump G4220B.

- introduction and specifications,
- installation,
- using and optimizing,
- troubleshooting and diagnose,
- maintenance,
- parts identification,
- hardware information,
- safety and related information.

© Agilent Technologies 2009, 2010-2012

Printed in Germany  
08/2012



G4220-90005  
Rev. C



**Agilent Technologies**