



CATALYST FOR SUCCESS

## ➔ ANALYSIS OF VOGLIBOSE IN PHARMACEUTICAL FORMULATIONS

### BY HPLC WITH POST-COLUMN DERIVATIZATION

Voglibose is an Alpha-Glucosidase inhibitor widely used for the treatment of diabetes. Alpha-glucosidase inhibitors are agents that delay the glucose absorption at the intestinal level and thereby prevent sudden surge of glucose after a meal. Voglibose is the safest and most effective drug of its class.

Since Voglibose has no UV chromophore, post-column derivatization is employed to produce a fluorescent derivative.

This abstract describes a very sensitive and robust analytical method for the analysis of Voglibose in pharmaceutical tablets. Simple sample preparation and fast analysis time allow for using this method in high throughput environments.

### METHOD

#### *Analytical Conditions*

*Column:* Amino column,  
4.6x250 mm,  
Catalog number 1446250

*Temperature:* 35 °C

*Flow Rate:* 0.6 mL/min

*Mobile Phase:* Sodium phosphate buffer, 20 mM  
pH 6.5 / Acetonitrile (37 : 63)

*Injection Volume:* 50 µL

**Note:** We strongly recommend that the column be flushed with acetonitrile water (80:20) for twenty minutes before making any injections.

#### *Sample Preparation*

Crush 5 tablets and mix with 25 mL of mobile phase. Sonicate for 10 min and filter liquid portion through 0.45 µm filter. Put in HPLC autosampler vial and inject 50 µL.

#### *Post-column Conditions*

*Post-Column System:* Pinnacle PCX

*Heated Reactor Volume:* 3.5 mL

*Temperature:* 100 °C

*Cooling Coil:* 0.15 mL (at room temperature)

*Reagent:* Taurine (6.25 g), Sodium Periodate (2.56 g)  
in 1000 mL of water


*Flow Rate:* 0.6 mL/min

*Detection:* FLD


$\lambda_{ex}$ : 350 nm,  $\lambda_{em}$ : 430 nm

#### *Repeatability studies for different concentration levels.*


	0.5 ppm	100 ppm
Average RT, min	21.25	21.26
RSD, %, N = 6	0.36	0.08
Average Peak Area	9.22	1,562.69
RSD, %, N = 6	1.48	0.79




*Calibration Curve for Analytical Range 0.5-200 ppm*



*Chromatogram of Voglibose Standard, 50 ppm, 50 µL Injection*



*Chromatogram of Voglibose Standard, 0.5 ppm, 50 uL Injection*



*Chromatogram of Voglibose Tablets (Volix™, 0.2 mg), 50 µL Injection*