



Carbohydrate Analysis

Prevail™ Carbohydrate ES HPLC Columns and ELSD



The Challenge of Separating and Detecting Carbohydrates



Silica-based amino columns are powerful sugar separators, but they degrade under conventional aqueous mobile phase conditions required for sugar analyses. The result is loss of retention, poor peak shape, and column voiding. The problem is made more apparent when using ELS and MS detectors because these detectors respond to non-volatile bonded phase loss as baseline noise and/or baseline shift.

Because they lack chromophores, carbohydrates present a particular detection challenge. The traditional approaches to carbohydrate detection are either differential refractive index (RI) or low-wavelength UV (190 – 205nm). The poor sensitivity and gradient incompatibility of RI detection make it less than ideal for carbohydrate analysis by HPLC. Low-wavelength UV also lacks sensitivity and places stringent requirements on the optical properties of the mobile phase. Both techniques are subject to severe baseline shifts and long equilibration times as the detector responds to spectroscopic changes in the solvent.


Benefits of Evaporative Light Scattering Detection (ELSD) for Carbohydrate Analysis

- More sensitive than RI and low-wavelength UV
- Gradient compatible
- Stable baselines, even with gradients
- Solvent choice not limited by UV-cutoff

The Evaporative Light Scattering Detector (ELSD) is a universal HPLC detector that detects any compound less volatile than the mobile phase to low nanogram levels. It has many advantages over RI and low-wavelength UV detection for carbohydrate analyses.

Since carbohydrates are not good UV absorbers, they must be monitored at low wavelengths. This can be hindered by the mobile phase solvent's UV cutoff and poor baseline stability. With ELSD, the mobile phase evaporates before detection, so it is blind to spectroscopic properties and the compositional changes in the mobile phase. The result is quiet, stable baselines, and better sensitivity than with low-wavelength UV detection.

Although refractive index (RI) detection is universal in nature, its lack of sensitivity and gradient incompatibility make it less than ideal for carbohydrate analysis. Evaporation of the mobile phase before detection makes the ELSD gradient compatible for improved resolution and faster separations. In addition, ELSDs can be 10 to 1000 times more sensitive than RI detectors.



Benefits of Prevail™ Carbohydrate ES Columns for Carbohydrate Analysis

- **Highly stable and rugged**
- **Stable baselines, even with gradients**
- **Outstanding resolution and peak shapes**

The Alltech® Prevail™ Carbohydrate ES Column is free from the problems that plague other carbohydrate columns. It is extremely rugged and maintains stable baselines with all detectors, including the ELSD. The polymeric Prevail™ Column shows none of the chemical degradation that limits the lifetime of silica-based amino columns.




Quiet Gradient Baselines

The Prevail™ Carbohydrate ES Column works very well for isocratic separations, especially when used with the ELSD. Gradients with Prevail™ Carbohydrate ES Columns show flat, quiet baselines, free from noise and drift. Not only are gradients possible, they're recommended because they offer powerful selectivity, reduced run times, and enhanced detector sensitivity as peak efficiencies are maximized. Finally, solvents may be chosen for their chromatography advantages instead of their spectral properties.

Exceptional Resolution and Peak Shape

Prevail™ Carbohydrate ES Columns are versatile enough for mono- and oligosaccharides and sugar alcohols. They quickly resolve fructose, glucose, lactose, maltose, and oligosaccharides with high peak efficiencies. Further, Prevail™ Carbohydrate ES Columns produce single peaks for reducing sugars like glucose at ambient temperatures, eliminating the need for column heating.

Enhance Detector Sensitivity With Stable Gradient Baselines



- | | |
|------------------|-----------------------|
| 1. Dextrose | 11. Maltoundecaose |
| 2. Maltose | 12. Maltododecaose |
| 3. Maltotriose | 13. Maltotridecaose |
| 4. Maltotetraose | 14. Maltotetradecaose |
| 5. Maltopentaose | 15. Maltopentadecaose |
| 6. Malthexaose | |
| 7. Maltoseptaose | |
| 8. Maltooctaose | |
| 9. Maltononaose | |
| 10. Maltodecaose | |

Column: 250 x 4.6mm as identified
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	15
%B:	35	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Reduce Run Times Using Alternative Solvents


CHROM 9350

Unmodified Mobile Phase

Column: Prevail™ Carbohydrate ES, 5µm
250 x 4.6mm (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	50
%B:	20	35

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient



1. Fructose
2. Glucose
3. Sucrose
4. Maltose
5. Maltotriose
6. Maltotetraose
7. Maltopentaose


CHROM 9360

Modified Mobile Phase with Acetone

Column: Prevail™ Carbohydrate ES, 5µm
250 x 4.6mm (Part No. 35101)
Mobile Phase: A: Acetonitrile:Acetone (75:25)
B: Water
Gradient:

Time:	0	18
%B:	24	50


Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient





Conquistador® Tequila

CHROM 10585



1. Fructose
2. Mannitol
3. Glucose
4. Sucrose
5. Unknown


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm (Part No. 35101)
Mobile Phase: Acetonitrile : Water (75:25)
Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Orange Juice Profiles

CHROM 9532

CHROM 9531

CHROM 9530




1. Fructose
2. Glucose
3. Sucrose

Column: Prevail™ Carbohydrate ES, 5µm, 53 x 7.0mm Rocket™ (Part No. 35104)
Mobile Phase: Acetonitrile:Water (75:25)
Flow Rate: 1.3mL/min
Detector: ELSD
Column Temp: Ambient

Beer

CHROM 9351



1. Fructose
2. Glucose
3. Maltotetraose
4. Maltotriose
5. Maltopentaose
6. Maltohexaose
7. Maltoheptaose


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	50
%B:	20	35

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Wine - Concord Grape

CHROM 9476



1. Sorbitol
2. Fructose
3. Maltose
4. Inositol


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	60
%B:	17	35

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Splenda® Sweetener

CHROM 10531



1. Sucralose
2. Glucose
3. Maltose
4. Maltodextrins


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	5	10	15	25
%B:	20	20	30	30	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Pepsi One® Cola

CHROM 10530



1. Sucralose
2. Unknown 1
3. Fructose
4. Glucose
5. Unknown 2
6. Sucrose
7. Maltose
8. Maltodextrins


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	5	10	15	25
%B:	20	20	30	30	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient



Kellogg's Frosted Flakes® Cereal Regular vs. 1/3 less sugar



1. Fructose
2. Glucose
3. Sucrose
4. Maltose
5. Oligosaccharides

CHROM
10528


CHROM
10527

Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	5	10	15	25
%B:	20	20	30	30	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Clover Honey



1. Fructose
2. Glucose
3. Sucrose


CHROM
9365

Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	15
%B:	25	40

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Infant Formula



1. Sucrose
2. Lactose


CHROM
10529

Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	5	10	15	25
%B:	20	20	30	30	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Corn Syrup Solids



1. Dextrose
2. Maltose
3. Maltotriose
4. Maltotetraose
5. Maltopentaose
6. Maltohexaose
7. Maltoheptaose
8. Maltooctaose
9. Maltononaose
10. Maltodecaose
11. Maltoundecaose
12. Maltododecaose
13. Maltotridecaose

CHROM
9354


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	15
%B:	35	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Flavor Sample Solution

Chromatogram courtesy of Charles Stewart, Firmenich Inc.



1. Fructose
2. Glucose
3. Sucrose


CHROM
10579

Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	30
%B:	20	29

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Vegetable Fermentation Broth Standards



1. iso-erythritol
2. Fructose
3. Sorbitol
4. Mannitol
5. Glucose
6. Inositol
7. Sucrose
8. Maltitol
9. Maltose
10. Raffinose
11. Maltotriose

CHROM
10004

Column: Prevail™ Carbohydrate ES, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: 0.04% NH₄OH
Gradient:


Time:	0	60
%B:	17	35

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient



Sugar Alcohols

CHROM
9352




1. iso-Erythritol
2. Xylitol
3. Mannitol
4. Maltitol

Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: Acetonitrile:Water (70:30)
Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Cough Syrup

CHROM
10556



1. Fructose
2. Glucose
3. Sucrose
4. Maltose


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	5	10	15	25
%B:	15	15	25	25	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Robitussin® Flu Syrup

CHROM
10557



1. Sucralose
2. Fructose
3. Sorbitol
4. Glucose
5. Sucrose
6. Maltose

Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:


Time:	0	5	10	15	25
%B:	15	15	25	25	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Spiked Urine

Chromatogram courtesy of Rosemary Rose,
Monash University, Dept. of Gastroenterology

CHROM
10584




1. Sucralose
2. Rhamnose
3. Cellobiose

Column: Prevail™ Carbohydrate ES, 5µm, 150 x 4.6mm, (Part No. 35102)
Mobile Phase: Acetonitrile:Water (70:30)
Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Tobacco

CHROM
10526



1. Fructose
2. Glucose
3. Sucrose
4. Unknown


Column: Prevail™ Carbohydrate ES, 5µm, 250 x 4.6mm, (Part No. 35101)
Mobile Phase: A: Acetonitrile B: Water
Gradient:

Time:	0	5	10	15	25
%B:	20	20	30	30	50

Flow Rate: 1.0mL/min
Detector: ELSD
Column Temp: Ambient

Low Level Detection - 10ng on Column

CHROM
10498



1. Fructose
2. Sucrose
3. Raffinose

Column: Prevail™ Carbohydrate ES, 5µm, 150 x 4.6mm, (Part No. 35102)
Mobile Phase: A: Water B: Acetonitrile
Gradient:


Time:	0	10	11	14
%B:	80	50	80	80

Flow Rate: 1.0 mL/min
Detector: ELSD
Column Temp: Ambient




Glucooligosaccharides in Maltodextrin

CHROM
9534



Partially Hydrolyzed Dextran

CHROM
9540



Prevail™ Carbohydrate ES HPLC Columns

Prevail™ Carbohydrate ES HPLC Columns are available in a variety of formats. Choose from analytical, solvent reducing, and preparative formats.

High-Speed Rocket™ Columns Reduce Analysis Time and Solvent Consumption

Separations normally performed on conventional 150mm x 4.6mm columns run faster on Rocket™ Columns and use less mobile phase, saving time and money. The column design allows higher flow rates with reduced back pressures while maintaining good plate counts. At a flow rate of 3mL/min, a simple sugar separation can be accomplished in less than 5 minutes on a 53 x 7mm Rocket™ Column.



Actual Size (33 x 7mm)

PREVAIL™ HPLC COLUMNS			
PACKING	FORMAT	I.D. x LENGTH	PART No.
Carbohydrate ES, 5µm	LC/MS	2.0 x 150mm	35116
	Analytical	4.6 x 150mm	35102
	Analytical	4.6 x 250mm	35101
	Rocket™	7.0 x 53mm	35104
	Prep	7.0 x 100mm	35103
	Prep	10 x 250mm	35115

PREVAIL™ ALL-GUARD™ CARTRIDGES*			
PACKING	I.D. x LENGTH	QTY	PART No.
Carbohydrate ES	4.6 x 7.5mm	3	96435
All-Guard™ Guard Cartridge Holder		1	80101

*All-Guard™ Holder required

PREVAIL™ PREP-GUARDS*			
PACKING	I.D. x LENGTH	QTY	PART No.
Carbohydrate ES	7.6 x 50mm	ea	35105
Direct-Connect™ Column Coupler		ea	28195

*Direct-Connect™ Column Coupler required



To learn more about Alltech® ELSDs, call 1-800-255-8324 or visit our website: www.discoverysciences.com/ELSD

ALLTECH® is a trademark, registered in the United States and/or other countries, of Alltech Associates, Inc. ALL-GUARD™, DIRECT-CONNECT™, PREVAIL™, and ROCKET™ are trademarks of Alltech Associates, Inc. GRACE® and GRACE DAVISON® are trademarks, registered in the United States and/or other countries, of W. R. Grace & Co.-Conn. GRACE DAVISON DISCOVERY SCIENCES™ is a trademark of W. R. Grace & Co.-Conn. CONQUISTADOR® is a registered trademark of Heaven Hill Distilleries, Inc. DOMINICK'S® is a registered trademark of Dominick's Finer Foods, Inc. KELLOGG'S FROSTED FLAKES® is a registered trademark of Kellogg Company. MINUTE MAID® is a registered trademark of Coca-Cola Company. PEPSI ONE® is a registered trademark of PepsCo, Inc. ROBITUSSIN® is a registered trademark of A.H. Robins Company, Incorporated. SPLEND® is a registered trademark of McNeil-PPC, Inc. This trademark list has been compiled using available published information as of the publication date of this brochure and may not accurately reflect current trademark ownership. Alltech Associates, Inc. is a wholly owned subsidiary of W. R. Grace & Co.-Conn. Grace Davison Discovery Sciences is a product group of W. R. Grace & Co.-Conn., which now includes all product lines formerly sold under the Alltech brand. © Copyright 2008 Alltech Associates, Inc. All rights reserved.



Grace Davison Discovery Sciences American Headquarters

2051 Waukegan Road • Deerfield, Illinois 60015-1899

Tel: 847.948.8600 • Fax: 847.948.1078 • Email: discoverysciences@grace.com • Web: www.discoverysciences.com

Other Regional Headquarters:

In Europe:
Brandstraat 12
B-9160 Lokeren, Belgium
Tel: +32 (0)9-340-65 65
Email: discoverysciences.BE@grace.com

In Asia:
19th Floor, K.Wah Center
1010 Hua Hai Zhong Road
Shanghai 200031 PRC
Tel: 86-21-5467 4678
Email: dsbiz.asia@grace.com

In Australia/New Zealand:
30 Brookhollow Avenue
Baulkham Hills
New South Wales 2153, Australia
Tel: +61 2 9865 7177
Email: custcare-syd.discoverysciences@grace.com