

OpenLAB CDS EZChrom Edition A.04.06

Tips and Tricks for GC Users



Jennifer McCulley

GC Software Product Manager

Problem / Solution

Problem:

I have purchase a new 7890B GC with OpenLAB CDS EZChrom A.04.06 replacing a 6890 GC. I want a way to transfer method to the 7890B GC.

Solution:

GC Method Transfer



Instrument Control Driver



Method Migration



6890 GC
Method



Open Method in
OpenLAB CDS
EZChrom Edition,
Download to the 7890B,
then upload to the
software. Save method



Agilent Technologies

OpenLAB CDS EZChrom A.04.06 for GC Users

May 2014

Problem / Solution

Problem:

I want to view data in a simple UI tailored to my needs without changing screen each time I start the software.


I want to view data in my personalized preferences while my colleagues have another preference.

Solution:

OpenLAB Data Analysis-Customized layouts

OpenLAB Data Analysis

- Customized layouts: User selects which windows and items they want to view. They are persisted even if the software is closed.
- Use “up” and “down” arrows to review all the samples in the result set in a few minutes.
- Peak detail: View each peak in the sample in a zoom-in.



Problem / Solution

Problem:

I need to customize labeling the peaks to determine if a peak is present or missing.

I need to view the groups sections within the chromatogram.



Solution:

OpenLAB Data Analysis-Groups Views



Enhanced UI for productive data review

- Display of compounds and groups in the chromatogram
- Display and editing of integration events in the chromatogram
- Customizable peak annotations
- Enhanced layout management
- Global shifting of expected RT's of all compounds
- Configuration of tables
 - Choose the columns you want to see
 - Define number precision to display
 - Re-order columns
 - Show identified and/or unidentified peaks in the injection results
- GC dual channel data review



Problem / Solution

Problem:

I need to use specific data analysis for my result calculation.

- I want to use published Relative Retention Times
- I need to report results in mass%.
- I need to add multipliers in the data analysis to obtain the result for my analysis.

Solution:

OpenLAB Data Analysis-Quantitation features



Compound ID & Quantitation features

- Compound identification by Relative Retention Time Groups of Groups
- Compound multipliers
 - Allows to correct for purity of standard compounds
- Mass% calculation
- Summarize amounts of timed groups and other compounds in a single named group

Compound Table General

☒ External standard ☐ Internal standard

Number of levels: 3

Curve calculation: From average per level

RF definition: Response per amount

Warning: If you change the RF definition, you need to clear your calibration curve otherwise your results will be wrong. Use the "Clear all calibration" RunType for the 1st standard in the injection list to remove the old calibration curve.

☐ Normalize to 100.00 %

Concentration calculation: Amount * Multipliers * Dil. factor

☐ Include ISTD amount

☒ Calculate mass %

Method: Valley2Valley

| # | Type | Name | Signal | Ex |
|---|--------|--------------|---------------------------|----|
| 1 | [Icon] | C10-C25 | Channel A -- Front Signal | |
| 2 | [Icon] | C25-C40 | Channel A -- Front Signal | |
| 3 | [Icon] | C40-Cend | Channel A -- Front Signal | |
| 4 | [Icon] | C10+C40toEnd | | |
| 5 | [Icon] | C10 | Channel A -- Front Signal | |
| 6 | [Icon] | C18 | Channel A -- Front Signal | |
| 7 | [Icon] | C25 | Channel A -- Front Signal | |
| 8 | [Icon] | C40 | Channel A -- Front Signal | |

| Name | RT |
|--|--------|
| <input type="checkbox"/> C10 | 4.285 |
| <input type="checkbox"/> C18 | 7.092 |
| <input type="checkbox"/> C40 | 12.682 |
| <input type="checkbox"/> C25 | 8.366 |
| <input checked="" type="checkbox"/> C10-C25 | 0.000 |
| <input type="checkbox"/> C25-C40 | 0.000 |
| <input checked="" type="checkbox"/> C40-Cend | 0.000 |



Problem / Solution

Problem:

I need to use specific data analysis for my result calculation.

I need to use a specific calibration curve and weighting types.

I need to use calibration curve functions.


Solution:

OpenLAB Data Analysis-Quantitation features



Calibration & Quantitation features

- New calibration curve and weighting types
 - quadratic, logarithmic, exponential, double-logarithmic - based on Area, Height, Area%, Height%
 - 1/Response and 1/Response² weighting
- Enhanced calibration curve functions
 - Enable/disable individual calibration points
 - display of last modification time of a calibration curve
 - display of last calibration time for a compound
 - curve calculation on average point per level or from all individual points



Compound Table General

☒ External standard ☐ Internal standard

Number of levels 3

Curve calculation From average per level

RF definition Response per amount

⚠ If you change the RF definition, you need to clear your calibration curve otherwise your results will be wrong.
 Use the "Clear all calibration" RunType for the 1st standard in the injection list to remove the old calibration curve.

☐ Normalize to 100.00 % ☐ Include ISTD amount

Concentration calculation Amount * Multipliers * Dil. factor ☐ Calculate mass %



Problem / Solution

Problem:

I need to report the peaks in groups of certain compounds.

I need to report peaks in groups of groups.

I need to report the instrument parameters on my sample report.

Solution:

OpenLAB Data Analysis-Intelligent Reporting



A chromatogram plot showing detector response over time. The x-axis is labeled 'Retention Time, Comp' and the y-axis is labeled 'Retention Time, Comp'. There are two distinct peaks. The first peak is smaller and occurs at an earlier retention time. The second peak is larger and occurs at a later retention time. A red arrow points to the baseline of the first peak.

- Baseline is red



Problem / Solution

Problem:

I need to see a trend chart of samples with upper and lower limits

Solution:


OpenLAB Data Analysis-Intelligent Reporting




OpenLAB Data Analysis-Intelligent Reporting

Trend Charting

- Upper and Lower control limits can be entered interactively
- Outliers can be flagged automatically
- 2- Sigma limits & 3- Sigma limits
- Applying Shewart- or Westgard- rules are possible
- x points above or below average
- x points showing up- or down-trend



Reference (average area): 309.8



Problem / Solution

Problem:

I need more flexibility in designing templates for the customized reports.


Solution:

OpenLAB Data Analysis-Intelligent Reporting










OpenLAB Data Analysis-Intelligent Reporting

- New “modified injection” flag (new DA) such as manually modified

| Sequence Summary Report | | | | | | |
|-------------------------|----------|-------------|------|---------|-----------------------|--|
| | | | | | |  Agilent Technologies |
| Line# | Location | Sample Name | Inj# | Cal Lvl | Inj.Date | Manually modified? |
| 8 | Vial 2 | sam2 | 2 | | 1/30/2009 11:48:45 AM | None |
| 11 | | sam3 | 1 | | 1/30/2009 12:33:36 PM | Man. Integ. |
| 12 | | sam4 | 1 | | 1/30/2009 12:42:32 PM | Man. Integ and Compound ID |

Single Sequence Summary [Rev. 0] Printed: 10/15/2013 11:29:57 AM Page 1 of 1

-  Clear Calculation Variable
-  Page No
-  Page No-Total Pages
-  Signatures
-  Total Pages

Font Format

Family : Arial

Aharoni
Algerian
Andalus
Angsana New
AngsanaUPC
Aparajita
Arabic Typesetting
Arial

Size : 14

Style : Normal

Weight : Normal

Decoration : Normal

Preview




AbcXyz

Column Properties




Charts and Customization

- Expression support for chart axis scaling
- Simple peak filter settings for chart control (like table/matrix)
- Snippets to insert page numbers into header or footer



In addition.....

OpenLAB CDS: Scalable in Storage, Lab Management and Administration



OpenLAB CDS: Networked Workstation Configuration

OpenLAB CDS Configurations: Networked workstation



- **Instrument Control:** Local on each workstation
- **Administration:** Central in OpenLAB Core Server Software
- **Storage:** Local on each workstation
- **Benefits:**
 - Central administration for all users, licenses and all user privileges
 - Status information in lab-at-a-glance view from all instruments connected

Fits well with:

- Laboratories with many instruments and few users, budget-controlled
- Laboratories looking for central lab monitoring without putting their instruments on the network



OpenLAB CDS Networked Workstation with Central Storage

OpenLAB CDS Configurations: Networked workstation with OpenLAB ECM



- **Instrument Control:** Local
- **Administration:** Central
- **Storage:** Central (OpenLAB ECM)
- **Benefits:**
 - Central administration (users, licenses, user privileges)
 - Instrument status information (lab-at-a-glance view)
 - Result data available from anywhere
 - Storage in central OpenLAB ECM with database storage

Fits well with:

- Laboratories who need GLP/GMP compliance and central storage
- Laboratories looking for central lab monitoring without putting their instruments on the network and have a need of database storage



Problem / Solution

Problem

- I cannot afford to duplicate work one of my colleagues has already done.
- I need to securely store my data for a set period of time.

Solution

OpenLAB Data Store allows you to centrally store your data, share it with colleagues, backup and archive.


OpenLAB Data Store – Central Data Storage for OpenLAB CDS

Efficient Data Storage and Retrieval

- ✓ Centralized Storage for OpenLAB CDS
- ✓ Efficient search and retrieval of files
- ✓ Local Language Support (Chinese & Japanese)
- ✓ Free of Charge PostgreSQL database

21 CFR Part 11 Compliant

- ✓ Built to support FDA regulations
 - ✓ Data Integrity and Traceability
 - ✓ Electronic Signatures
 - ✓ Archival Capabilities



For More Information ...



For more information, check the Agilent web site
or contact your Agilent sales representative.