

**ASMS 2015**

**WP 071**

Trace Analysis of Dioxins  
and Dioxin-Like PCBs  
Utilizing GC/MS/MS with a  
New High Efficiency  
Source

Jessica Westland; Agilent  
Technologies, Wilmington, DE

## Introduction

Polychlorinated dibenzo-p-dioxins (PCDDs) and polychlorinated dibenzofurans (PCDFs) and polychlorinated biphenyls (PCBs) are highly toxic Persistent Organic Pollutants (POPs) with properties that are detrimental to human health and have been linked to cancer, endocrine disruption, and reproductive disorders. Many of these compounds are not intentionally produced, but are created as by-products of industrial processes, pesticide manufacturing, combustion processes, and other sources. Even though these toxic compounds are not deliberately produced; their chemical properties impart their stability in the environment. These compounds are lipophilic chemicals that accumulate in the fatty tissues of animals that form part of the food chain.

Since these toxic compounds enter and remain in the human food chain, they are monitored and regulated by Food Researchers and Environmental Agencies (such as the

European Commission (EC), US EPA, and WHO).

As of June 2014 the European Union (EU) has instituted new regulations governing the levels of PCDDs, PCDFs, dioxin-like (DL; including non-ortho (NO) and mono-ortho (MO) substituted) PCBs and non-dioxin-like (NDL) PCBs in food and feed. In a new amendment to EU legislations *No. 589/2014* and *No. 709/2014*, the use of GC/MS/MS systems has been accepted as a *confirmatory* technique for checking compliance with maximum levels (ML). Previously, the use of a High Res Mass Spectrometry (HRMS) was needed to confirm and quantify dioxins, due to the ability to identify, confirm, and quantitate the trace levels of dioxins. With the use of a new high-efficiency Electron Ionization (EI) source incorporated in the 7010 Triple Quadrupole GC/MS system, the GC/MS/MS can now confidently detect and quantitate dioxins and dioxin-like PCBs at ultra-trace levels.

## Experimental

### Sample Preparation

The most frequently used methods for the determination of PCDD/PCDF and DL-PCB in foodstuffs and animal feed combine fat extraction (i.e. Soxhlet) with cleanup steps using different column chromatographies (i.e. silica gel coated with sulphuric acid, florisil, alumina, and active carbon)

Manual dioxin sample preparation is tedious and comprehensive; multicolumn automated systems have been made to automate dioxin sample extraction to reduce analysis times and to attempt to reduce costs

The final extract is collected as two separate fractions:

- 1) PCDDs/PCDFs and the NO-PCB congeners
- 2) MO-PCB congeners and the NDL-PCB congeners

Along with the analysis of the native compounds: 13C-ISTDs are required for each individual standard for accurate identification and quantitation; Surrogates are also added prior to extraction to correct for analyte recovery

### Simple straight forward configuration

Analysis utilizes:

Agilent 7010 MS/MS – Increased sensitivity  
Consistent GC method – Increased productivity  
MMI inlet – Flexibility with injection techniques/volumes  
GC column – Validated for dioxin analysis  
RT locked MRM transitions – Optimized & validated

### High-Efficiency EI Source (Agilent 7010 Triple Quadrupole GC/MS)

MS sensitivity depends on the number of ions measured

This new ultra-efficient EI source maximizes the number of ions that are created and transferred out of the source body and into the quadrupole analyzer

GC Conditions			
Column	DB 5MSUI 60 m x 0.25 mmID x 0.25 µm		
Injection port liner	2mm id dimpled splitless liner, UI		
Injection mode	Cold-splitless (compressed air/CO <sub>2</sub> cooled MMI)		
Injection volume	1 µL		
Inlet temperature program	60 °C	0.31 min	
	600 °C/min	330 °C	5 min
	Carrier gas		
	He, constant flow 1.00 mL/min		
Oven program	60 °C	1 min	
	30 °C/min	270 °C	1 min
	2 °C/min	310 °C	0 min
	5 °C/min	350 °C	0.5 min
MS transfer line temperature	350 °C		

Advantages:

- Increased response and better precision at all levels
- Lower detection limits
- More precise ion ratios and better qualitative information

Unit mass resolution allows for sufficient resolution to separate two peaks one mass unit apart

As well as the ability to minimize possible interferences on the analytes of interest

MS/MS parameters	
Electron Energy	70 eV
Tune	eih.s.tune.xml
EM gain	10
MS1 resolution	Unit
MS2 resolution	Unit
Collision Cell	1.5 mL/min N <sub>2</sub> 4 mL/min He
Quant/Qual transitions	Fraction Specific
Dwell times	Fraction Specific
Collision energies	Optimized
Source temperature	350 °C
Quad temperatures	150 °C

gram

aration of the difficult hexa-dioxin/

Peak to Peak Valley = 4.6%      Peak to Peak Valley = 14.3%



between the difficult MO

PCB 123/118 = 2/10 ppb,  $^{13}\text{C}_{12}$  = 10 ppb


$$t = \frac{\bar{x} - \mu_0}{\text{BSD} \times c}$$

---

### A 'Performance LOQ'

icate injections @ lowest calibrat

$$= (\text{Student's } t \text{ value at } 99\% \text{ confidence level}) \times S$$

				DDB	GILP	DDB
--	--	--	--	-----	------	-----

100

PCB - 189	1.095	3.13	8.8
-----------	-------	------	-----

---

PCB-189	1.095	3.13	9	32
---------	-------	------	---	----

100

PCB - 189	0.64	0.02	4
-----------	------	------	---

(TCDD spike at 0.79pg-TEQ/g fat)

PER 1	PER 2	PER 3	PER 4	PER 5
PER 1	PER 2	PER 3	PER 4	PER 5

