

VMA, HVA, 5-HIAA LC-MS/MS Analysis Kit in Human Urine Samples



Good NEWS:

You don't have to deal with
the problems of ECD anymore !

Robust, sensitive and
quantitative analysis
of VMA, HVA, 5-HIAA
NOW AVAILABLE
in 15 minutes

Suitable for any
LC-MS/MS System

Main methods and procedures that have been
selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit

- 15 minutes analysis time
- Sensitive separation with Zivak® VMA HPLC column
- Deuterated internal standards
- No derivatisation, evaporation or SPE process
- Approximately 1000 analysis with one column

VMA, HVA, 5-HIAA




LC-MS/MS Analysis Kit

in Human Urine Samples



Sample Chromatograms and Validation Results

Sample chromatograms are for urine control level 2 of Zivak® VMA LC-MS/MS Analysis Kit



VMA

LOD	($\mu\text{mol/L}$)	0.3
LOQ	($\mu\text{mol/L}$)	0.9
Accuracy	(%)	97.5
Precision Intra-Assay	(%CV)	2.9
Precision Inter-Assay	(%CV)	4.3
Linearity	(R^2)	0.9990

HVA

LOD	($\mu\text{mol/L}$)	0.2
LOQ	($\mu\text{mol/L}$)	0.6
Accuracy	(%)	98.7
Precision Intra-Assay	(%CV)	2.3
Precision Inter-Assay	(%CV)	3.5
Linearity	(R^2)	0.9992

5-HIAA

LOD	($\mu\text{mol/L}$)	0.2
LOQ	($\mu\text{mol/L}$)	0.6
Accuracy	(%)	96.3
Precision Intra-Assay	(%CV)	3.3
Precision Inter-Assay	(%CV)	4.9
Linearity	(R^2)	0.9989

The test results have been validated with an AB Sciex API 4000 LC-MS/MS instrument in Nov. 2014.

Psychoactive Drugs LC-MS/MS Analysis Kit

in Human Urine Samples



Reliable, sensitive and
quantitative detection of
most common drugs
in 12 minutes

Suitable for any
LC-MS/MS System

Main methods and procedures that have been
selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit

- Also *available* with Zivak® Multitasker Full Automated Sample Preparation AND Injection System
- 12 minutes analysis time
- Sensitive separation with Zivak® Psychoactive Drugs HPLC column
- Deuterated internal standards
- Approximately 1000 analysis with one column
- No evaporation, derivatisation or SPE process

Psychoactive Drugs

LC-MS/MS Analysis Kit

in Human Urine Samples



Validation Results


No	Analyte	LOD ($\mu\text{mol/L}$)	LOQ ($\mu\text{mol/L}$)	Accuracy (%)	Precision Intra-Assay (%CV)	Precision Inter-Assay (%CV)	Linearity (R^2)
1	Amphetamine	0.27	0.90	99.24	2.6	2.9	0.9994
2	Methamphetamine	0.31	1.02	100.9	4.8	5.3	0.9998
3	MDA	0.08	0.26	95.32	2.6	2.9	0.9976
4	MDMA	0.12	0.39	98.54	2.9	3.3	0.9991
5	MDEA	0.07	0.23	98.74	3.4	3.7	0.9998
6	Diazepam	0.10	0.33	95.86	2.4	2.6	0.9999
7	Morphine	0.08	0.26	100.5	6.7	6.9	0.9999
8	Codeine	0.31	1.02	92.10	4.4	4.7	0.9969
9	Cocaine	0.09	0.30	96.84	2.8	3.1	0.9989
10	Alprazolam	0.09	0.30	96.96	2.8	3.2	0.9996
11	Delta-9-THC	0.50	1.65	98.82	4.1	4.5	0.9998
12	Clonazepam	0.16	0.53	97.42	2.9	3.3	0.9998
13	Flunitrazepam	0.13	0.43	98.76	4.4	4.7	0.9996
14	6-acetylmorphine	0.08	0.26	99.12	2.6	2.8	0.9986
15	Heroin	0.12	0.40	98.10	6.3	6.7	0.9987

The test results have been validated with an AB Sciex API 4000 LC-MS/MS instrument in Nov. 2014.

Total Homocysteine LC-MS/MS Analysis Kit

in Human Plasma Samples

METABOLIC DISEASES



2 MINUTES

**NOW *cost effective,*
sensitive and
quantitative detection**

Suitable for any
LC-MS/MS System

Main methods and procedures that have been selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit


- 2 minutes analysis time
- Sensitive separation with Zivak® Total Homocysteine HPLC column
- Deuterated internal standards
- No evaporation or SPE process
- Approximately 1000 analysis with one column

Total Homocysteine LC-MS/MS Analysis Kit in Human Plasma Samples



Sample Chromatogram and Validation Results

Sample chromatogram is for plasma control level 1 of Zivak® Total Homocysteine LC-MS/MS Analysis Kit



Homocysteine		
LOD	($\mu\text{mol/L}$)	0.06
LOQ	($\mu\text{mol/L}$)	0.18
Accuracy	(%)	98.5
Precision Intra-Assay	(%CV)	3.5
Precision Inter-Assay	(%CV)	3.8
Linearity	(R^2)	0.997

The test results have been validated with an AB Sciex API 4000 LC-MS/MS instrument in Nov. 2014.

LC-MS/MS Analysis Kit for the in-vitro Diagnostic Quantitative Determination of Phenylketonuria in Human Serum & Plasma and Dried Blood Samples

METABOLIC DISEASES



2 MINUTES

NOW cost effective,
accurate, sensitive
and quantitative
detection

Suitable for any
LC-MS/MS System

Main methods and procedures that have been
selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit


- 2 minutes analysis time
- Minimum sample preparation process
- Sensitive separation with Zivak® PKU HPLC column
- Approximately 1000 analysis with one column

LC-MS/MS Analysis Kit for the in-vitro Diagnostic Quantitative Determination of Phenylketonuria in Human Serum & Plasma and Dried Blood Samples

METABOLIC DISEASES


Sample Chromatograms and Validation Results

Sample chromatograms are for DBS calibrator of Zivak® PKU LC-MS/MS Analysis Kit




Phenylalanine (PHE)

LOD	($\mu\text{mol/L}$)	0.18
LOQ	($\mu\text{g/L}$)	0.54
Accuracy	(%)	95
Precision Intra-Assay	(%CV)	2.3
Precision Inter-Assay	(%CV)	2.6
Linearity	(R^2)	0.998



Tryptophan (TRP)

LOD	($\mu\text{mol/L}$)	0.20
LOQ	($\mu\text{g/L}$)	0.60
Accuracy	(%)	94
Precision Intra-Assay	(%CV)	2.1
Precision Inter-Assay	(%CV)	2.5
Linearity	(R^2)	0.996



Tyrosine (TYR)

LOD	($\mu\text{mol/L}$)	0.21
LOQ	($\mu\text{g/L}$)	0.63
Accuracy	(%)	96
Precision Intra-Assay	(%CV)	1.9
Precision Inter-Assay	(%CV)	2.3
Linearity	(R^2)	0.997

The test results have been validated with an AB Sciex API 4000 LC-MS/MS instrument in Nov. 2014.



HPLC-UV ANALYSIS KITS

25-OH Vitamin D2/D3
in Human Serum Samples

Quantitative Determination of PKU
in Human Serum & Plasma and
Dried Blood Samples

VD-200 Vitamin D2/D3 Analyser

World only fully automated analyser
with a barcode reader,
centrifuge, vortex,
dual parallel LC and more



3.5 MINUTES



25-OH Vitamin D2/D3 HPLC-UV Analysis Kit

in Human Serum Samples



NOW *cost effective*
HPLC analysis kit
with reliable accuracy
in 3.5 *minutes* available

Suitable for any
HPLC System



Main methods and procedures that have been selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit

- Also *available* with Zivak® Multitasker Full Automated Sample Preparation AND Injection System
- 3.5 minutes analysis time
- Sensitive separation with Zivak® 25-OH Vitamin D HPLC column
- Approximately 1000 analysis with one column


*Reference upon request

25-OH Vitamin D2/D3 HPLC-UV Analysis Kit in Human Serum Samples



Sample Chromatogram and Validation Results

Sample chromatogram is for serum calibrator of Zivak® Vitamin D2/D3 HPLC Analysis Kit



25-OH Vitamin D3		
LOD	(µg/L)	0.8
LOQ	(µg/L)	2.4
Accuracy	(%)	98
Precision Intra-Assay	(%CV)	2.7
Precision Inter-Assay	(%CV)	3.1
Linearity (R ²)	(0-200 µg/L)	0.9995

25-OH Vitamin D2		
LOD	(µg/L)	0.8
LOQ	(µg/L)	2.4
Accuracy	(%)	97
Precision Intra-Assay	(%CV)	2.4
Precision Inter-Assay	(%CV)	2.8
Linearity (R ²)	(0-200 µg/L)	0.9989

The results have been validated with a Zivak® COMET-4100 HPLC-UV system in Nov. 2014.

Tested by **NIST*** with an accuracy of 96%

HPLC-UV Analysis Kit for the in-vitro Diagnostic Quantitative Determination of Phenylketonuria

in Human Serum & Plasma and Dried Blood Samples



NOW *cost effective*,
accurate, sensitive
and quantitative
HPLC analysis
in 8 minutes

Suitable for any
HPLC System

Main methods and procedures that have been
selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit


- 8 minutes analysis time
- Minimum sample preparation process
- Sensitive separation with Zivak® PKU HPLC column
- Approximately 1000 analysis with one column

HPLC-UV Analysis Kit for the in-vitro Diagnostic Quantitative Determination of Phenylketonuria in Human Serum & Plasma and Dried Blood Samples

METABOLIC DISEASES

Sample Chromatogram and Validation Results

Sample chromatogram is for DBS calibrator of Zivak® PKU HPLC-UV Analysis Kit



Phenylalanine (PHE)		
LOD	($\mu\text{mol/L}$)	0.35
LOQ	($\mu\text{mol/L}$)	1.05
Accuracy	(%)	97
Precision Intra-Assay	(%CV)	3.9
Precision Inter-Assay	(%CV)	4.2
Linearity	(R^2)	0.997

Tryptophan (TRP)		
LOD	($\mu\text{mol/L}$)	0.60
LOQ	($\mu\text{mol/L}$)	1.80
Accuracy	(%)	98
Precision Intra-Assay	(%CV)	3.4
Precision Inter-Assay	(%CV)	3.8
Linearity	(R^2)	0.996

Tyrosine (TYR)		
LOD	($\mu\text{mol/L}$)	0.85
LOQ	($\mu\text{mol/L}$)	2.55
Accuracy	(%)	97
Precision Intra-Assay	(%CV)	4.9
Precision Inter-Assay	(%CV)	5.3
Linearity	(R^2)	0.998

The results have been validated with a Zivak® COMET- 4100 HPLC-UV system in Nov. 2014.



HAEMOGLOBIN ANALYSIS KITS & ANALYSER

Haemoglobin A1c
in Whole Blood Samples

Haemoglobin Variants
in Whole Blood Samples

HB-100 HbA1c/Variant Analyser

1,6 analysis run time with sample inversion,
rotational barcode reading and
200 sample capacity.



Haemoglobin A1c HPLC-UV Analysis Kit

in Whole Blood Samples

HAEMOGLOBIN AND DIABETES



1.6 MINUTES

analysis of
HbA1c with short program
and minimum sample
preparation process

*Automated sample preparation
AND injection
from whole blood samples*

Main methods and procedures that have been
selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit


- Available with Zivak® Haemoglobin Analyser
- 1.6 minutes analysis time
- Minimum sample preparation process
- Sensitive separation with Zivak® Haemoglobin HPLC column
- Approximately 1000 analysis with one column

Haemoglobin A1c HPLC-UV Analysis Kit in Whole Blood Samples

HAEMOGLOBIN AND DIABETES

Sample Chromatogram and Validation Results

Sample chromatogram is for whole blood control level 2 of Zivak® Hb A1c HPLC-UV Analysis Kit



Component	Retention time
Hb ₀	1.347

Component	Retention time
HbF	0.623

Component	Retention time
HbA _{1c}	0.788

(%)

The test results have been validated with a Zivak® Haemoglobin Analyser in Nov. 2014.

Haemoglobin Variants HPLC-UV Analysis Kit

in Whole Blood Samples



5 minutes analysis of
A1c, A0, A2, F, D, E, S, C variants
with minimum
sample preparation process

*Automated sample preparation
AND injection
from whole blood samples*

Main methods and procedures that have been
selected are based on EN ISO 13485 and 98/79/EC.

Properties of the analysis kit


- Available with Zivak® Haemoglobin Analyser
- Minimum sample preparation process
- 5 minutes analysis time
- Sensitive separation with Zivak® Haemoglobin HPLC column
- Approximately 1000 analysis with one column

Haemoglobin Variants HPLC-UV Analysis Kit in Whole Blood Samples



Sample Chromatogram and Validation Results

Sample chromatogram for whole blood control level 2 of Zivak® Hb Variants HPLC Analysis Kit



Component	Retention time	Component	Retention time	Component	Retention time	Component	Retention time
HbA1c	2.4	HbF	3.4	HbA ₀	3.9	HbAE	4.2

Component	Retention time	Component	Retention time	Component	Retention time	Component	Retention time
HbA1 ₂	4.5	HbD	4.7	HbS	5.535	HbC	2.2

The test results have been validated with a Zivak® Haemoglobin Analyser in Nov. 2014.



ZIVA[®]
TECHNOLOGIES