## **CONFIRM OIL INTEGRITY** AND PREVENT CATASTROPHIC FAILURE

The Measure of Confidence

### Agilent Transformer Oil Gas Analyzer (TOGA)

The performance and longevity of electrical transformers depend upon the stability of the oil used as an insulator and coolant. Normal transformer operation subjects the oil to electrical and mechanical stresses, causing aging, oxidation, vaporization, electrolytic action, and decomposition. These can change the oil's chemical properties and result in gas formation.

Analyzing these dissolved gases provides crucial diagnostic information about the transformer's current and future stability — and helps you determine whether a transformer should be decommissioned.

# Reliably characterize oil composition and integrity immediately after installation

Based on Agilent 7890B GC system, **Agilent Transformer Oil Gas Analyzers** are factory-configured and chemically tested to help you detect individual gas components – and measure their ratios – to predict and prevent transformer failure.



### Agilent Transformer Oil Gas Analyzers include innovative technology and reflect our stringent quality control process. Systems include:

#### Factory

- · System setup and leak testing
- · Instrument checkout
- · Installation of appropriate columns
- · Factory-run checkout method using application checkout mix

#### Delivery

- · Instrument manual for running the method
- CD-ROM with method parameters and checkout data files for easy out-of-the-box operation
- Application related consumables included no separate ordering required
- · Easy consumables re-ordering information

#### Installation

- Duplicate factory checkout with checkout sample onsite by factory-trained support engineer
- · Optional application startup assistance



### **Agilent Technologies**

## Perform fast, unattended analysis of transformer oil gas using these built-in features:

- A robust method for measuring decomposition products, including H<sub>2</sub>, O<sub>2</sub>, CO, CO<sub>2</sub>, and C<sub>1</sub>-C<sub>3</sub>.
- Compliance with ASTM D3612 including customizable data output reports.
- **Rigorous, highly efficient packing technology,** used in Agilent J&W Packed GC columns, ensures column-to-column reproducibility and separation efficiency.
- **Unattended operation:** Analyzers are configured around the Agilent 7890B GC, and interfaced to the Model 7697A Headspace Sampler.

#### **Ordering information:**

Part Number	Configured per	Target Analytes
G3445 Series #571	ASTM D3612-C	$H_2,O_2,N_2,CO,CO_2,CH_4,C_2$ (ethane, ethylene, acetylene), $C_3$ (propane, propylene), and $C_4$ (1-butene)
7890-0047	ASTM D3612-A	$\rm H_2,~O_2,~N_2,~CO,~CO_2,~CH_4,~C_2$ (ethane, ethylene, acetylene), $\rm C_3$ (propane, propylene), and $\rm C_4$ (1-butene)
7890-0552*	ASTM D3612-C	$\rm H_2,O_2,N_2,CO,CO_2,CH_4,C_2$ (ethane, ethylene, acetylene), $\rm C_3$ (propane, propylene), and $\rm C_4$ (1-butene)

\*With methanizer bypass to improve low end performance.

# Ensure operational stability and generate reproducible data... day in and day out



This transformer gas analysis was completed in just 10 minutes. Trace levels of C0 and  $CO_2$  were analyzed by conversion to  $CH_4$ , followed by FID detection.



It took just 15 minutes to complete this transformer gas analysis per ASTM D3612-A. The analyzer was also configured with data reporting macros per the ASTM method. (Microsoft® Excel required.)

#### Put your applications on the fast track

Contact your local Agilent Representative or Agilent Authorized Distributor at **agilent.com/chem/contactus** 

Or call 800-227-9770 (in the U.S. or Canada)

Visit **agilent.com/chem/appkits** for a description of available Analyzers and Application Kits

This information is subject to change without notice.

© Agilent Technologies, Inc. 2014 Printed in U.S.A., September 25, 2014 5991-5045EN



**Agilent Technologies**