

UHPLC Expert eSeminar Curriculum

Agilent Liquid Phase Separation News from Pittcon 2010!



Christian Gotenfels
Product Manager UHPLC

- 1290 Infinity LC**
- 10 x Higher Sensitivity
 - Higher sample capacity
 - Total cost of ownership
 - Third party Software Control



Michael Frank
Product Manager HPLC

- Ultra-High Pressure Valve Solutions



Martin Vollmer
Product Manager SFC

- Supercritical Fluid Chromatography
Superior Sensitivity with
green chemistry



The world's most sensitive UV detector!



1290 Infinity Diode Array Detector




Agilent Technologies

History of Sensitivity Gain (DAD)

- *The last 30 years*



Sensitivity
Gain



Detector Model (Intro Year)




Agilent Technologies

What is Sensitivity?



Limit of detection (LOD):

The lowest concentration, or smallest mass flow, which can be distinguished from the noise by a certain predefined probability (Signal/Noise).




Limit of detection (LOD):

$$h_{\text{Signal}} = 3(2) \times h_{\text{Noise}}$$




Effects of path length increase

Conventional flow cells:




10 mm pathlength
13 μL geom. volume




60 mm pathlength
78 μL geom. volume

Peak Dispersion



- Loss of resolution
- Loss of signal height
- Increase of Noise


Max-Light High Sensitivity cell:



Optofluidic waveguides
(total internal reflection)




60 mm pathlength
4 μL σ_v dispersion volume



Agilent Technologies

Optofluidic waveguides – Agilent Max-Light cartridge cell utilize total-internal reflection in a non-coated fused silica fiber




[Check DAD video on the website](#)



11.4 x Higher Sensitivity

-1290 Infinity DAD compare to 1200 Series DAD SL



Columns: 150 x 4.6mm Zorbax SB C18, 5µm

Sample: Anthracene: 835 pg/µL

Mobile phase: A: Water, B: Acetonitrile

Elution: isocratic 80 % B

Injection volume: 5 µL

Flow: 1.5 mL/min


DAD: 251/4nm, Ref= 450/80nm,
2.5Hz, slit width 4nm

	1290 DAD 60 mm	1200 DAD SL 10mm
Height [mAU]	28.87579	4.93845
Noise [mAU]	0.009806	0.01908
Signal/ noise	2944	259
Sensitivity increase	+11.4	



11.6 x Higher Sensitivity

- 1290 Infinity DAD compare to 1200 Series VWD A/B



Columns: 150 x 4.6mm Zorbax SB C18, 5µm

Sample: Anthracene: 835 pg/µL

Mobile phase: A: Water, B: Acetonitrile

Elution: isocratic 80 % B

Injection volume: 5 µL

Flow: 1.5 mL/min


DAD: 251/4nm, Ref= 450/80nm,
2.5Hz, slit width 4nm

VWD: 251 nm, 2.5 Hz

	1290 DAD 60 mm	1200 VWD 10mm
Height [mAU]	28.87579	4.79998
Noise [mAU]	0.009806	0.01894
Signal/ noise	2944	253
Factor	+11.6	

1290 Infinity Diode Array Detector


- *Baseline robustness due to significant reduced RI effects*



Agilent Technologies

Comparison with other vendors

- Noise Specifications




¹ Dry flow cell

² 8 nm slit



Comparison with other vendors

- *Signal/Noise*



¹ Dry flow cell

² 8 nm slit



10 x Higher Sensitivity

- *Customer benefits*

- Higher data quality for more confidence
- Detection of lower levels of impurities (genotoxins)
- Simplified sample preparation methods (10 x less sample is required for extraction and clean up)
- Prepared to meet more stringent regulative requirements



1290 Infinity DAD Features & Benefits

- Ultra sensitivity through revolutionary Agilent Max-Light cartridge cell with 60 mm optical path length (typically noise: $< \pm 0.6 \mu\text{AU}/\text{cm}$)
- Universal Agilent Max-Light cartridge standard cell with 10 mm optical path length provide high sensitivity (Noise: $< \pm 3 \mu\text{AU}$) and lowest peak dispersion for 2.1, 3 and 4.6 mm ID columns
- More reliable and robust peak integration process due through less baseline drift
- Multiple wavelength and full spectral detection at high sampling rate of 160 Hz, keeping pace with the analysis speed of ultra fast LC
- Programmable slit from 1 to 8 nm provides optimum incident light conditions for rapid optimization of sensitivity, linearity and spectral resolution
- RFID technology for flow cells and lamp provide new levels of data traceability
- Electronic temperature control (ETC) provides maximum baseline stability and practical sensitivity under fluctuating ambient temperature and humidity conditions





Agilent 1290 Infinity LC Injector HTS/HTC

Make your high throughput
LC/MS analysis more robust



Agilent Technologies

1290 Infinity Sample Injection Choices



1290 Infinity Autosampler (G4226A)



- Up to 2 well plates cooled
- 40 ppm or 0.004% typical carryover
- Injection volume 40 μ L and up to 100 μ L possible*

1290 Infinity LC Injector HTC (G4278A) **NEW**



- 50 cm wide
- Up to 12 well plates cooled
- ~40 ppm or 0.004% typical carryover
- Injection volume up to 100 μ L with DLW, and up to 5 mL** with Active or Fast wash

1290 Infinity LC Injector HTS (G4277A) **NEW**



- 80 cm wide
- Up to 24 well plates cooled
- ~40 ppm or 0.004% typical carryover
- Injection volume up to 100 μ L with DLW, and up to 5 mL** with Active or Fast wash

* May 2010

** Separate HW required



Fast Dynamic Load & Wash (DLW)



DLW wash station

DLW pump (more robust performance)





Syringe




Fast Dynamic Load & Wash Cycle Demonstration



Cycle Start




Carryover Test UV



Injection of 1200 ng Chlorhexidine on column; followed by 2 blank injections using the new DLW option. Less than 0.003% (30ppm) of carryover could be detected.




Carryover, 1.0 μ M (+/-) Propranolol (260.2 pg on column) (+MRM: 260.2 \rightarrow 116.1)



Dwell time = 20 ms
Interchannel delay = 3.5 ms
Cycle time = 50 ms

Flow rate = 0.8 mL/min



Total cost of ownership

“All costs, including direct and indirect costs, associated with owning capital assets required to support your business operations”.

Total cost of ownership =

Process Factor

» Administrative & Management

Product Factor

» Product Cost

- Instrument price
- Operating cost
- Service cost
- Uptime
- etc

Productivity Factor

» UHPLC Productivity

- Speed
- Resolution
- Sensitivity



Agilent Technologies

UHPLC Productivity with HPLC service costs

- *Lowest total cost of ownership*

“UHPLC systems are seen as less reliable and more complex to maintain, resulting in higher service costs compare to HPLC systems”.

This has changed with 1290 Infinity LC, because of

1290 Infinity Design for supportability:

- Instrument monitoring
- Enhanced maintenance procedures
- Superior Diagnostic (onsite/remote)
- Increased lifetime of parts
- Easy and fast exchange of parts



Agilent Technologies




Agilent 1290 Infinity LC

**Third party
Software Control**

Instrument Control Framework


The Ultimate 3rd Party Connectivity Solution



Agilent Technologies

Instrument Control Framework

The Ultimate 3rd Party Connectivity Solution



Advantages for Customers

Full Control

- All actual and future modules and systems supported – 1120, 1200, 1290, ...
- All features supported – complete Agilent RC.net driver functionality as Agilent SW

Fast, Easy, Trouble-free Control

- Immediate and easy access new releases – via Plug and Play of new/updated Agilent RC.net Driver

Future Proof Control

- All future modules, systems and updates supported



Agilent Technologies

Dionex – Agilent Joined Statement

Dionex & Agilent Technologies
Joint Statement
on
Dionex Chromeleon Control of the Agilent 1290 Infinity LC

10th November, 2009

Dear Agilent 1290 Infinity LC and Dionex Chromeleon customer,

We have received a number of enquiries from joint customers regarding Dionex Chromeleon control of the Agilent 1290 Infinity LC.

Agilent and Dionex have a long standing agreement on cooperation in the development of software control of each other's instruments. The agreement provides Dionex access to the full set of Agilent 1200LC Series instrument control codes and/or drivers with assurance of continued updating of information should the codes/drivers change in any way and technical support from Agilent for instrument driver development and troubleshooting.

For the Agilent 1290 Infinity LC, the two companies are actively collaborating on developing Chromeleon control and early prototypes already exist. According to the current project plan, instrument control for the Agilent 1290 Infinity LC is planned to be available middle of 2010.

Regards



Dr. Stefan Schütte

Sr. Marketing Director Liquid Chromatography
Liquid Phase Separations Business/LPAD-LSG
Agilent Technologies R&D and Marketing GmbH & Co. KG
Hewlett-Packard-Strasse 8
76337 Waldbronn
Germany



Fraser McLeod

Director of Product Marketing
Life Sciences Business Unit
Dionex Softron GmbH
Domierstrasse 4
82110 Germering
Germany



Agilent Technologies

Agilent Customer Letter on ICF

Agilent Technologies

Status on use of ICF by Waters

Agilent is fully cooperating with Waters to enable 1290 Infinity LC control in Waters Empower®. Already in January 2009, a prototype ICF was delivered to Waters for evaluation together with the 1200 driver set that is currently used by Agilent ChemStation and EZChrom. In August 2009, an updated and backwards-compatible version of the ICF prototype, including 1290 support, was provided. Apart from the ICF activity, Agilent has also provided Waters with the 1290 Infinity LC Instrument low level control codes to enable optional extension of the current ICS 1200 LC driver set to control the 1290 Infinity LC. However, it is our understanding that Waters supports the ICF approach to provide 1290 control in Empower as it minimizes both their R+D and support investment and allows any future update on Agilent instrument drivers to be seamless and automatically integrated into Waters Empower® chromatographic data system at the same point in time when Agilent releases new or updated instruments.

Waters has explicitly requested that for information on the timing of availability of 1290 Infinity LC control in Waters Empower® users should contact their responsible Waters sales representative.

As control of 1290 by third-party software is work in progress we will provide regular updates on this letter to keep you fully informed.

Please let us know if additional clarification is required.



Dr. Stefan Schütte
Sr. Marketing Director Liquid Chromatography
Liquid Phase Separations Business/LPAD-LSG
Agilent Technologies R&D and Marketing GmbH & Co. KG
Hewlett-Packard-Strasse 8
76337 Waldbronn
Germany

Jan 2009: Waters received ICF from Agilent

Aug 8, 2009: Waters received backwards compatible 2nd revision of ICF, including 1290 support.

Aug 8, 2009: Waters received 1290 control codes (allows traditional ICS driver development)

2009: Agilent receives enhancement requests from Waters (e.g. Import of Empower plate formats)

Agilent's understanding is that implementation method of choice for 1290 is ICF, feasibility proven.



Agilent Technologies

Thanks for your attention!

