

KIT CONTENT

Order No.	Volume	Symbol	Component
ZV-3061-02S1-10	1 x 3 ml	STD I	Plasma Calibrator, lyophilised
ZV-3061-02K1-10	1 x 3 ml	CONT I	Plasma Control, lyophilised
ZV-3061-02K2-10	1 x 3 ml	CONT II	Plasma Control, lyophilised
ZV-3061-02R1-10	1 x 4 ml	R1	Reagent 1, Contains deuterated internal standard, lyophilised
ZV-3061-02R2-10	1 x 4 ml	R2	Reagent 2, Contains reducing agent
ZV-3061-02R3-10	1 x 40 ml	R3	Reagent 3, Contains base
ZV-3061-02R4-10	1 x 5 ml	R4	Reagent 4, Contains derivatisation reagent and organic modifier
ZV-3061-02R5-10	1 x 40 ml	R5	Reagent 5, Contains aqueous diluent
ZV-3061-02MP-10	1 x 1000 ml	MP	Mobile Phase A, Contains organic solvent
ZV-3061-02WB-10	1 x 500 ml	WB	Washing solution, Contains organic solvent
ZV-3061-KK-En_Rev01	1 x 1 pcs.	KK	User Guide