QUANT-MY-WAY: SIMPLIFY AND PERSONALIZE YOUR MASSHUNTER QUANTITATIVE ANALYSIS EXPERIENCE

A Faster Route to Insight

What is Quant-My-Way?

Quant-My-Way is a set of tools to designed to assist you in exposing only the capabilities of interest so you can create your own streamlined version of MassHunter Quantitative Analysis for each of your specific assay(s). You can customize the new Ribbon with just the actions you need, and you can hide as much or as little as your assay requires.

E = +	Agilent MassHunter Quant-My-Way Fl	avors Setu	ıp - Vanilla		? – 8 ×
Reload Setting: Robot Setting:					
Flavor Settings 《 General Settings Each: Table Settings Column Labels Mibbons Sample Types Outliers	Ribbons Close commands frame: A Adox Caunitative Analysis Access Asiasima Curve Actions Access Asiasima Curve Actions Act Calination Add Samples Add/more Columns Add/more Columns	• * * *	Add >> << Remove	Constraint the Ribbox: Ribbon Table Method Method Entit Method from Acquired MMID Esta Method from Acquired MMID Esta Method from Acquired MMID Esta Method from Acquired Adm Ara New Method from Acqu	
Flavor Root	Apply Target RTs to Qualifier			Update New Group New Button New Menu Properties.	~

Customize your Ribbon with only the features you need.

MORE DATA MORE POSSIBILITIES

TURN DATA INTO INSIGHT WITH MASSHUNTER SOFTWARE

Learn more: www.agilent.com/chem/ masshunterinsights



With Quant-My-Way, you can choose what features of MassHunter Quantitative Analysis are displayed.

	lent MassHu	nter Quantitative Ana Home	lysis (Vanilla) - Dri View	ugsOfAbuse - Quan	t-My-Way.batch.bi	n – 🗆 ×
File File New	Open Batch	Save Analyze		ienerate Report		
Batch Tab	le	✓ Sample Type: Sam			• 🔰 ISTD: Meth-	* #× 15 ∎ III
	Sampl	e	Amp Results	Cocaine Results	MDMA Results	Meth Results
Name	Approved	DA Method File	Final Conc.	Final Conc.	Final Conc.	Final Conc.
Sample-1		Quant-My-Way.m			6.4064	6.7843
Sample-2		Quant-My-Way.m	4.6270	4.3835	5.1292	5.8008
Sample-3		Quant-My-Way.m	13.6529	10.9739	10.7484	14.2957
+ MRM (150.1 -> \$\$ x10 ³ 0 1.2- 1.1- 1- 0.9-				2.239 min.		
0.8- 0.7- 0.6- 0.5- 0.4- 0.3- 0.2- 0.1-						

As part of Quant-My-Way, the user interface has been updated and streamlined. Of course, MassHunter Quantitative Analysis will continue to contain all the features you already know and love.



Updated and Streamlined User Interface

If you prefer the look and feel of earlier versions of MassHunter Quantitative Analysis, you can easily switch to the classic view.



The familiar navigation, colors and icons from earlier versions of MassHunter Quantitative Analysis will continue to be available.

Faster, easier, and more productive MS analysis from Agilent MassHunter www.agilent.com/chem/masshunterinsights

> For Research Use Only. Not for use in diagnostic procedures. Information, descriptions and specifications in this publication are subject to change without notice.

> > © Agilent Technologies, Inc. 2017 Published in USA, May 31, 2017 5991-8192EN

