

Agilent Technologies

<http://www.chem.agilent.com/layouts/agilent/downloadFirmware.aspx?whid=69761>

Firmware files for Agilent 1100/1200/1260/1290/1120/1220 series HPLC modules. Depending on the revision, the set may contain

- A.xx.xx (for modules without on-board LAN)
- B.xx.xx (for modules with on-board LAN)
- B.02.xx (for G4208A Instant Pilot) - for use with A.06.xx/B.01.xx/B.06.xx
- A.05.xx (for G4208A Instant Pilot) - for use with A.05.xx
- B.04.02 (for G1323B Control Module) - for use with A.06.xx/B.01.xx/B.06.xx
- B.03.22 (for G1323B Control Module) - for use with A.05.xx and below
- A.02.10 (for G1323A Control Module) - for use with A.05.xx and below

and corresponding documents (revision change information and compatibility matrices).

In general, the latest firmware is backward compatible series modules. Limitations are based on

- new features vs. control software and
- hardware (board type)

Please refer to the ***FWUpdate\_Infinity.PDF*** or ***FWUPDATE.PDF*** (part of ***Agilent HPLC Modules Firmware Update Tools & Procedures***).

The firmware update tools are not available in this set. Return to the overview page for LC modules firmware and select Agilent ***HPLC Modules Firmware Update Tools & Procedures*** for more information on the tools and procedures.

Select the required firmware set and download the required/all files and follow the instructions described in the ***FWUpdate\_Infinity.PDF*** or ***FWUPDATE.PDF***.

Use always the latest FW-Update Tools provided on this page.